

Product datasheet for **SC303009**

ECM2 (NM_001393) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	ECM2 (NM_001393) Human Untagged Clone
Tag:	Tag Free
Symbol:	ECM2
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

Fully Sequenced ORF: >SC303009 representing NM_001393.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGAAGATTGCAGTTTTGTTTTGTTTTTCTGCTTATCATTTTTCAAACGACTTTGAAAAAATGAA
GAAATTCCTAGGAAGCAAAGGAGGAAGATCTACCACAGAAGGTTGAGGAAAAGTTCAACCTCACACAAG
CACAGATCAAACAGACAGCTTGGAATTCAGCAAAACAGTTTTTACACCAGTAGCAAGACTTCCTATT
GTAACTTTGATTATAGCATGGAGGAAAAGTTGAATCCTTTTCAAGTTTTCTGGAGTAGAATCAAGT
TATAATGTGTACCAGGAAAAGAGGACACTGTTTGGTAAAGGGCATAACCATGTACAACAAAGCTGTG
TGGTCGCCTGAGCCCTGCACTACCTGCCTCTGCTCAGATGGAAGAGTTCTTTGTGATGAAACCATGTGC
CATCCCCAGAGGTGCCCCAAACAGTTATACCTGAAGGGGAATGCTGCCGGTCTGCTCCGCTACTGTC
TCCTATTCTCTACTCAGTGGTATAGCATTAAATGATAGAAATGAATTTTCTGGTGATTCTTCAGAACAA
AGAGAACCTACCAATTTACTTCATAAGCAACTGCCACCTCCTCAGGTGGGAATGGACCGAATAGTAAGA
AAAGAAGCACTTCAATCTGAGGAGGATGAAGAAGTGAAGAAGAAGATACAGAGCAAAAGAGAGAGACC
CCTGAATCTAGAAATCAGGGGCACTTTACAGTGAGGGGGACAGCAGAGGAGGAGACAGAAAGCAGAGG
CCTGGAGAGGAGAGGAGGCTGGCACACCAGCAACAACGCCAAGGAAGGGAGGAGGAGGAGGATGAGGAG
GAGGAGGGTGAGGAGGGTGAGGAGGATGAGGAGGACGAGGAGGACCCGGTAAGAGGAGATATGTTCCGA
ATGCCCTCTCGATCCCCGCTTCTGCTCCTCCCAGAGGCACACTGCGCCTGCCAAGCGGGTGCTCTCTG
TCCTACAGGACCATCAGCTGCATCAACGCCATGCTTACCAGATACCACCGCTGACAGCACCACAGATA
ACAAGTCTGGAGCTCACTGGCAATTCATCGCCTCCATCCCAGATGAAGCATTAAATGGATTACCAAT
TTGGAAGGCTTGATCTGAGTAAAAATAATCACTTCTTCAGGCATAGGTCCAAAGCATTCAAGCTT
CTGAAGAAGTTAATGCGTTTGAATATGGATGAAATAATTTGATACAGATTCCTTCACAATTGCCATCT
ACATTAGAAGAACTTAAAGTCAATGAGAACAATCTTCAGGCTATCGATGAAGAAAGTTTATCAGACTTA
AATCAGTTGGTCACCTTAGAATTGGAAGGAAACAATCTCAGTGAAGCCAATGTCAATCCTTTAGCTTTC
AAACCTTTGAAGAGCCTAGCCTACTTGCGTCTGGGAAAAAATAAATTTAGAATTATACCGCAGGGTCTT
CCTGGTTCTATTGAGGAATTATACCTAGAAAATAACCAAATTGAAGAAATAACTGAAATTTGTTTCAAT
CATACCAGAAAGATCAATGTCTTGTACTACGTTATAACAAAATTGAAGAAAATAGGATTGCTCCTTTA
GCCTGGATAAATCAAGAAAATCTAGAATCCATTGATCTCTCTACAACAAGCTCTATCACGTCCCGTCC
TATCTACCAAGTCCTTGCTGCACCTAGTACTCCTTGGGAACCAGATTGAACGGATCCCTGGCTATGTG
TTTGGCCACATGGAACCAGGCCTGGAATACTTGTACCTGTCATTTAACAACTTGCTGATGATGGCATG
GACCGTGTCTCCTTCTATGGGGCATATCATTCTCTGAGAGAATTTTCTGGATCACAATGACTTAAAA
TCTATACCACCTGGGATACAAGAAATGAAAGCACTACATTTTCTGAGGCTGAACAACAACAGATACGG
AACATTCTTCCAGAAGAAATTTGCAATGCTGAAGAGGATGATGACTCAAACTCTGGAACATCTTCACTT
GAAACAAATTATATTAATTAAGAGAAATACCATCTTACACATTTTATGCATAAGATCATACTCAAGT
ATCGTTCTTAACCACAAAACATCAAGTAA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
  
```

Restriction Sites: SgfI-MluI

ACCN: NM_001393

Insert Size: 2100 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001393.3</u>
RefSeq Size:	3307 bp
RefSeq ORF:	2100 bp
Locus ID:	1842
UniProt ID:	<u>O94769</u>
Cytogenetics:	9q22.31
Protein Families:	Secreted Protein
MW:	79.8 kDa
Gene Summary:	ECM2 encodes extracellular matrix protein 2, so named because it shares extensive similarity with known extracellular matrix proteins. Three transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Oct 2010] Transcript Variant: This variant (1) represents the longest transcript and encodes the longest isoform (1).