

Product datasheet for SC303211

NPTX1 (NM_002522) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	NPTX1 (NM_002522) Human Untagged Clone
Tag:	Tag Free
Symbol:	NPTX1
Synonyms:	NP1
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>SC303211 representing NM_002522. Blue=Insert sequence Red=Cloning site Green=Tag(s)


```
GCTCGTTTGTAGTAACCGTCAGAATTTTGTAAACGACTACTATAGGGCGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGCCGGCCGGCCGCGCGCGCGCACCTGTGCGCTGCTCGCCCTCTGCCTCCTGGGCGCCGGGGCCAG
GATTTTCGGGCCGACGCGCTTCATCTGCACCTCGGTGCCGTGGACGCCGACATGTGCGCCGCGTCCGTG
GCCCGCCGGCGCCGAGGAGCTCCGGAGCAGCGTGTGTCAGCTCCGCGAGACGGTGTGTCAGCAGAAG
GAGACCATCTGAGCCAGAAGGAGACCATCCGCGAGCTGACCGCAAGCTGGGCCGCTGCGAGAGCCAG
AGCAGCTGGACCCGGAGCCGGCAGGGCCGGCGGGCGGCCGCAAGCAGCCCGGCTCGGGCAAG
AACACCATGGGCGACCTGTCCCGGACACCGGCCCGGAGACGCTCAGCCAACTCGGGCAAACCTTTGCAA
TCGCTCAAACCCGCTGGAGAACCTCGAGCAGTACAGCCGCTCAATTCCTCCAGCCAGACCAACAGC
CTCAAGGATCTGCTGAGAGCAAGATCGATGAGCTGGAGAGGCAGGTGCTGTCCCGGGTGAACACCTG
GAGGAGGGCAAGGGGGGCCCCAGGAACGACACCGAGGAGGGTCAAGATCGAGACCGCCCTGACCTCC
CTGACACCGGATCAGCGAGCTCGAGAAAGGTCAGAAAGACAACCGCCCTGGAGACAAGTTCAGCTC
ACATTCCTGCGGACCAACTATATGTATGCCAAGGTGAAGAAGAGCCTGCCAGAGATGTACGCCCTTC
ACTGTCTGCATGTGGCTCAAGTCCAGCGCCACGCCAGGTGTGGGCACGCCCTTCTCCTACGCTGTGCC
GGCCAGGCAACGAGCTGGTCTCATTGAGTGGGGCAACAACCCCATGGAGATCCTCATCAATGACAAG
GTGGCCAAAGTTGCCTTTTGTCAATGATGGCAAGTGGCACCACATCTGTGTCACCTGGACACCCGG
GACGGGGTCTGGGAGCCCTACCAGGATGGCACGCAGGGTGGCAGTGGCGAGAATTGGCGCCCTATCAC
CCCATCAAGCCCCAGGGCGTGTGGTGTGGGCCAGGAGCAGGACACTCTGGGTGGTGGGTTTGATGCC
ACCCAGGCATTTGTGGGTGAGCTGGCCACTTCAACATCTGGACCGCAAGCTGACCCCCGGGGAGGTG
TACAACCTGGCCACCTGCAGCACCAAGGCTCTGTCCGGCAATGTCATCGCTGGGCTGAATCCCACATC
GAGATCTACGGAGGGGCCACCAAGTGGACCTTCGAGGCCGTGCGCCAGATCAACTGA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
```



[View online »](#)

Restriction Sites: SgfI-MluI

Plasmid Map:



ACCN: NM_002522

Insert Size: 1299 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_002522.3](#)

RefSeq Size: 5437 bp

RefSeq ORF: 1299 bp

Locus ID: 4884

UniProt ID: [Q15818](#)

Cytogenetics: 17q25.3

Protein Families: Secreted Protein

MW: 47.1 kDa

Gene Summary: NPTX1 is a member of the neuronal pentraxin gene family. Neuronal pentraxin 1 is similar to the rat NP1 gene which encodes a binding protein for the snake venom toxin taipoxin. Human NPTX1 mRNA is exclusively localized to the nervous system. [provided by RefSeq, Jul 2008]