

Product datasheet for **SC303303**

ZNF37A (NM_003421) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	ZNF37A (NM_003421) Human Untagged Clone
Tag:	Tag Free
Symbol:	ZNF37A
Synonyms:	KOX21; ZNF37
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)




[View online »](#)

Fully Sequenced ORF: >SC303303 representing NM_003421.
Blue=Insert sequence Red=Cloning site Green=Tag(s)

```
GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGATCACATCCCAGGGATCAGTGTGTTTTAGGGATGTGACTGTGGGCTTCACTCAAGAGGAGTGGCAG
CATCTGGACCCTGCTCAGAGGACCCTGTACAGGGATGTGATGCTGGAGAACTACAGCCACCTTGTCTCA
GTAGGGTATTGCATTCTAAACCAGAAGTGATTCTCAAGTTGGAGAAAGGCGAGGAGCCATGGATATTA
GAGGAAAAATTTCCAAGCCAGAGTCATCTGGAATTAATTAATACCAGTAGAAAATTTCAATAATGAAG
TTCAATGAGTTTAAACAAAGGTGGAAAAATGTTTCTGTGATGAAAAGCATGAAAATAATTCATTCTGAAGAG
GAACCTTCTGAATATAATAAAAAATGGGAACAGCTTCTGGCTGAATGAAGACCTCATTGGCATCAGAAA
ATTAATAATTGGGAACAATCTTTTGAATACAATGAATGTGGGAAAGCTTTCCTGAGAATCACTCTTC
CTTGTACATAAGAGAGGTTACACAGGACAGAAAACCTGCAATATACTGAACATGGGAAAACCTGTGAT
ATGTCATTTTTTCATCACTCATCAGCAAACATCCAAGAGAAAACCACTATGGTAATGAATGTGGAGAA
AATATCTTTGAGGAATCCATTCTCCTTGAACATCAGAGTGTACCATTAGCCAGAAGTTAAATCTC
ACTCCAATTCAGAGAACCCTCAATTAACAATATTATTGAATATAATGAGTGTGGAACATTTTTCACTG
GAAAAATTAGTCCTTCATTTACAACAGAGAACACATACAGGAGAAAAACCTTATGAATGTGATGAATGT
GGAAAAACCTTACCCAGAAAGTCAGCCACACAAGACATCAGAGAACACACAGGGGGAAAAACCTAT
GAATGTGATGAATGTGGGAAAGACCTTCTATAAGAAATTCAGACCTCATTAAACATCAAAGAATTCATAA
GGGAGAGACCTTATGGATGTCATGAATGTGGGAAATCCTTCAGTGAAAAGTCAACCTTACTCAACAT
CAAAGAACGCACACAGGGGAGAAACCATATGAATGTGATGAATGTGGGAAACCTTCTCATTAAAGTCA
GTCCTTACTGTGCATCAGAAAACACACAGGGGAGAAAGCCCTATGAATGCTATGCATGTGGGAAAGCC
TTTCTCAGAAAAATCAGACCTCATTAAACATCAAAGAATACACACAGGTGAAAACCTTATGAATGTAAT
GAATGTGGGAAAGTCATTCTCTGAGAAGTCAACCTTACTAAACATCTAAGAATCACAACAGGTGAGAAA
CCTTATGAATGTATTAGTGTGGGAAATTTTTCTGCTACTACTCCGTTTTACAGAACATCTGAGAAGA
CACACAGGGGAGAAACCTTTGGATGTAATGAATGTGGGAAACCTTCCGTGAGAAGTCAGCCCTAATT
GTTACCAGAGAACTCATATAAGACAGAAACCTATGGATGTAATCAATGTGGGAAATCATTCTGTGTG
AAGTCAAACCTCATTGCACATCATAGAACACACAGGGGAGAAACCTATGAATGTAATGTTTGTGGA
AAATCATTCTATGTTAAGTCAAACCTAACTGTACATCAGAGAATACACTGGGGAGAAACCTATAAAT
GTAGTAAACGAGGAAATTAATCTGGGTGA
ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGCGC
```

Restriction Sites: SgfI-MluI

Plasmid Map:


ACCN: NM_003421

Insert Size: 1686 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_003421.2](#)

RefSeq Size: 6584 bp

RefSeq ORF: 1686 bp

Locus ID: 7587

UniProt ID: [P17032](#)

Cytogenetics: 10p11.1

Protein Families: Transcription Factors

MW: 65.4 kDa

Gene Summary: May be involved in transcriptional regulation.[UniProtKB/Swiss-Prot Function]
Transcript Variant: This variant (2) lacks an internal segment in the 5' UTR, compared to variant 1. Variants 1-10 all encode the same isoform (a). Sequence Note: This RefSeq record was created from transcript and genomic sequence data to make the sequence consistent with the reference genome assembly. The genomic coordinates used for the transcript record were based on transcript alignments.