

Product datasheet for **SC303308**

ZNF189 (NM_003452) Human Untagged Clone

Product data:

| | |
|---------------------------|---|
| Product Type: | Expression Plasmids |
| Product Name: | ZNF189 (NM_003452) Human Untagged Clone |
| Tag: | Tag Free |
| Symbol: | ZNF189 |
| Mammalian Cell Selection: | Neomycin |
| Vector: | pCMV6-Entry (PS100001) |
| E. coli Selection: | Kanamycin (25 ug/mL) |




[View online »](#)

Fully Sequenced ORF: >SC303308 representing NM_003452.
Blue=Insert sequence Red=Cloning site Green=Tag(s)

```
GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGGCTTCCCCGAGCCCCCGCGGAGTCGAAGGGTTGCTGACATTTGAGGATGTGGCTGTGTTTTT
ACCCAGGAGGAGTGGGATTATCTGGACCCAGCTCAGAGAAGCCTGTATAAAGATGTCATGATGGAGAAT
TATGGAAACCTGGTCTCACTGGATGTTTTGAACAGAGATAAGGATGAGGAGCCAACGTAAAAACAAGAG
ATTGAAGAAATTGAGGAAGAAGTGAACACAGGGTGAATAGTTACAAGAATCAAAAGTAAAATTGAC
CAGGATCCTATGGGTAGAGAAACATTTGAACCTGTTGGTAGGTTAGATAAAACAAAGAGGGATCTTCTTA
TGGGAAATACCAAGGGAATCTTTGACCCAGGAACAGAGAATGTTTCAAGAAAACTAACATTATCCGT
AAAAGACCAAACCTCAGAAGAGAAATGCCATAAATGTGAAGAATGTGAAAAGGGTTTTGTCCGAAGGCC
CATTTTCATTCAACATCAAAGGGTCCATACTGGTGAGAAACCTTTTCAGTGAATGAATGTGGAAAAAGT
TTTAGTCGAGTTCATTTGTTATTGAACATCAGAGAATTCACACTGGGAAAAGGCCCTATGAGTGAAT
TACTGTGAAAAACCTTTAGTGTGAGCTCAACCCTTATTAGACATCAGAGAATCCCACTGGAGAAAAGA
CCCTATCAGTGAATCAGTGTAACAGAGCTTCAGCCAGAGAAGGAGCCTTGTTAAACATCAAAGGATT
CATAACAGGTGAGAAACCCATAAATGTAGTGACTGTGGAAAAGCCTTCAGTTGAAAATCACACCTTATT
GAGCATCAAAGAACTCACACTGGTGAAGAACCTTATCACTGTACCAAATGTAAGAAGAGCTTTAGTCGA
AATTCAATTGCTTGTGAGCATCAAAGAATTCACACTGGGAAAAGACCCCAATAAATGTGGTGAATGTGGG
AAAGCCTTTCGATTAAGCACATACCTTATAACAACCAAAAAATTCACACTGGCGAGAAGCCTTTTCTT
TGTATTGAGTGTGAAAAAGTTTCAGTCGGAGCTCATTCTTATTGAACATCAGAGGATCCATACTGGT
GAAAGACCTTATCAGTGCAAAGAGTGTGGAAAAGTTTCAGTCAGCTTTGCAACCTTACTCGTCATCAG
AGAATTCACACAGGAGACAAGCCCAATAAATGTGAGGAATGTGAAAAGCCTTTAGTAGAAGCTCAGGT
CTTATTACAGCATCAGAGAATTCACACCAGGAGAAGACTTATCCATACAATGAAACTAAGGAAAAGTTTT
GATCCAAATTGCAGTCTTGTATACAGCAGGAAGTCTACCCTAAGGAGAAATCTTATAAATGTGATGAA
TGTGGGAAAACCTTTAGTGTTAGTGCTCATCTGTACAACATCAAAGAATCCCACTGGTAAAAGCCC
TATCTATGTACTGTCTGTGGAAAAGCTTCAGCCGGAGCTCATTCTTATTGAACATCAGAGAATCCAC
ACTGGTGAGAGACCCTATCTGTGCAGACAGTGTGGAAAAGCTTTAGTCAGCTTTGTAATCTTATTCGA
CATCAGGGTGTTCACACAGGTAATAAACCATAAATGTGATGAATGTGAAAAGGCCCTTAGCCGGAAC
TCGGGTCTTATTCAGCATCAGAGAATACACACAGGAGAGAAAACCTTATAAGTGTGAGAAGTGCACAAA
AGTTTCAGTCAACAGCGAGTCTTGTCAACCATCAGAAGATCCATGCAGAGGTGAAAACCAAGAAACC
CATGAATGTGACGCTTGTGGTGAAGCCTTAAATGCCGATTTCTCTTATTACAGCATCAGAAATTGCAC
ACAGCATGGATGCAATAA
ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
```

Restriction Sites: SgfI-MluI

Plasmid Map:


ACCN: NM_003452

Insert Size: 1881 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_003452.3](#)

RefSeq Size: 3213 bp

RefSeq ORF: 1881 bp

Locus ID: 7743

UniProt ID: [O75820](#)

| | |
|--------------------------|--|
| Cytogenetics: | 9q31.1 |
| Domains: | KRAB, zf-C2H2 |
| Protein Families: | Transcription Factors |
| MW: | 73 kDa |
| Gene Summary: | <p>Kruppel-like zinc finger proteins such as ZNF189 contain a conserved stretch of 7 amino acids that connects a variable number of DNA-binding zinc finger repeats of the cys(2)his(2) (C2H2) type (summarized by Odeberg et al., 1998 [PubMed 9653648]). Approximately 30% of human Kruppel-like zinc finger proteins contain an N-terminal Kruppel-associated box (KRAB) domain. The KRAB domain consists of approximately 75 amino acids that may be subdivided into an A box, which is present in every KRAB domain and is essential for transcriptional repression, and a B box, which is not always present.[supplied by OMIM, May 2010]</p> <p>Transcript Variant: This variant (1) encodes the longest isoform (1).</p> |