

Product datasheet for **SC303316**

Synapsin III (SYN3) (NM_003490) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Synapsin III (SYN3) (NM_003490) Human Untagged Clone
Tag:	Tag Free
Symbol:	SYN3
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)




[View online »](#)

Fully Sequenced ORF: >SC303316 representing NM_003490.
Blue=Insert sequence Red=Cloning site Green=Tag(s)

```
GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGAATTTCTCCGGCGACGTCTCTCTGACAGCAGCTTCATGGCCAACCTGCCTAATGGCTATATGACG
GACCTGCAACGCCAGATAGCTCCACCAGCTCACCTGCTTCCCCGCCATGGAGAGGAGGCACCCCCAG
CCCTGGCTGCCCTTCTCTCTCCAGGATCCAGCCTTTTATAGCTCCCTCTCCAGTGCCATGAAGCAG
GCCCTCAGGCCACCTCAGGACTGATGGAGCCTCCAGGTCCTCCACGCCATTGTTCAAAGACCCAGG
ATCCTGTTGGTATCGATGATGCCATACAGACTGGTCAAGTATTTCCATGGGAAGAAGGTGAATGGA
GAGATTGAGATCCGAGTGGAGCAGGCTGAATTCTCAGAGTTGAACCTAGCTGCCTATGTGACCGGGGC
TGCATGGTGGACATGCAGGTCGTGAGAAATGGACCAAAGTGGTGAAGCAGATCCTCAAGCCAGACTTC
ATCCTGGTCCGCCAGCATGCCTACAGCATGGCCCTGGGGGAAGACTACCGCAGCCTGGTCATCGGCCTG
CAGTATGGAGGGCTGCCTGCTCAACTCTCTACTCCGTCTACAACCTCTGCAGCAAGCCTGGGTG
TTCTCTCAGCTCATTAGATCTTCCATTCCCTGGGTCTGAGAAGTTCCCGCTTGTGGAGCAAACATTT
TTCCCCAACCATAGCCAATGGTCACAGCCCCACACTTCCCGGTGGTAGCAAGCTGGGACATGCCAC
GCTGGAATGGGAAAGATCAAAGTGGAAAACAGCTTGACTTCCAGGACATCACCAGCGTGGTCGCCATG
GCCAAAACCTACGCCACCACCGAGGCCCTTCATCGACTCCAAGTACGACATCCGCATCCAGAAAATTGGA
TCCAACACTACAAGGCTTACATGAGAACCTCCATCTCTGGGAACTGGAAGGCCAACACAGGCTCTGCCATG
CTGGAGCAGGTGGCCATGACAGAGAGGTACAGGCTGTGGGTGGACAGCTGCTCGGAAATGTTTGGCGGC
CTGGACATCTGTCCGTCAAGGCTGTCCACAGCAAGGATGGCAGAGATTACATCATCGAGGTAATGGAC
AGCTCAATGCCGCTGATTGGAGAGCATGTGGAAGAGGACAGACAGCTGATGGCCGACCTTGTGTCTCC
AAAATGAGCCAGCTCCCGATGCCAGGAGGCACAGCGCCCTCCCCCTCAGACCTTGGGCTCCACAGATT
AAATCAGCGAAATCCCCAGGGCAAGCCAGCTGGGGCCTCAGCTAGGCCAGCCCCAGCCACGCCACCT
CCGCAAGGAGGCCCTCGCCAAGCTCAGTCTCTCAGCCCCAGAGATCTGGAAGCCCTCCCAACAGAGG
CTCTCCCCACAAGGCCAGCAGCCCTGAGCCCCAGTCCGGATCTCCACAGCAGCAAAGGTCAACAGGC
TCTCCGCAGCTATCCCGGCATCCAGTGGCAGCTCCCCAAACCAGGCTCCAAGCCAGGTGCCACCCTC
GCCTCACAGCCCCGGCCCCCTGTGCAGGGCCGTAGTACCTCCAGCAGGGTGAAGAGTCCAAGAAGCCA
GCACCACCCATCCGCATCTCAACAAATCTCAGTCCCTGACTAACAGCCTCAGCACATCCGACACCTCC
CAGCGTGGGACCCCAAGTGAAGACGAGGCCAAGGCTGAAACCATCCGCAACCTGAGGAAGTCTTTTGCC
AGCCTGTTCTCTGACTAA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
```

Restriction Sites: SgfI-MluI

Plasmid Map:


ACCN: NM_003490

Insert Size: 1743 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_003490.3](#)

RefSeq Size: 3047 bp

RefSeq ORF: 1743 bp

Locus ID: 8224

UniProt ID: [O14994](#)

Cytogenetics: 22q12.3

Protein Families:	Secreted Protein
MW:	63.3 kDa
Gene Summary:	<p>This gene is a member of the synapsin gene family. Synapsins encode neuronal phosphoproteins which associate with the cytoplasmic surface of synaptic vesicles. Family members are characterized by common protein domains, and they are implicated in synaptogenesis and the modulation of neurotransmitter release, suggesting a potential role in several neuropsychiatric diseases. The protein encoded by this gene shares the synapsin family domain model, with domains A, C, and E exhibiting the highest degree of conservation. The protein contains a unique domain J, located between domains C and E. Based on this gene's localization to 22q12.3, a possible schizophrenia susceptibility locus, and the established neurobiological roles of the synapsins, this family member may represent a candidate gene for schizophrenia. The TIMP3 gene is located within an intron of this gene and is transcribed in the opposite direction. Alternative splicing of this gene results in multiple splice variants that encode different isoforms. [provided by RefSeq, Oct 2008]</p> <p>Transcript Variant: This variant (IIIa) encodes the full-length isoform IIIa.</p>