

## Product datasheet for **SC303332**

### H3FL (HIST1H3B) (NM\_003537) Human Untagged Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	H3FL (HIST1H3B) (NM_003537) Human Untagged Clone
Tag:	Tag Free
Symbol:	H3C2
Synonyms:	H3/l; H3C1; H3C3; H3C4; H3C6; H3C7; H3C8; H3C10; H3C11; H3C12; H3FL; HIST1H3B
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>SC303332 representing NM_003537. Blue=Insert sequence Red=Cloning site Green=Tag(s)


```
GCTCGTTTAGTGAACCGTCAGAATTTGTAAACGACTCACTATAGGGCGCCGGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGGCTCGTACTAAACAGACAGCTCGGAAATCCACCGCGGTAAAGCGCCACGCAAGCAGCTGGCTACC
AAGGCTGCTCGCAAGAGCGCGCCGGCTACCGCGCGGTGAAAAGCCTCACCGTTACCGCCGGGCACT
GTGGCTCTGCGCGAGATCCGCCGTACCAAAAGTCGACCGAGTTGCTGATTGGAAGCTGCCGTTCCAG
CGCCTGGTGCAGAAATCGCCAAAGACTCAAGACCGATCTTCGTTCCAGAGCTCTGCCGTGATGGCG
CTGCAGGAGGCTTGTGAGGCTACTTGGTAGGGCTCTTTGAGGACACAACTTTGCCCATCCATGCT
AAGCGAGTGACTATTATGCCAAAGACATCCAGCTCGCTCGCCGATTTCGCGGAGAAAGAGCGTAA
ACGCGTACGCGCGCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCCGC
```

Restriction Sites: Sgfl-MluI



[View online »](#)

## Plasmid Map:



ACCN: NM\_003537

Insert Size: 411 bp

**OTI Disclaimer:** Due to the inherent nature of this plasmid, standard methods to replicate additional amounts of DNA in *E. coli* are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at [custsupport@origene.com](mailto:custsupport@origene.com) or by calling 301.340.3188 option 3 for pricing and delivery.

The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

<b>OTI Annotation:</b>	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<a href="#">NM_003537.3</a>
<b>RefSeq Size:</b>	472 bp
<b>RefSeq ORF:</b>	411 bp
<b>Locus ID:</b>	8358
<b>UniProt ID:</b>	<a href="#">P68431</a>
<b>Cytogenetics:</b>	6p22.2
<b>Protein Pathways:</b>	Systemic lupus erythematosus
<b>MW:</b>	15.4 kDa
<b>Gene Summary:</b>	Histones are basic nuclear proteins that are responsible for the nucleosome structure of the chromosomal fiber in eukaryotes. This structure consists of approximately 146 bp of DNA wrapped around a nucleosome, an octamer composed of pairs of each of the four core histones (H2A, H2B, H3, and H4). The chromatin fiber is further compacted through the interaction of a linker histone, H1, with the DNA between the nucleosomes to form higher order chromatin structures. This gene is intronless and encodes a replication-dependent histone that is a member of the histone H3 family. Transcripts from this gene lack polyA tails; instead, they contain a palindromic termination element. This gene is found in the large histone gene cluster on chromosome 6p22-p21.3. [provided by RefSeq, Aug 2015]