

Product datasheet for **SC303349**

BFSP2 (NM_003571) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	BFSP2 (NM_003571) Human Untagged Clone
Tag:	Tag Free
Symbol:	BFSP2
Synonyms:	CP47; CP49; CTRCT12; LIFL-L; PHAKOSIN
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>SC303349 representing NM_003571. Blue=Insert sequence Red=Cloning site Green=Tag(s)


```
GCTCGTTTGTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGCCGGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGAGTGAGAGGCGAGTGGTAGTGGACTTGCCACCAGTGCCAGCTCCAGCATGCCCTCCAGAGGCGC
AGGGCGTCCTTCAGGGGGCCACGGTCATCATCCTCCCTGGAGAGCCCCCAGCCTCCAGGACCAATGCC
ATGAGTGGCCTTGTCCGAGCACCCGGGGTCTATGTAGGAACAGCACCCAGTGGGTGCATAGTGGCTTG
GGTGCCCGTGTGACCCGCGGGCCCTCGGCATCAGCAGTGTCTTCTTCAGGGCTGCGGAGCTCAGGC
CTGGCCACCGTCCGGCTCCAGGTTTGGAGAGGACCATGGTGTCTTGAGGACCTAGGGGGCTGCCTG
GTGGAATATATGGCCAAAGTGCACGCCCTTGAGCAAGTCAGTCAGGAGCTGAAACACAACCTGCGGATG
CACCTGGAGAGCAAAGCCACACGCTCGGAAACTGGGGTGCCCTACGGGCTTCTGGGCCAGCAGCTGC
CAGCAGGTGGGTGAGGCAGTCTTGGAAAATGCCCGGCTCATGCTGCAGACAGAACTATCCAGGCCGGA
GCAGATGACTTTAAAGAGAGATATGAAAATGAGCAGCCATTTGAAAGGCAGCAGAAAGAGAAATTAAC
TCTCTGTATAAAGTCATTGATGAGGCTAATTTGACTAAAATGGACCTGGAGAGTCAAATAGAAAAGTCTG
AAAGAAGAACTTGGCTCTCTATCAAGAACTATGAAGAGGATGTGAAGCTGCTGCACAAACAGTTGGCA
GGGTGTGAGCTGGAACAAATGGATGCTCCCATTTGGCACTGGTCTGGACGACATCCTTGAGACGATCAGA
ATTCAGTGGGAGAGAGATGTTGAAAAGAACCAGGTGGAGGAGGAGCCCTGCTCCAAGCTAAGCAACAG
GCGGAGGTGGCCACATGTCCCAGACCCAGGAGGAGAAGCTGGCAGCTGCCCTCAGGGTGGAGTTACAC
AACACTTCGTGCCAAGTCCAGAGCTCCAGGCTGAGACAGAATCCTTACGTGCCCTGAAACGAGGCCCTG
GAGAACACCTTGACAGTGCACAGCAAGCACTGGCATGACATGGAGCTCCAGAACCTGGGCGCTGTGGTCGGC
CGGCTGGAGGCGGAGCTCAGGAAATCCGAGCGGAGGCGGAGCAGCAGCAACAGGAGCGCGCATCTG
CTGGCCCCAAGTGCCAGCTGCAGAAGGACGTGGCGTCTACCACGCCCTGCTGGACAGGGAGGAGAGC
GGCTGA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
```



[View online >](#)

Restriction Sites: SgfI-MluI

Plasmid Map:



ACCN: NM_003571

Insert Size: 1248 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_003571.3](#)

RefSeq Size: 1595 bp

RefSeq ORF: 1248 bp

Locus ID: 8419

UniProt ID: [Q13515](#)

Cytogenetics: 3q22.1

MW: 45.9 kDa

Gene Summary: More than 99% of the vertebrate ocular lens is comprised of terminally differentiated lens fiber cells. Two lens-specific intermediate filament-like proteins, the protein product of this gene (phakinin), and filensin, are expressed only after fiber cell differentiation has begun. Both proteins are found in a structurally unique cytoskeletal element that is referred to as the beaded filament (BF). Mutations in this gene have been associated with juvenile-onset, progressive cataracts and Dowling-Meara epidermolysis bullosa simplex. [provided by RefSeq, Jun 2009]