

## Product datasheet for **SC303376**

### Keratin 36 (KRT36) (NM\_003771) Human Untagged Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	Keratin 36 (KRT36) (NM_003771) Human Untagged Clone
Tag:	Tag Free
Symbol:	KRT36
Synonyms:	HA6; hHa6; KRTHA6
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

**Fully Sequenced ORF:** >SC303376 representing NM\_003771.  
 Blue=Insert sequence Red=Cloning site Green=Tag(s)


```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGGCCACCCAGACCTGCACCCCTACCTTCTCCACTGGGTCTATCAAGGGCCTCTGTGGCACAGCAGGC
GGCATCTCTCGGGTGTCTCCATCCGTTCTGTGGGCTCCTGCAGGGTCCCAGTCTCGCCGGTGCTGCA
GGGTACATCTCTTCTGCTAGGTCGGGCCTCTGGCCTTGGGAGCTGCTTGCCTGGCTCCTACCTGTCT
TCTGAGTGCCACACCTCTGGCTTTGTGGGAGCGGGGGCTGTTCTGCGAGGGCTCCTTCAACGGCAGC
GAGAAGGAGACTATGCAGTTCCTGAACGACCGCCTGGCCAACCTGGAGAAGGTGCGTCAGCTGGAG
CGGGAGAACGCGGAGCTGGAGAGCCGCATCCAGGAGTGGTACGAGTTTTCAGATCCCATACATCTGCCCA
GACTACCAGTCTACTTCAAGACCATCGAAGATTTCCAGCAGAAGATCCTGCTGACTAAGTCTGAGAT
GCCAGGCTGGTCTGCAGATTGATAATGCCAAGCTGGCTGCTGACGACTCCGGACCAAGTATGAGACA
GAGCTGTCTCTCGGCAGCTAGTGGAGGCCGACATCAACGGCCTGCGTAGGATCCTGGATGAGCTGACC
CTGTGCAAGGCTGACCTGGAGGCTCAGGTGGAGTCCCTGAAGGAGGAGCTGATGTGCCTCAAGAAGAT
CACGAGGAGGAAGTCAGTGTACTCCGTTGCCAACTTGGGACCGACTGAAATGTGGAGGTGGACGCTGCT
CCCCCAGTGGATCTCAACAAGATCCTGGAGGATATGAGATGCCAGTACGAGGCCCTGGTGGAGAATAAC
CGCAGAGATGTGGAGGCCTGTTCAACACCCAGACTGAGGAGCTGAACCAGCAGGTGGTGTCCAGCTCG
GAGCAGCTGCAGTGTGCCAGACGGAGATCATCGAGCTGAGACGTACGGTCAACGCGCTAGAGATTGAG
CTGCAGGCTCAGCACAGCATGCGGAATTCCTTGGAAATCCACCTGGCCGAAACCGAGGCCCGCTACAGC
TCCCAGCTGGCCAGATGCAGTGCCTGATCAGCAACGTGGAGGCCAGCTGTCTGAGATCCGCTGCGAC
CTGGAGCGGCAGAACCAGGAGTACCAGGTGTTACTGGACGTCAAGGCCCGGCTGGAGGGCAGATCGCT
ACCTACCGCCACCTGCTGGAGGGAGAGGACTGCAAGCTTCTCCCAACCTTGTGCCACGGCATGCAAG
CCTGTTATTAGAGTTCCTTCTGTCCCCCGGTGCCCTGTGTCCCCTCTGTGCCCTGCACCCCGGCTCCC
CAGGTTGGCACTCAGATCCGCACCATCACCGAGGAGATCAGAGATGGGAAAGTCATCTCCTCCAGGGAG
CACGTGCAGTCCCGCCCGCTGTGA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
  
```

**Restriction Sites:**

SgfI-MluI

**Plasmid Map:**



**ACCN:** NM\_003771

**Insert Size:** 1404 bp

<b>OTI Disclaimer:</b>	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
<b>OTI Annotation:</b>	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"> <li>1. Centrifuge at 5,000xg for 5min.</li> <li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li> <li>3. Close the tube and incubate for 10 minutes at room temperature.</li> <li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li> <li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li> </ol>
<b>RefSeq:</b>	<u><a href="#">NM_003771.4</a></u>
<b>RefSeq Size:</b>	1688 bp
<b>RefSeq ORF:</b>	1404 bp
<b>Locus ID:</b>	8689
<b>UniProt ID:</b>	<u><a href="#">O76013</a></u>
<b>Cytogenetics:</b>	17q21.2
<b>MW:</b>	52.2 kDa
<b>Gene Summary:</b>	The protein encoded by this gene is a member of the keratin gene family. This type I hair keratin is an acidic protein which heterodimerizes with type II keratins to form hair and nails. The type I hair keratins are clustered in a region of chromosome 17q12-q21 and have the same direction of transcription. [provided by RefSeq, Jul 2008]