

Product datasheet for **SC303384**

ADAM7 (NM_003817) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	ADAM7 (NM_003817) Human Untagged Clone
Tag:	Tag Free
Symbol:	ADAM7
Synonyms:	ADAM-7; ADAM 7; EAPI; GP-83; GP83
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

Fully Sequenced ORF: >SC303384 representing NM_003817.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCCGCATCGCC
ATGCTTCCCGGGTGTATATTCTTGATGATTTTACTCATTCTCAGGTTAAAGAAAAGTTCATCCTTGGG
GTAGAGGGTCAACAACCTGGTTCGTCCTAAAAAGCTTCTCTGATACAGAAGCGAGATACTGGACACACC
CATGATGATGACATACTGAAAACGTATGAAGAAGAATTGTTGTATGAAATAAACTAAATAGAAAAACC
TTAGTCTTCATCTTCTAAGATCCAGGGAGTTCCTAGGCTCAAATTACAGTGAAACATTCTACTCCATG
AAAGGAGAAGCGTTCACCAGGCATCCTCAGATCATGGATCATTGTTTTTACCAAGGATCCATAGTACAC
GAATATGATTCAGCTGCCAGTATCAGTACGTGTAATGGTCTAAGGGGATTCTTCAGAATAAACGACCAA
AGATACCTCATTGAACCACTGAAATACTCAGATGAGGGAGAACATTTGGTGTCAAATAAACCTGAGG
GTGCCGTATGGTGCCAATTATTCCTGTACAGAGCTTAATTTTACCAGAAAACTGTTCCAGGGGATAAT
GAATCTGAAGAAGACTCCAAAATAAAAGGCATCCATGATGAAAAGTATGTTGAATTGTTCTATTGTTGCT
GATGATACTGTGTATCGCAGAAATGGTCATCCTACAATAAACTAAGGAACCGAATTTGGGGAATGGTC
AATTTTGTCAACATGATTTATAAAAACCTTAAACATCCATGTGACGTTGGTTGGCATTGAAATATGGACA
CATGAAGATAAAATAGAATAATTCAAATATAGAACTACCTTATTGCGTTTTTCAATTTGGCAAGAA
AAGATCCTTAAACACGGGAAGGATTTTGTATCATGTTGTATTACTCAGTGGGAAGTGGCTCTACTACAT
GTGCAAGGAATTTCTTATCCAGGGGGTATGTGCCTGCCCTATTATTCCACCAGTATCATTAAAGGATCTT
TTACCTGACACAAAACATAAATGCAAACAGAATGGCACATCAACTGGGGCATAACCTTGGGATGCAGCAT
GACGAGTTCCTATGCACCTGTCTTCAGGAAAAATGCGTGATGGACAGTATGGAAGCATTCTGCACTG
AAATTCAGTAAATGCAGCCAAAACCAATACCACCAGTACTGAAGGATTATAAGCCAACATGCATGCTC
AACATTCATTTCTTACAATTTTCATGATTTCCAATTTTGGGAAACAAGAAGTTGGATGAGGGTGAA
GAGTGTGACTGTGCCCTGTCCAGGAGTACTAATCCTTGCTGTGATGCACACACATGTGTAAGG
CCAGGATTTACTGTGCAGAAGGAGAATGCTGTGAATCTTGTCAGATAAAAAAAGCAGGGTCCATATGC
AGACCGCGCAAGATGAATGTGATTTTCTGAGATGTGCACTGGCCACTCGCCTGCCTGTCTAAGGAC
CAGTTCAGGGTCAATGGATTTCTTGAAGAACTCAGAAGGCTACTGTTTCATGGGAAATGTCCAAC
CGTGAGGATCAGTGTCTGAACTATTTGATGATGAGGCAATAGAGAGTATGATATCTGCTACAAGATG
AATACAAAAGGAAATAAATTTGGATACTGCAAAAACAAGGAAAACAGATTTCTCCCTGTGAGGAGAAA
GATGTCAGATGTGAAAGATCTACTGCACTGGAGGGGAGCTTCTCTCTCTTGGAGAAGACAAGACT
TATCACCTAAGGATCCCAAGAAGAATGCTACTGTCAAATGCAAACTATTTTTTATACCATGATTCT
ACAGACATTGGCCTGGTGGCGTCAGGAACAAAATGTGGAGAGGGAATGGTGTGCAACAATGGTGAATGT
CTAAACATGGAAGGTCTATATCTCAACCAATTGCCCTCTCAGTGCAATGAAAATCCTGTGGATGGC
CACGGACTCCAGTGCCACTGTGAGGAAGGACAGGCACCTGTAGCCTGTGAAGAAACCTTACATGTTACC
AATATCACCATCTTGGTTGTTGTGCTTGTCTGTTATTGTCGGTATCGGAGTCTTATACTATTAGTT
CGTTACCGAAAATGTATCAAGTTGAAGCAAGTTCAGAGCCACCTACAGAAACCTGGGAGTGGAGAAC
AAAGGATACTTTGGTGTGAGCAGCAGATAAGGACTGAGCCAATCCTGCCAGAAATTCATTTCTAAAT
AAACCTGCAAGTAAAGATTCAAGAGGAATCGCAGATCCCAATCAAAGTCCAAGTGA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
  
```

Restriction Sites: SgfI-MluI

ACCN: NM_003817

Insert Size: 2265 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_003817.3
RefSeq Size:	3397 bp
RefSeq ORF:	2265 bp
Locus ID:	8756
UniProt ID:	Q9H2U9
Cytogenetics:	8p21.2
Protein Families:	Druggable Genome, Transmembrane
MW:	85.7 kDa
Gene Summary:	This gene encodes a member of the ADAMs family of zinc proteases. These transmembrane proteins play roles in multiple processes including cell signaling, adhesion and migration. The encoded protein lacks protease activity and may play roles in protein-protein interactions and cell adhesion processes including sperm-egg fusion. Mutations in this gene may be involved in the progression of melanoma. [provided by RefSeq, Oct 2011]