

Product datasheet for **SC303385**

Matrilin 4 (MATN4) (NM_003833) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Matrilin 4 (MATN4) (NM_003833) Human Untagged Clone
Tag:	Tag Free
Symbol:	MATN4
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)




[View online »](#)

Fully Sequenced ORF: >SC303385 representing NM_003833.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCCGGATCGCC
ATGAGAGGCCTTCTTTGCTGGCCGTGTTGCTGCTCCTTCTTCAGCCCTGGGAAACCCAGCTCCAGTTG
ACAGGTCCCAGGTGTACACTGGGCCCTGGATCTGGTGTTCGTGATTGACAGCTCCCGCAGCGTGCGC
CCTTTCGAGTTCGAGACCATGCGGCAGTTCCTCATGGCCCTCCTCCGAGGCCTGAACGTGGTCCCAAC
GCCACGCGCGTTGGCGTGATCCAGTATTCGAGTCAAGTGACAGAGCGTCTTCCCTCTCCGCGCGTTCTCT
CGCCGCGAGGACATGGAGCGCGCATCCGCGACCTGGTGCCTCTGGCGCAAGGCACCATGACGGGACTG
GCAATCCAGTACGCCATGAACGTGGCCTTCAGTGTGGCCGAGGGCGCGGACCGCCAGAGGAGCGCGTG
CCGCGTGTGCTGTATCGTGACAGACGGGCGGCCCCAGGACCGCGTGGCCGAGGTGGCGGCACAGGCG
CGCGCCCGCGGATTGAAATTTACGCGTGGGGTGCAGCGCGGACGTGGGCTCCCTGCGGCCATG
GCATCGCCCCGCTAGACGAGCAGTCTTCTCGTAGAGTCTTCGACCTCATCCAGGAGTTCGGCCTG
CAGTTCAGAGCCGGCTGTGTGCCATTGATCTGTGTGCTGAAGGGACCCATGGATGTGAGCACCCTGC
GTCAATCCCCAGGCTCCTATTTCTGCTACTGCCAAGTTGGCTTTGTACTCCAGCAGGACCAGAGGAGC
TGACGGGCCATTGACTACTGCAGCTTTGGGAACCATAGCTGTCAGCATGAGTGTGTTAGCACCCCTGGT
GGGCCACGGTGCCACTGCAGAGAGGGCCATGACTTGCAGCCTGATGGGAGGAGCTGTGAGTCCGGGAC
CTTTGCAATGGCGTGGACCATGGCTGTGAGTTCAGTGTGTGAGCGAGGGCCTCTCTACCGCTGCCTG
TGCCCCGAGGGGCGGCAACTTCAGGCAGATGGCAAGAGCTGCAACCGGTGCCGGGAAGGCCACGTGGAC
CTTGTCTGCTGGTTGATGGCTCCAAGAGCGTGCCTCCACAAAATTTCGAGCTAGTGAAGCGCTTCGTG
AACCAGATTGTGGACTTCTAGATGTGTCCCCGAGGGCACGCGGTGGGGTGGTGCAGTTCCTCGAGC
CGCGTGCACCGAGTTCCTCTGGTGCCTACGGCACCGCAGCCGAGGTGAAGCAGGCGGTCTTGCC
GTGGAGTACATGGAACGCGGCACCATGACAGGGCTGGCGTTGCGGCACATGGTGGAGCACAGCTTCTCC
GAGGCGCAGGGTGCACGGCCCCGTGCCCTAACGTGCCTCGTGTGGCCTGGTCTTCACGGATGGCCGC
TCCCAGGATGACATCTCGGTGTGGGCAGCGCGCCAAGGAGGAAGGCATCGTCATGTACCCGTGGGC
GTGGGCAAGGCGGTGGAGGCGGAGCTGCGCGAGATCGCCTCGGAGCCAGCGGAACTGCACGTGTCTAT
GCCCGGACTTCGGCACCATGACGCACCTGCTGGAGAACCTCAGAGGCAGCATCTGTCCAGAGGAGGGC
ATCAGCGCAGGGACAGAGCTTCGGAGCCATGCGAATGCGAAAGCCTCGTGGAGTTCAGGGCCGCAGC
CTGGGGGCGCTCGAGAGCCTGACGCTGAACCTGGCCAGCTGACGGCGCGCTGGAGGATCTGGAGAAC
CAGCTGGCCAACCAAGTGA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTAAACGGCCGGC
  
```

Restriction Sites: SgfI-MluI

Plasmid Map:


ACCN: NM_003833

Insert Size: 1746 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_003833.4](#)

RefSeq Size: 2335 bp

RefSeq ORF: 1746 bp

Locus ID: 8785

Cytogenetics: 20q13.12

Protein Families: Secreted Protein

MW: 64.1 kDa

Gene Summary: This gene encodes a member of von Willebrand factor A domain-containing protein family. The proteins of this family are thought to be involved in the formation of filamentous networks in the extracellular matrices of various tissues. This family member is thought to be play a role in reorganizing and regenerating the corneal matrix in granular and lattice type I dystrophies. It may also be involved in wound healing in the dentin-pulp complex. Alternative splicing results in multiple transcript variants. [provided by RefSeq, May 2013]
Transcript Variant: This variant (1) represents the longest transcript and encodes the longest isoform (1).