

Product datasheet for SC303439

KRT33A (NM_004138) Human Untagged Clone

Product data:


Product Type:	Expression Plasmids
Product Name:	KRT33A (NM_004138) Human Untagged Clone
Tag:	Tag Free
Symbol:	KRT33A
Synonyms:	Ha-3I; HA3I; hHa3-I; K33A; Krt1-3; KRTHA3A
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>SC303439 representing NM_004138. Blue=Insert sequence Red=Cloning site Green=Tag(s)

```
GCTCGTTTGTAGTAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGTCTTACAGTTGTGGCCTGCCAGCCTGAGCTGCCGCACCAGCTGCTCCTCCCGGCCCTGTGTGCC
CCCAGCTGCCACGGCTGCACCCTGCCGGGGCCTGCAACATCCCCGCAATGTGAGCAACTGCAACTGG
TTCTGTGAGGGCTCCTTCAATGGCAGTGAGAAGGAGACCATGCAGTTCTGAACGACCGCTGGCCAGC
TACCTGGAGAAGGTGCGTCAGCTGGAGCGGGACAACGCGGAGCTGGAGAACCTCATCCGGGAGCGGTCA
CAGCAGCAGGAGCCCTTGGTGTGTGCCAGCTACCAGTCTACTTCAAGACCATTGAGGAGCTCCAGCAG
AAGATCCTGTGCAGCAAGTCTGAGAATGCCAGGCTTGTGGTGCAGATCGACAATGCCAAGCTGGCCTCA
GATGACTTCAGGACCAAATATGAGACCGAGCTGTCCCTGCCGCAGCTGGTGGAGTCGGACATCAATGGC
CTGCGCAGGATCCTGGATGAGCTGACCCTGTGCAGGTCTGACCTGGAGGCCAGGTGGAGTCCCTGAAG
GAGGAGCTGCTGTGCCTCAAGCAGAACCATGAGCAGGAGGTTAACACCCTGCGCTGCCAGCTTGGAGAC
CGCCTCAACGTGGAGGTGGACGCTGCTCCCACTGTGGACCTGAACCAGGTCTGAATGAGACCAGGAGT
CAGTATGAGGCCCTGGTGGAAACCAACCGCAGGGAAGTGGAGCAATGGTTCGCCACGCAGACCGAGGAG
CTGAACAAGCAGGTGGTATCCAGCTCGGAGCAGCTGCAGTCTACCAGGCGGAGATCATCGAGCTGAGA
CGCACGGTCAATGCCCTGGAGATCGAGCTGCAGGCCAGCACAACCTGCCGAGACTCTCTGGAAAACAG
CTGACAGAGAGCGAGGCCCGCTACAGCTCCAGCTGCCAGTGTCCAGGTGCAGAGACTGATCACCACGTGGAG
TCCCAGCTGGCGGAGATCCGCAAGTACCTGGAGCGGCAGAACAGGAGATCAGGTGCTGCTGGAGCTG
CGGGCGCGGCTGGAGTGTGAGATCAACACGTACCGGAGCCTGCTGGAGAGCGAGGACTGCAAGCTCCCC
TCCAACCCTGCGCCACAACCAATGCATGTGACAAGTCCACTGGGCCCTGTATCTCTAATCCCTGTGGC
CTACGTGCTCGGTGTGGCCTTGAACACATTTGGGTACTAG
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
```

Restriction Sites: Sgfl-Mlul



[View online >](#)

Plasmid Map:


ACCN: NM_004138

Insert Size: 1215 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_004138.3](#)

RefSeq Size: 1287 bp

RefSeq ORF: 1215 bp

Locus ID: 3883

UniProt ID: [O76009](#)

Cytogenetics: 17q21.2

MW: 45.9 kDa

Gene Summary: This gene encodes a member of the keratin gene family. This gene is one of multiple type I hair keratin genes that are clustered in a region of chromosome 17q12-q21 and have the same direction of transcription. As a type I hair keratin, the encoded protein is an acidic protein which heterodimerizes with type II keratins to form hair and nails. There are two isoforms of this protein, encoded by two separate genes, keratin 33A and keratin 33B. [provided by RefSeq, May 2012]