

Product datasheet for **SC303487**

FOXD2 (NM_004474) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	FOXD2 (NM_004474) Human Untagged Clone
Tag:	Tag Free
Symbol:	FOXD2
Synonyms:	FKHL17; FREAC-9; FREAC9
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)


Fully Sequenced ORF: >SC303487 representing NM_004474.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAATACGACTCACTATAGGGCGGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGACCCTGGGCAGCTGCTGCTGCGAGATCATGTCTCCGAGAGCTCCCCGGCCGCGCTGTCCGAGGCC
GACGCAGACATAGACGTGGTGGCGCGCGCAGCGGGGGGAGCTCCCAGCTCGCTCCGGCCCCCGC
GCCCCCGGGACGTGCTCCCCACGGCCACGAGCCTCCCGCGGAGGAAGCCGAGGCAGACTTAGCCGAG
GACGAGGAGGAGTCTGGTGGCTGCTCGGACGGCGAGCCCCGCGCTCTGGCGTCCCGGGGGCGCGGCC
GCAGCGGGGAGCCCGGGGCCAGGCAGCGCGCGCGGCCCGCGGCGCAGCGGGGCCGGGGCCGGGACCCTG
TCGGGGGGCGCGGCGACGCGGAGCCCGTGGTGAAGCCGCCCTACTCGTACATCGCGCTCATCACCATG
GCCATCTGCAGAGCCCCAAGAAGCGGCTGACGTTGAGCGAGATCTGCGAGTTCATCAGCGGCCGCTTC
CCCTACTACCGGGAGAAGTTCGCCCTGGCAGAACAGCATCCGCCACAACCTCTCTCAACGACTGC
TTCGTAAGATCCCCCGGAGCCGGCAACCCGGGCAAGGGCAACTACTGGACGCTGGACCCGGAGTCG
GCCGACATGTTTCGACAACGGCAGCTTCTCGCGGCTCGAAGCGCTTCAAGCGGACGCCCTGCCGCG
CCGACCCACACCCGACCCCTACCCGGAGCTGCTGCTGCGTGGCGGGCCGCGCGGCGGGGGATCCC
GGCGCTTTCCTGCCCGGCTTCGCTGCCTACGGCGCTACGGCTACGGCTACGGGCTGGCTCTCCCGCC
TACGGCGCACCCCGCGGGGGCCGGCCCCGCATCCGCACCCGCACCCGCACGCCTTCGCTTTCGCCGCG
GCAGCCGCGCGCTCCTTGCCAGCTGTGCGTACCCCCAGGCCGCGCCGCGCGCTCCACCCGGACCT
CCGACGGCCTCGGTGTTTCGACGGCGCGGATCGGCCCCAGCTCCTGCGCCTGCCTCAGGCTCGGGCCG
GGCCCGGGCCCCGAGCCCTGCCCGCTTCTGGGCGCGAGCTGGGTGCGCCAAAGCCTTCTACGCG
GCGTCCCTGAGTCTCCCGCAGCCGGCACCCGCGCGGTCTGCCACCCGACTTCTGCGCCAGGCCTC
AAGACGGACCGGGCGGTGGTGCAGGCGCGGGGGCCGGGGCAGGGCAGAGGCCTTCTCTCTATA
GACCACATCATGGGCCACGGTGGCGCGGGGACGACCCCGGGCGCGCGGAGGCTCTCCGGGACCG
CCATTGCGGCGAGCCGCGGTCTGGGGCCAAAGCCAGGTCTTGCCATGCTGACTGCTCCGGCCCTG
GCTCCGTTGCTGGCCACATTGCCTCTCGCATCCCGGGACGCGCTGCTGCTCAGGGTCCCGGTTT
GCCAGCAAAGTCGCCGCTTAGTGGCTGCCACTTCTGA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGCG
  
```

Restriction Sites: SgfI-MluI

Plasmid Map:



ACCN: NM_004474

Insert Size:	1488 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_004474.3
RefSeq Size:	4675 bp
RefSeq ORF:	1488 bp
Locus ID:	2306
UniProt ID:	O60548
Cytogenetics:	1p33
Protein Families:	Druggable Genome, Transcription Factors
MW:	48.7 kDa
Gene Summary:	This gene belongs to the forkhead family of transcription factors which is characterized by a distinct forkhead domain. The specific function of this gene has not yet been determined. [provided by RefSeq, Jul 2008]