

Product datasheet for **SC303552**

dystrophia myotonica containing WD repeat motif (DMWD) (NM_004943) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	dystrophia myotonica containing WD repeat motif (DMWD) (NM_004943) Human Untagged Clone
Tag:	Tag Free
Symbol:	dystrophia myotonica containing WD repeat motif
Synonyms:	D19S593E; DMR-N9; DMRN9; gene59
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

Fully Sequenced ORF: >SC303552 representing NM_004943.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCCGGATCGCC
ATGGCGCGCGGCGCGGAGGGCGGCTCGGGCCCCGGCGCCCATGGGGGACTGCGCGGAGATTAAG
TCGCAATTCGCGACGCGGAGGGTTTCTACAAGCTACTCCCGGGGACGCGCGCCGCTCGCAGGTCGGGT
CCGGCTTCGCGCCAGACTCCGGTGCCGCCTCAGCCACCGCAGCCCCCGCGGCCCTGCCTCCGCCTCC
GGTCCCGGCGCTGCAGGCCCGCGCTCGTCCCGCGCGCCGACGGCCCGGACCCGGGCGCCCTGCC
GCCGTGCGCTCAGCCTCGTGCGCTCGGGGAGCCGGACAGCGCGGGGCGGGGAGCCGCCCGCCACG
CCCGCGGGGCTGGGCTCGGGGGAGACCGCTGCTTCAACTTGGGCGTGAGCTCTATTTCTACCCA
GGCTGCTGTCGTGGGAGCCAACGGTCCATTGACCTCAACAAGCCAATTGACAAGCGGATCTACAAG
GGCACCCAGCCACCTGCCACGATTTCAACCAGTTCAGTCTGCCACCGAGACCATCTCGTCTGGTG
GGCTTCTCAGCGGTCAAGTGCAGTACCTGGATCTCATAAAAAGGACACCAGCAAGCTGTTCAATGAG
GAGCGGTTGATCGACAAGACCAAGTGACATATCTGAAGTGGCTGCCTGAGTCGGAGAGCCTGTTCTG
GCATCACACGCCAGTGGCCACCTGTACCTGTACAACGTACGCCACCCCTGCGCCTCGGCCCGCCCCAG
TACAGCCTGCTGAAGCAGGGCGAGGGCTTCTGTCTATGCTGCCAAGAGCAAGGCACCCCGCAACCCG
CTGGCCAAGTGGGCGGTGGTGAGGGGCCCTCAACGAGTTTCGCCTTCTCGCCGATGGCCGGCACCTG
GCCTGTGTGAGCCAGGATGGTGCCTGCGCGTCTTCCACTTCGACTCCATGCTCCTGCGTGGGCTCATG
AAGAGCTACTTTGGGGCCTGCTGTGTGTGTGCTGGAGCCCTGACGGCCGCTACGTGGTGACGGGTGGC
GAAGATGACCTGGTCACCGTGTGGTCTTACCGAGGGCCGCGTGGTGGTTCGAGGCCATGGCCACAAG
TCCTGGGTCAACGCTGTGGCCTTTGACCCTACACCACAAGGGCAGAGGAGGCGGCGACAGCAGCCGGT
GCTGATGGGGAGCGGAGCGCGAAGAGGAGGAGGAGGCCAGGCTGCGGGCACAGGCTCGGCCGGG
GGCGCCCGCTCTCTCCACTGCCAAGGCTGGCTCCATTACTTACCGCTTTGGTCTCGGCGGGCCAGGAC
ACGCAGTCTGCCTGTGGACCTCACTGAAGACGTGCTTACCCGCACCCCCCTGGCCCGCACCCGC
ACCCTCCCTGGCACACCTGGCACACGCCACCGCCGCCAGCAGCTCGAGGGGTGGCGAGCCTGGCCCA
GGCCCCCTGCCTCGCTCGCTGTCCCGCTCAACAGTCTCCCGCACCCAGCTGGCGGGGCAAGCGGGC
GGCCCGGGTGTGGCGGCAGAGCCTGGCACACCATTAGCATTGGCCGCTTCGCCACGCTCACACTGCAG
GAGCGGCGGACCGGGGGCAGAGAAGGAGCACAAGCGCTACCACAGCCTGGGCAACATCAGCCGGGT
GGCAGTGGCGGCAGTGGCAGTGGTGGGAGAAGCCAGCGCCCTGTTCCCGCAGCCGCTGGACCCC
GCCAAGGTGCTGGGCACTGCGCTGTGCCGCGCATCCACGAGTGCCTGCTGGAGCCCTTGTGTGC
AAGAAGATCGCCAGGAGCGGCTCACAGTCTCCTGTTCTGAGGACTGCATCATCACTGCCTGCCAG
GAGGGCCTCATCTGCACCTGGGCCCGGGCGGGCAAGGCGTTCACAGACGAGGAGACCGAGGCCAGACA
GGGAAGGAAGTTGGCCAGGTACCCAGCAAGTCAGTGGTAGAGGGCATCTCTCCAACCAGGCAAC
TCCCCGAGTGGCACAGTGGTGTA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
  
```

Restriction Sites: SgfI-MluI

ACCN: NM_004943

Insert Size: 2025 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_004943.1</u>
RefSeq Size:	3308 bp
RefSeq ORF:	2025 bp
Locus ID:	1762
UniProt ID:	<u>Q09019</u>
Cytogenetics:	19q13.32
MW:	70.4 kDa