

Product datasheet for **SC303570**

RBM1A1 (NM_005058) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	RBM1A1 (NM_005058) Human Untagged Clone
Tag:	Tag Free
Symbol:	RBM1A1
Synonyms:	RBM; RBM1; RBM2; RBM1A; RBM1B; RBM1C; YRRM1; YRRM2
Vector:	<u>pCMV6 series</u>
Fully Sequenced ORF:	>NCBI ORF sequence for NM_005058, the custom clone sequence may differ by one or more nucleotides ATGGTAGAAGCAGATCATCCTGGCAAGCTTTTCATTGGTGGCCTCAATAGAGAAACCAAT GAGAAGATGCTTAAAGCAGTATTTGGGAAACATGGTCCCATATCAGAAGTTCTTTTGATA AAGGATCGAACCAGCAAATCCAGAGGCTTTGCATTTATTACTTTTGAGAACCCTGCAGAT GCTAAGAATGCTGCCAAAGATATGAATGGAAAGTCTTTGCATGGAAAAGCAATAAAAGTA GAACAAGCCAAGAAACCATCTTTTCAAAGTGGTGGTAGGCGGAGACCACCAGCTTCTTCG AGAAACAGAAGCCCTTCAGGAAGTCTGAGATCTGCAAGAGGAAGCCGTGGAGGAACAAGA GGGTGGCTTCCCTCACATGAAGGGCACCTGGATGATGGTGGATACACTCCTGATCTCAAG ATGAGTTATTCTAGGGGACTCATTCCAGTAAAAGAGGTCCATCTTCAAGAAGTGGAGGT CCTCCTCCGAAAAATCTGCTCCTTCTGCTGTGGCAAGAAGCAATAGTTGGATGGGAAGC CAAGGTCCCATGTCAAAAGAAGAGAGAATTATGGAGTTCCTCCACGCAGAGCGACAATA TCTTCTGGAGAAATGATCGCATGTCAACAAGACATGATGGTTATGCAACTAACGATGGA AATCATCCAAGTTGCCAAGAAACGAGGGATTATGCTCCACCATCTAGAGGCTATGCATAC CGTGATAATGGTCATTCTAATCGGGATGAACATTCCTCTAGAGGATATAGAAATCATCGA AGTTCCCGAGAAACTAGGGATTATGCTCCACCATCTAGAGGCCATGCATACCGTGATTAT GGTCATTCTCGTCGGGATGAAAGTTATTCTAGAGGATACAGAAATCGTCGAAGTTCCCGA GAAACTAGGGAGTATGCTCCACCATCTAGAGGCCATGGATACCGTGATTATGGTCATTCT CGTCGACATGAAAGTTATTCTAGAGGATATAGAAATCATCCAAGTTCCTCGAGAAACCAGG GATTATGCTCCACCACATAGAGACTATGCATACCGTGATTATGGTCATTCTAGTTGGGAT GAACATTCCTCTAGAGGATATAGTTATCATGATGGCTACGGTGAGGCCCTTGGTAGAGAT CATTCTGAACATCTAAGTGAAGTTCTTATAGAGATGCACTTCAGAGATACGGGACCTCT CATGGTGCACCACCTGCAAGAGGGCCTCGGATGTCTTATGGTGAAGCACCTGCCACGCA TATAGTAATACACGAGATAGATATGGCAGAAGTTGGGAGAGTTACTCGAGCTGTGGTGAT TTTCATTATTGTGATCGTGAGCATGTTTGCAGAAAAGACCAAAGGAATCCGCCTTCTCTG GGTAGGGTGTCCCTGATCCTCGTGAAGCATGTGGTAGCTCAAGTTATGTGGCATCTATA GTAGATGGTGGGAGAGTCGATCTGAAAAAGGAGACTCGAGCAGATATTA
Restriction Sites:	Please inquire
ACCN:	NM_005058



[View online »](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_005058.2 , NP_005049.1
RefSeq Size:	1900 bp
RefSeq ORF:	1491 bp
Locus ID:	5940
UniProt ID:	P0DJD3
Cytogenetics:	Yq11.223
Gene Summary:	This gene encodes a protein containing an RNA-binding motif in the N-terminus and four SRGY (serine, arginine, glycine, tyrosine) boxes in the C-terminus. This protein is thought to function as a splicing regulator during spermatogenesis. Multiple closely related paralogs of this gene are found in a gene cluster in the AZFb azoospermia factor region of chromosome Y. Most of these related copies are thought to be pseudogenes, though several likely encode functional proteins. [provided by RefSeq, Mar 2016]