

Product datasheet for **SC303581**

TNFSF18 (NM_005092) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	TNFSF18 (NM_005092) Human Untagged Clone
Tag:	Tag Free
Symbol:	TNFSF18
Synonyms:	AITRL; GITRL; hGITRL; TL6; TNLG2A
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>SC303581 representing NM_005092. Blue=Insert sequence Red=Cloning site Green=Tag(s)

```
GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGCCGGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGACATTGCATCCTTCACCCATCACTTGTGAATTTTGTTCACAGCTCTCATTCTCCAAAATG
TGTTTGAGCCACTTGAAAAATAGCCTTAAGCCATTCAAGAAGCTCAAGGAGCTCAGAGATCCTCTGG
AAGCTGTGGCTCTTTGCTCAATAGTTATGTTGCTATTTCTTTGCTCCTTCAGTTGGCTAATCTTTATT
TTTCTCCAATTAGAGACTGCTAAGGAGCCCTGTATGGCTAAGTTTGGACCATTACCCTCAAATGGCAA
ATGGCATCTTCTGAACCTCCTTGCCTGAATAAGGTGTCTGACTGGAAGCTGGAGATACTTCAGAATGGC
TTATATTTAATTTATGGCCAAGTGGCTCCAATGCAAACATAATGATGTAGCTCCTTTTGGAGTGGG
CTGTATAAAAACAAAGACATGATACAAACTCTAACAACAAATCTAAAATCCAAAATGTAGGAGGGACT
TATGAATTGCATGTTGGGACACCATAGACTTGATATTCAACTCTGAGCATCAGTTCTAAAAAATAAT
ACATACTGGGGTATCATTCTACTAGCAAATCCCAATTCATCTCCTAG
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
```

Restriction Sites:	Sgfl-MluI
ACCN:	NM_005092
Insert Size:	600 bp



[View online »](#)

OTI Disclaimer: Due to the inherent nature of this plasmid, standard methods to replicate additional amounts of DNA in E. coli are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at custsupport@origene.com or by calling 301.340.3188 option 3 for pricing and delivery.

The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_005092.3](#)

RefSeq Size: 748 bp

RefSeq ORF: 600 bp

Locus ID: 8995

UniProt ID: [Q9UNG2](#)

Cytogenetics: 1q25.1

Protein Families: Druggable Genome, Transmembrane

Protein Pathways: Cytokine-cytokine receptor interaction

MW: 22.7 kDa

Gene Summary: The protein encoded by this gene is a cytokine that belongs to the tumor necrosis factor (TNF) ligand family. This cytokine is a ligand for receptor TNFRSF18/AITR/GITR. It has been shown to modulate T lymphocyte survival in peripheral tissues. This cytokine is also found to be expressed in endothelial cells, and is thought to be important for interaction between T lymphocytes and endothelial cells. [provided by RefSeq, Jul 2008]