

Product datasheet for **SC303618**

Histone H1t (HIST1H1T) (NM_005323) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Histone H1t (HIST1H1T) (NM_005323) Human Untagged Clone
Tag:	Tag Free
Symbol:	Histone H1t
Synonyms:	dj221C16.2; H1.6; H1FT; H1t; HIST1H1T
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>NCBI ORF sequence for NM_005323, the custom clone sequence may differ by one or more nucleotides

```
ATGTCTGAAACCGTGCCTGCAGCTTCTGCCAGTGTGGTGTAGCCGCTATGGAGAACTCCAACCAAGA  
AGCGAGGGAGGAAGCCGGCTGGCTTGATAAGTGAAGTCGCAAAGTGCCGAACCTCTGTGTCCAAGTT  
GATCACCGAGGCCCTTTCAGTGTACAGGAACGAGTAGGTATGTCTTTGGTTGCGCTCAAGAAGGCATTG  
GCCGCTGCTGGCTACGACGTAGAGAAGAATAACAGCCGCATCAAAGTGTCCCTCAAGAGCTTAGTGAACA  
AGGGAATCCTGGTGCAAACAGGGTACTGGTGTCCGGTTCCTTTAAGCTTAGTAAGAAGGTGATTCC  
TAAATCTACCAGAAGCAAGGCTAAAAAGTCAGTTTCTGCCAAGACCAAGAAGCTGGTTTTATCCAGGGAC  
TCCAAGTCACCAAAGACTGCTAAAACCAATAAGAGAGCCAAGAAGCCGAGAGCGACAACCTCTAAAAGT  
TTAGGAGCGGGAGAAAGGCTAAAGGAGCCAAGGTAAGCAACAGCAGAAGAGCCAGTGAAGGCAAGGGC  
TTCGAAGTCAAAATTGACCCAACATCATGAAGTTAATGTTAGAAAGCCACATCTAAGAAGTAA
```

Restriction Sites:	Please inquire
ACCN:	NM_005323



[View online »](#)

OTI Disclaimer:	<p>Due to the inherent nature of this plasmid, standard methods to replicate additional amounts of DNA in E. coli are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at custsupport@origene.com or by calling 301.340.3188 option 3 for pricing and delivery.</p> <p>The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info</p>
OTI Annotation:	<p>This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.</p>
Components:	<p>The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).</p>
Reconstitution Method:	<ol style="list-style-type: none"> 1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<p>NM_005323.3, NP_005314.2</p>
RefSeq Size:	<p>725 bp</p>
RefSeq ORF:	<p>624 bp</p>
Locus ID:	<p>3010</p>
UniProt ID:	<p>P22492</p>
Cytogenetics:	<p>6p22.2</p>
Gene Summary:	<p>Histones are basic nuclear proteins responsible for nucleosome structure of the chromosomal fiber in eukaryotes. Two molecules of each of the four core histones (H2A, H2B, H3, and H4) form an octamer, around which approximately 146 bp of DNA is wrapped in repeating units, called nucleosomes. The linker histone, H1, interacts with linker DNA between nucleosomes and functions in the compaction of chromatin into higher order structures. This gene is intronless and encodes a replication-dependent histone that is a member of the histone H1 family. Transcripts from this gene lack polyA tails but instead contain a palindromic termination element. This gene is found in the large histone gene cluster on chromosome 6. [provided by RefSeq, Aug 2015]</p>