

Product datasheet for **SC303686**

PDCD7 (NM_005707) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	PDCD7 (NM_005707) Human Untagged Clone
Tag:	Tag Free
Symbol:	PDCD7
Synonyms:	59K; ES18; HES18
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

Fully Sequenced ORF: >SC303686 representing NM_005707.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)


```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAATACGACTCACTATAGGGCGGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGGCCCTGCCACCATTTCTCGGCCAGGGTCGCCACAGCCACC GCCCCCCGAGCCGCCCTCCTGCT
CCTTTTCGGCTGTCCGCCACCGCCGCTGCCCTCCCGGCTTTCCCGCCGCTCTCCCCAGGGCCCGGC
CCTTTTCCGGGGCCTCCGCCCCCTTCTTCAGCCTCCGCTGGCTCTGCAGCCCCGAGCCTCCGCGGAG
GCCTCCCGCGGGGAGGGCGGCTGGCGCCTTCTACCCGGTGCCACCACC GCCCGCTGCCTCCTCCGCCG
CCCCAGTGTGGCCCTTCCCGGGGACCGACGCCGGCGAGCGGCCGCGGCCACC GCCTCCCGGCCGGG
CCGCCCTGGAGCCCGGGTGGCTGAGGCGCCGCCGCCGCCGCCGCGGACGTGCTCGGGGATGCGGCCCTC
CAACGCCTGCGGACCGGCAGTGGCTGGAGGCGGTGTTCCGGACCCCGCGGGGAGGCTGTCCGGTG
CCCCAGCGCACGATGCCGGGCCAGCCTTGCCGAAGTGC GCGCGGATTGCTCCGGGCTGTGCGCCTG
GTGCGGGGCTGCGCGGCTGAGCCAGGCCCTGCGCGAGGCCGAAGCCGACGGCGGCCCTGGGTCTG
CTGTACTCCAGACCGCGCCGCTGCGCGGGAAGTGGCCGAGCGGCTACAGCCGTTGACCCAGGCTGCC
TATGTGGGCGAGGCGCGGAGGAGGCTGGAGAGGGTCCGGCGCCGCCGGCTGCGGGTTCGCGAGAGGGCC
CGGAACGCGAGGCCGAGCGGGAGGCAGAGGCCGCGGGCAGTGAACGCGAGCAGGAGATTGACCGC
TGGAGGGTGAAGTGTGTGCAGGAGTGGAGGAGAAGAAGCGGGAGCAGGAACTCAAAGCAGCCGCTGAT
GGCGTACTATCTGAAGTGAAGAAAAACAAGCAGATACCAAAGAATGGTGGACATTCTACGGGCTTTG
GAGAAATTGAGGAACTGAGGAAAGAGGCTGCAGCGAGGAAAGGGGTCTGTCTCCAGCCTCAGCAGAT
GAGACTTTTACGCATCATCTTCAGCGACTGAGAAACTCATTAAAAAGCGCTCTGAACTGTATGAAGT
GAAGAGAGAGCCCTCAGAGTTATGCTAGAAGGAGAACAAGAGGAAGAGAGGAAAGAGAATTAGAAAAG
AAACAAAGAAAAGAAAAGAGAAAATTTACTTCAGAAACGTGAAATTGAGTCCAAGTTGTTGGGGAT
CCAGATGAGTTCCCACTTGCTCACCTCTTGAGCCTTTCGACAGTATTATCTCCAAGCCGAGCACTCC
CTGCCAGCGCTCATCCAGATCAGGCATGATTGGGATCAGTACCTGGTGCCATCCGATCATCCCAAAGGC
AACTTCGTTCCCAAGGATGGGTCTTCCCGGCTCCCGACACGACATCTGGGCAACTGCTGTTAAG
CTGCATTAG
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
  
```

Restriction Sites:

SgfI-MluI

Plasmid Map:



ACCN:

NM_005707

Insert Size:	1458 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_005707.1</u>
RefSeq Size:	2866 bp
RefSeq ORF:	1458 bp
Locus ID:	10081
UniProt ID:	<u>Q8N8D1</u>
Cytogenetics:	15q22.31
MW:	54.7 kDa
Gene Summary:	This gene encodes a 59 kDa protein that is associated with the U11 small nuclear ribonucleoprotein (snRNP), which is a component of the minor U12-type spliceosome responsible for catalyzing pre-mRNA splicing of U12-type introns. [provided by RefSeq, Dec 2010]