

Product datasheet for **SC303690**

CLEC3A (NM_005752) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	CLEC3A (NM_005752) Human Untagged Clone
Tag:	Tag Free
Symbol:	CLEC3A
Synonyms:	CLECSF1
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>SC303690 representing NM_005752. Blue=Insert sequence Red=Cloning site Green=Tag(s)

GCTCGTTTAGTGAACCGTCAGAATTTGTAAATACGACTCACTATAGGGCGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCC**CGGATCGCC**
ATGGCTCAGCAGGCTTGCCCCAGAGCCATGGCAAAGAATGGACTTGAATTTGCATCCTGGTGATCACC
TTACTCCTGGACCAGACCACCAGCCACACATCCAGATTAAGCCAGGAAGCACAGCAAACGTCGAGTG
AGAGACAAGGATGGAGATCTGAAGACTCAAATTGAAAAGCTCTGGACAGAAGTCAATGCCTTGAAGGAA
ATCAAGCCCTGCAGACAGTCTGTCTCCGAGGCACTAAAGTTCACAAGAAATGCTACCTTGTTCAGAA
GGTTTGAAGCATTCCATGAGGCCAATGAAGACTGCATTTCAAAGGAGGAATCCTGGTTATCCCCAGG
AACTCCGACGAAATCAACGCCCTCAAGACTATGGTAAAAGGAGCCTGCCAGGTGCAATGACTTTTGG
CTGGGCATCAATGACATGGTCACGGAAGCAAGTTTGTGACGTCAACGGAATCGCTATCTCCTTCCTC
AACTGGGACCGTGCACAGCCTAACGGTGGCAAGCGAGAAAACGTGTCTCCTGTTCTCCCAATCAGCTCAG
GGCAAGTGGAGTGATGAGGCCTGTCGCAGCAGCAAGAGATACATATGCGAGTTACCACATCCCTCAA**TAG**

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC

Restriction Sites:	Sgfl-MluI
ACCN:	NM_005752
Insert Size:	621 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).



[View online »](#)

OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_005752.4
RefSeq Size:	1952 bp
RefSeq ORF:	621 bp
Locus ID:	10143
UniProt ID:	O75596
Cytogenetics:	16q23.1
MW:	23.2 kDa
Gene Summary:	Promotes cell adhesion to laminin-332 and fibronectin.[UniProtKB/Swiss-Prot Function] Transcript Variant: This variant (1) represents the longer transcript and encodes the longer isoform (1).