

Product datasheet for **SC303693**

PRMT3 (NM_005788) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	PRMT3 (NM_005788) Human Untagged Clone
Tag:	Tag Free
Symbol:	PRMT3
Synonyms:	HRMT1L3
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

Fully Sequenced ORF: >SC303693 representing NM_005788.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)


```

GCTCGTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGTGCTCGTTAGCGTCAGGCGCTACCGCGCGCCGGGGCGCTGTGGAGAATGAGGAGACCTGCCAGAA
CTGTCCGACAGCGGGGACGAGGCCGCTGGGAGGATGAGGACGATGCAGATCTCCCCACGGCAAGCAG
CAGACCCCTGCCTGTTCTGTAACAGGTTATTCACATCTGCTGAAGAAACATTTTCACACTGTAAGTCT
GAGCATCAGTTTAAATATTGACAGCATGGTTCATAAACATGGACTTGAATTTTATGGATACATTAAGCTA
ATAAATTTTATTAGACTTAAGAATCCTACAGTTGAGTACATGAATTCCATATAACAACCCAGTGCCTTGG
GAGAAAGAAGAGTATTTGAAGCCAGTATTAGAAGATGACCTTTTACTTCAATTTGATGTAGAAGATCTT
TATGAACCGGTGTCAGTACCCTTCTCATACCCCAATGGACTCAGTAAAAACATCTGTTGTTGAAAAA
TTGAAACATATGGAAGCCAGGGCACTGTCTGCTGAAGCCGATTGGCCAGAGCACGTGAGGATCTGCAA
AAAATGAAACAATTTGCTCAGGATTTTGTGATGCACACAGATGCAGAACCTGCTCGTCATCTACTAGT
GTCATTGCGGACCTCCAGGAGGATGAGGATGGTGTGTTATTTTCAGCTCATACGGGCATTATGGGATACAT
GAAGAAATGCTAAAGGACAAAATACGAACAGAAAGCTACCGAGATTTTCATATAACAAAATCCACATATC
TTCAAAGACAAGGTAGTTTTGGATGTTGGGTGTTGGAATCTCTCTATGTTTGTCTAAAGCT
GGGGCGAAGAAGGTTCTTGGAGTTGATCAATCTGAAATACTTTACCAGGCAATGGATATTATAAGACTA
AATAAACTTGAAGATACTATTACACTAATTAAGGAAAGATTGAAGAAGTTCATCTTCTGTAGAAAAA
GTAGATGTTATCATATCTGAGTGGATGGCTATTTTCTTCTGTTTGTGATCTATGTTAGATTCTGCCTT
TATGCAAAGAACAATACTTGGCAAAAGGAGGCTCGGTCTACCCTGACATTTGCACTATCAGCCTTGTA
GCAGTGAGTGATGTAATAACATGCTGATAGAATTGCTTTTTGGGATGATGCTATGGCTTCAAGATG
TCCTGCATGAAGAAAGCAGTTATCCAGAAGCTTTGTGGAAGTTTATAGATCCGAAGACTTTATTTCA
GAACCTTGTGGTATTAAGCATATAGATTGCCATACGACGCTATCTCAGATTTGGAATTTTCATCAGAT
TTTACCCTGAAAAACACAAGGACATCCATGTGCACGGCAATTGCTGGCTACTTTGATATATATTTTGG
AAGAATTGCCACAACAGGGTCGTGTTCTCTACGGGCCCTCAGAGCACAAAAACACTGGAACAAAACA
GTATTTCTACTGGAAAAACCATTTTTCAGTTAAAGCAGGTGAAGCCTTGAAGGAAAGGTCACAGTTCAC
AAGAATAAGAAAGATCCACGTTCTCTACCCTGACCCCTCAGTTGAATAATCAACTCAAACCTATGGT
CTCCAGTGA
ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGCGC
  
```

Restriction Sites:

Sgfl-MluI

Plasmid Map:



ACCN:	NM_005788
Insert Size:	1596 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_005788.3</u>
RefSeq Size:	2743 bp
RefSeq ORF:	1596 bp
Locus ID:	10196
UniProt ID:	<u>O60678</u>
Cytogenetics:	11p15.1
Protein Families:	Druggable Genome
MW:	59.9 kDa
Gene Summary:	<p>This gene belongs to the protein arginine methyltransferase (PRMT) family. The encoded enzyme catalyzes the methylation of guanidino nitrogens of arginyl residues of proteins. The enzyme acts on 40S ribosomal protein S2 (rpS2), which is its major in-vivo substrate, and is involved in the proper maturation of the 80S ribosome. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Aug 2013]</p> <p>Transcript Variant: This variant (1) encodes the longest isoform (1).</p>