

## Product datasheet for SC303715

### SPRR2A (NM\_005988) Human Untagged Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	SPRR2A (NM_005988) Human Untagged Clone
Tag:	Tag Free
Symbol:	SPRR2A
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>SC303715 representing NM_005988. Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGTCTTATCAACAGCAGCAGTGAAGCAGCCCTGCCAGCCACCTCCTGTGTGCCCCACGCCAAAGTGC
CCAGAGCCATGTCCACCCCCGAAGTGCCCTGAGCCCTGCCACCACCAAAGTGCCACAGCCCTGCCCA
CCTCAGCAGTGCCAGCAGAAATATCCTCCTGTGACACCTTCCCCACCCTGCCAGTCAAAGTATCCACCG
AAGAGCAAGTAA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
  
```

Restriction Sites:	SgfI-MluI
ACCN:	NM_005988
Insert Size:	219 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).


[View online »](#)



**Reconstitution Method:**

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

**RefSeq:** [NM\\_005988.2](#)

**RefSeq Size:** 692 bp

**RefSeq ORF:** 219 bp

**Locus ID:** 6700

**UniProt ID:** [P35326](#)

**Cytogenetics:** 1q21.3

**MW:** 8 kDa

**Gene Summary:** Cross-linked envelope protein of keratinocytes. It is a keratinocyte protein that first appears in the cell cytosol, but ultimately becomes cross-linked to membrane proteins by transglutaminase. All that results in the formation of an insoluble envelope beneath the plasma membrane.[UniProtKB/Swiss-Prot Function]