

Product datasheet for **SC303743**

NTN3 (NM_006181) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	NTN3 (NM_006181) Human Untagged Clone
Tag:	Tag Free
Symbol:	NTN3
Synonyms:	NTN2L
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)




[View online »](#)

Fully Sequenced ORF: >SC303743 representing NM_006181.
Blue=Insert sequence Red=Cloning site Green=Tag(s)

```
GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCCGCATCGCC
ATGCCTGGCTGGCCCTGGGGCTGCTGCTGACGGCAGGCAGCTCTTCGCCGCCCTGAGTCCTGGGCCG
CCGGCGCCCGCCGACCCCTGCCACGATGAGGGGGTGCGCCCGCGGCTGCGTGCCAGGACTGGTGAAC
GCCGCCCTGGCCGCGAGGTGCTGGCTTCCAGCACGTGCGGGCGGCCGGCCACTCGGGCCTGCGACGCC
TCCGACCCGCGACGGGCACACTCCCCGCCCTCCTTACTTCCCCAGGGGGCACGGCCAGCCCTCTGTGC
TGGCGCTCGGAGTCCCTGCCTCGGGCGCCCTCAACGTGACTCTCACGGTGCCCTGGGCAAGGCTTTT
GAGCTGGTCTTCGTGAGCCTGCGCTTCTGCTCAGCTCCCCAGCCTCCGTGGCCCTGCTCAAGTCTCAG
GACCATGGCCGACGTGGGCCCGCTGGGCTTCTTCTCTCCACTGTGACTGGACTATGGCCGTCTG
CCTGCCCTGCCAATGGCCAGCTGGCCAGGGCCTGAGGCCCTGTGCTTCCCCGACCCCTGGCCAG
CCTGATGCGAGCGCCTTCTGGCCTTCAAGTGCAGGACAGCAGCCCCCAGGCCCTGGACCTGGACAGC
AGCCAGTGCTCCAAGACTGGGTGACCGCCACCGAGTCCGTGTAGTGCTCACAAGCCCTAGCACGGCA
GGTGACCCAGGGACATGGAGGCCGTGTCCTTACTCCTACGCAGCCACCGACCTCCAGGTGGGGGGG
CGCTGCAAGTGAATGGACATGCCTCACGGTGCCTGCTGGACACACAGGGCCACCTGATCTGCGACTGT
CGGCATGGCACCGAGGGCCCTGACTGCGGCCGTGCAAGCCCTTACTGCGACAGGCCATGGCAGCGG
GCCACTGCCCGGGAATCCACGCCTGCCTCGCTTGTCTGCAACGGCCATGCCCGCCGTGCCGCTTC
AACATGGAGTGTACCGACTGTCCGGCCCGCAGCGGGGGTGTCTGTCTCAACTGCCGGCACAAACC
GCCGGCCGCACTGCCACTACTGCCGGGAGGGCTTCTATCGAGACCCTGGCCGTGCCCTGAGTGACCGT
CGGGCTTGACAGGGCCTGCGACTGTACCCGGTTGGTGTGCTGGCAAGACCTGCAACCAGACCACAGGC
CAGTGTCCCTGCAAGGATGGCGTCACTGGCCTCACCTGCAACCGCTGCGCCCTGGCTTCCAGCAAAGC
CGCTCCCCAGTGGCGCCCTGTGTTAAGACCCCTATCCCTGGACCCACTGAGGACAGCAGCCCTGTGCAG
CCCCAGGACTGTGACTCGCACTGCAAACTGCCCGTGGCAGCTACCGCATCAGCCTAAAGAAGTTCTGC
AAGAAGGACTATGCGGTGCAGGTGGCGGTGGGTGCGCGCGGCGAGGGCGCGCGCGCGTGGACACGCTTC
CCGGTGGCGGTGCTCGCCGTGTTCCGGAGCGGAGAGGAGCGCGCGCGCGGGAGTAGCGCGCTGTGG
GTGCCCGCGGGGATGCGGCTGCGGCTGCCCGCCTGCTCCCCGGCCGCGCTACCTCCTGCTGGGG
GGCGGGCCTGGAGCCGCGGCTGGGGGCGCGGGGGCCGGGGCCCGGGCTCATCGCCGCCCGGGAAGC
CTCGTGTACCTGGAGGACCGGTGGACGCGGCCCTGCGGAGGCTGCAGCGACCGAACGGCGGGGG
CGCTGCAGCGCCGCTGA
```

AGCGGACCGACGGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGAT
ATCCTGGATTACAAGGATGACGACGATAAGGTTTAA

Restriction Sites: SgfI-RsrII

Plasmid Map:


ACCN: NM_006181

Insert Size: 1743 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_006181.2](#)

RefSeq Size: 1986 bp

RefSeq ORF: 1743 bp

Locus ID: 4917

UniProt ID: [O00634](#)

Cytogenetics: 16p13.3

Protein Pathways: Axon guidance

MW: 61.5 kDa

Gene Summary: Netrins control guidance of CNS commissural axons and peripheral motor axons.
[UniProtKB/Swiss-Prot Function]