

Product datasheet for **SC303997**

PCDHB1 (NM_013340) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	PCDHB1 (NM_013340) Human Untagged Clone
Tag:	Tag Free
Symbol:	PCDHB1
Synonyms:	PCDH-BETA1
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

Fully Sequenced ORF: >SC303997 representing NM_013340.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCCGCATCGCC
ATGGCGGGTACGCGCAGAAAAATCTTTGCAAAACAGGCAAGTGGGATCTTCTCATTTTTCTGTGCATA
TCTGTGGGGGATGCGACAACATCCGCTATTCAGTGGCAGAGGAAATGGAGAGCGGCTCGTTTGTGGCC
AACGTAGCTAAGGACCTAGGACTGGAGGTAGGGAAGCTGGCTGCGCGCGGGGCGCGGCTGGTTTCCGAG
GGCAACAAAAATGCATTTCCGGCTCCACCGCAAGACGGGAGATTTGTTTGTGAAGGAGAACTGGATCGG
GAGTCACTTTGTGGCAAAGCCGACCCCTTGTTCTGCACTTTGAAGTAGTCCTGGTGGAGCCGCTGCAG
TCCTTCCGGGCGGAGGTCAAGGTATTTGATATCAATGACAATGCCCAAGTTTTCATAACAAGGAGCCG
CTTTTAAAGATTCGGGAGACACCCCTTTGGGTTACGTTTTCTCTGCAGAGCGCCAGGATCTGGAC
GTGGGCCTTAACGGTCTCCAGAACTACACCCTGAGTGCCAATGGGTATTTCCACCTGCACACCCGCTTC
TGCAGCCACGGGCTAAATATGCTGAGCTGGTGTGAACAAACCCTGGACCGAGAGGAGCAGCCTGAA
GTCAACTTGACAATTACGGCGGTGGACGGCGGGTCCCCGCCTAAGTCTGGCACAGCTCACATCCACGTG
GTGGTTCTGGATGTCAACGACCAGTGCACCGAGTTCTCGCGACTGGTGTACAGAGCCAGGTATCAGAG
AACAGCCCAATGGCTCTTTGGTGGCCACGGTGACTGCCGTGGACCTAGACGAGGGCACCAACAAAGCG
ATAACTTACTCTTTAGCTCAAAACCCAGAAGCAATTCTCAAGACGTTTCAGATTGACCCTCAAAATGGA
GAAGTTCGACTAAGAGGACCCCTCGATTTTGAAGCCATTGAAACATACGACATTGACATTCAGCTACA
GATGGTGGAGGCCTCTCTGCCACAGCAAAGTCTGGTAGAAGTGGTGGATGTGAATGACAATCCTCCC
GAAGTGTGGTCTCCTCTGTGTCCAGCCACTCCCTGAAGACTCACCCACAGACAGTAGTAGCCCTT
TTCATATCAGAGACCGGACATTCGAGTGGGAGGAAAAGTCACTGCTTCTCAGAGAAGACCTTCCC
TTTGTAAATCAAACCTACATTTGGGAATTTACTACTGGTCACTGACAGAAGCTTGGATCGGGAGGAG
GTCTCAGGCTATAATATACCATTGTTGCCATGGATACTGGACCACCTAGCTTGTCTGCCGAGACTATG
ATAGAGGTGCTAATATCCGACGTTAATGACAATCCTCCAATATTTCCGGGAAGATTCTATATCTTGACT
GTTCCGAGAAAACAACAGTCTGCGGTTTTTATTGGCAAAGTCCATGCTGAGGATCTTGATTTGGGTGAG
AATGCCCAAATAACATATCTCTGTTGCCTCAAAAAACGGAGATCTTTCAGTCTTTGCTTACATATCC
ATAAATTCAGGCAATGGGAAGCTCTACGCGCTGAGAACCATGGATTATGAGGCCATTCAAGATTTCAA
TTTGTGGTAAAGGCAACTGATGGGGCTTCTGTCACTGAGTAGCCAAGTACTGTCAGAGTGGTTGTC
CTAGATGACAATGACAATCGTCCAATGATCTTATACCCACTGCAGAACGGCACCTTGCCCTGCAATGAC
CTGGTGCCAGGTCTGCAGAGGCAGGCTACCTAGTGACCAAAGTGGTGGCTGTGGATGGTGACTCAGGT
CAGAAATCTTGGCTTTCATATCATCTACTTAAGGCCACTGACCTTGGGTTATTTCTGTTCAAAGACAA
AATGGAGAAATCCATACATTAAGGCAGATATCTGAGAGAGACCCATGATGCAGAAATTGATCATTCTT
GTTCAGGATCACGGCCAACCAGCTTTTCCACTACTGTCTCACTCAACATCCTGTGGTAGATGGCTTT
TCAGAGCCCTACCTGCAGTTCAGGATCCAACCAAGCATTCTAGAAAGGTAATCCATCCACTAAATAT
TTGGTCATTTCTGGTCATCCTTTCTTTCTCTTCTCTCTCTGTGATAGTATCTTATTATACAT
GTCTACCAAAAAGATTAATATAGAGAAAAGTTCACAATTCAAGAGCATTCTATGATGACTGTAATTTCC
TCTAACAACTGGTACAAGGACAAGGCAATGGATCCTTATCTCGGCTTGTCCATATGAAATGTGTTCA
GCCACTGGCACTGGTAAATAGTGAGTTTCGCTTTCTTAAGCGTTTTATGCCCAACTCCCTTTCCCTCAT
GCCACTGGGGAGATAAAAATGGAGGCTGGCTCCAGTTTGCCCAAAATTCGATAGGAATAAGTCTCAG
AGATTAGAGGGCCATGACCAGGTATCTGATGACTATATGTAG
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
  
```

Restriction Sites: SgfI-MluI
ACCN: NM_013340
Insert Size: 2457 bp

OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_013340.3
RefSeq Size:	2587 bp
RefSeq ORF:	2457 bp
Locus ID:	29930
UniProt ID:	Q9Y5F3
Cytogenetics:	5q31.3
Protein Families:	Druggable Genome, Transmembrane
MW:	90.5 kDa
Gene Summary:	<p>This gene is a member of the protocadherin beta gene cluster, one of three related gene clusters tandemly linked on chromosome five. The gene clusters demonstrate an unusual genomic organization similar to that of B-cell and T-cell receptor gene clusters. The beta cluster contains 16 genes and 3 pseudogenes, each encoding 6 extracellular cadherin domains and a cytoplasmic tail that deviates from others in the cadherin superfamily. The extracellular domains interact in a homophilic manner to specify differential cell-cell connections. Unlike the alpha and gamma clusters, the transcripts from these genes are made up of only one large exon, not sharing common 3' exons as expected. These neural cadherin-like cell adhesion proteins are integral plasma membrane proteins. Their specific functions are unknown but they most likely play a critical role in the establishment and function of specific cell-cell neural connections. [provided by RefSeq, Jul 2008]</p>