

Product datasheet for **SC304075**

ADAM29 (NM_014269) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	ADAM29 (NM_014269) Human Untagged Clone
Tag:	Tag Free
Symbol:	ADAM29
Synonyms:	CT73; svph1
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)




[View online »](#)

Fully Sequenced ORF: >SC304075 representing NM_014269.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCCGCATCGCC
ATGAAGATGTTACTCTGCTGCATTGCCTTGGGGTGTTCCTGTCCTGTTCTGGACACATCCAGGATGAG
CACCCCAATATCACAGCCCTCCGGATGTGGTGATTCTGTGAGGATAACTGGCACCACCAGAGGCATG
ACACCTCCAGGCTGGCTCCTATATCCTGCCCTTGGAGGCCAGAAACACATTATCCACATAAAGGTC
AAGAAGCTTTTGTTCCAAACCTCCCTGTGTTACCTACACAGACCAGGGTGTATCCTTGAGGAC
CAGCCATTTGCCAGAATAACTGCTACTATCATGGTTATGTGGAAGGGGACCAGAATCCCTGGTTTCC
CTCAGTACCTGTTTTGGGGTTTTCAAGGAATATTACAGATAAATGACTTTGCTTATGAAATCAAGCCC
CTAGCATTTTCTACCACGTTTGAACATCTGGTATAACAAGATGGACAGTGAGGAGAAAACATTTTCAACC
ATGAGATCCGGATTTATGCAAAATGAAATAACATGCCAATGGAATTTGAAGAAATTGATAATCCACT
CAGAAGCAAAGTCTTATGTGGGCTGGTGGATCCATTTTAGGATTGTTGAAATTGATGCTGCATTGAT
AATTATCTGTACATTCGTTATGAAAGGAACGACTCAAAGTTGCTGGAGGATCTATATGTTATTGTTAAT
ATAGTGGATTCCATTTTGGATGTCATTGGTGTTAAGGTGTTATTATTTGGTTTGGAGATCTGGACCAAT
AAAAACCTCATTGTAGTAGATGATGTAAGGAAATCTGTGCACCTGTATTGCAAGTGAAGTCCGGAGAAC
ATTACGCCCGGATGCAACATGACACCTCACATCTTTTCACAACTCTAGGATTAAGAGGGTTAAGTGGC
ATAGGAGCTTTTAGAGGAATGTGTACACCACACCGTAGTTGTGCAATTGTTACTTTCATGAACAAAAC
TTGGGCACTTTTTCAATTGCAGTGGCTCATCATCTAGGTCATAATTTGGGCATGAACCATGATGAGGAT
ACATGTCGTTGTTCAACCTAGATGCATAATGCATGAAGGCAACCCACCAATAACTAAATTTAGCAAT
TGTAGTTATGGTGATTTTTGGGAATATACTGTAGAGAGGACAAAGTGTGTTGCTTGAACAGTACACACA
AAGGACATCTTTAATGTGAAGCGCTGTGGGAATGGTGTGTTGGAAGAAGGAGAAGAGTGTGACTGTGGA
CCTTTAAAGCATTGTGCAAAAGATCCCTGCTGTCTGCAAAATGCACTCTGACTGATGGTTCTACTTGT
GCTTTTGGGCTTTGTTGCAAAGACTGCAAGTCCCTACCATCAGGGAAAAGTGTGTAGAAAGGAGGTCAAT
GAATGTGATCTTCCAGAGTGGTGAATGGTACTTCCCATAAGTGCCAGATGACTTTTATGTGGAAGAT
GGAATTCCTGTAAGGAGAGGGGCTACTGCTATGAAAAGAGCTGTGATGACCGCAATGAACAGTGTAGG
AGGATTTTTGGTGCAGGCGCAAATACTGCAAGTGAGACTTGCTACAAAGAATTGAACACCTTAGGTGAC
CGTGTGGTCACTGTGGTATCAAAAATGCTACATATAAAGTGAATATCTCAGATGTCCAGTGTGGA
AGAATTCAGTGTGAGAATGTGACAGAAATCCCAATATGAGTGATCATACTACTGTGATTGGGCTCGC
TTCAATGACATAATGTGCTGGAGTACTGATTACCATTTGGGGATGAAGGGACCTGATTTGGTGAAGTG
AAAGATGGAACAGAGTGTGGGATAGATCATATATGCATCCACAGGCACTGTGTCATATAACCATCTTG
AATAGTAATTGCTCACCTGCATTTTGAACAAGAGGGGCATCTGCAACAATAAACATCACTGCCATTGC
AATTATCTGTGGGACCCTCCCAACTGCCTGATAAAAGGCTATGGAGGTAGTGTGACAGTGGCCACCC
CCTAAGAGAAAGAAGAAAAAGAAGTTCTGTTATCTGTGTATATTGTTGCTTATTGTTTTGTTATTTTA
TTATGTTGCTTTTATCGACTTTGTAAAAAAGTAAACCAATAAAAAAGCAGCAAGATGTTCAAATCCA
TCTGCAAAAGAAGAGGAAAAAATTCAGCGTCGACCTCATGAGTTACCTCCCAGAGTCAACCTTGGGTG
ATGCCTCCCAGAGTCAACCTCCTGTGACGCCTCCCAGAGTCATCCTCAGGTGATGCCTCCCAGAGT
CAACCTCCTGTGACACCCTCCCAGAGTCAACCTCGGGTGTGCTTCTCAGAGTCAACCTCCTGTGATG
CCTTCCCAGAGTCATCCTCAGTTGACGCCTTCCCAGAGTCAACCTCCTGTGACACCCTCCCAGAGGCAA
CCTCAGTTGATGCCTTCCCAGAGTCAACCTCCTGTGACGCCTCCTAG
ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
  
```

Restriction Sites: Sgfl-Mlul

Plasmid Map:


ACCN: NM_014269

Insert Size: 2463 bp

OTI Disclaimer: Due to the inherent nature of this plasmid, standard methods to replicate additional amounts of DNA in E. coli are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at custsupport@origene.com or by calling 301.340.3188 option 3 for pricing and delivery.

The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

- Reconstitution Method:**
1. Centrifuge at 5,000xg for 5min.
 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
 3. Close the tube and incubate for 10 minutes at room temperature.
 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq:	NM_014269.4
RefSeq Size:	3325 bp
RefSeq ORF:	2463 bp
Locus ID:	11086
UniProt ID:	Q9UKF5
Cytogenetics:	4q34.1
Protein Families:	Druggable Genome, Transmembrane
MW:	92.8 kDa
Gene Summary:	<p>This gene encodes a member of the ADAM (a disintegrin and metalloprotease domain) family. Members of this family are membrane-anchored proteins structurally related to snake venom disintegrins, and have been implicated in a variety of biological processes involving cell-cell and cell-matrix interactions, including fertilization, muscle development, and neurogenesis. The protein encoded by this gene is highly expressed in testis and may be involved in human spermatogenesis. Alternative splicing results in multiple transcript variants that encode the same protein. [provided by RefSeq, Jul 2008]</p> <p>Transcript Variant: This variant (1) encodes the functional protein. Variants 1, 2, 3, 4, 5, 6, and 7 encode the same protein.</p>