

## Product datasheet for **SC304106**

### IL36 beta (IL36B) (NM\_014438) Human Untagged Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	IL36 beta (IL36B) (NM_014438) Human Untagged Clone
Tag:	Tag Free
Symbol:	IL36 beta
Synonyms:	FIL1; FIL1-(ETA); FIL1H; FILI-(ETA); IL-1F8; IL-1H2; IL1-ETA; IL1F8; IL1H2
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>NCBI ORF sequence for NM_014438, the custom clone sequence may differ by one or more nucleotides

```
ATGAACCCACAACGGGAGGCAGCACCCAAATCCTATGCTATTCGTGATTCTCGACAGATGGTGTGGGTCC
TGAGTGGAAATTCCTTTAATAGCAGCTCCTCTTAGCCGAGCATTAAAGCCTGTCACTCTTCATTTAATAGC
CTGTAGAGACACAGAATTCAGTGACAAGGAAAAGGTAATATGGTTTACCTGGGAATCAAGGGAAAAGAT
CTCTGTCTCTTCTGTGCAGAAATTCAGGGCAAGCCTACTTTGCAGCTTAAGCTTCAGGGCTCCCAAGATA
ACATAGGGAAGGACACTTGCTGGAACTAGTTGGAATTCACACATGCATAAACCTGGATGTGAGAGAGAG
CTGCTTCATGGGAACCTTGACCAATGGGAATAGGAGTGGGTAGAAAGAAGTGGGAAGAGTTCCTTTCAA
CATCACCATCTCAGGAAGAAGGACAAAGATTTCTCATCCATGCGGACCAACATAGGAATGCCAGGAAGGA
TGTAG
```

Restriction Sites:	Please inquire
ACCN:	NM_014438
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.



[View online »](#)

<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<u>NM_014438.3, NP_055253.2</u>
<b>RefSeq Size:</b>	1182 bp
<b>RefSeq ORF:</b>	495 bp
<b>Locus ID:</b>	27177
<b>UniProt ID:</b>	<u>Q9NZH7</u>
<b>Cytogenetics:</b>	2q14.1
<b>Protein Families:</b>	Druggable Genome, Secreted Protein
<b>Gene Summary:</b>	<p>The protein encoded by this gene is a member of the interleukin 1 cytokine family. Protein structure modeling indicated that this cytokine may contain a 12-stranded beta-trefoil structure that is conserved between IL1A (IL-A alpha) and IL1B (IL-1 beta). This gene and eight other interleukin 1 family genes form a cytokine gene cluster on chromosome 2. Two alternatively spliced transcript variants encoding distinct isoforms have been reported. [provided by RefSeq, Jul 2008]</p> <p>Transcript Variant: This variant (1) encodes a longer isoform (1). Sequence Note: This RefSeq record was created from transcript and genomic sequence data to make the sequence consistent with the reference genome assembly. The genomic coordinates used for the transcript record were based on transcript alignments.</p>