

Product datasheet for **SC304124**

CRTR1 (TFCP2L1) (NM_014553) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	CRTR1 (TFCP2L1) (NM_014553) Human Untagged Clone
Tag:	Tag Free
Symbol:	TFCP2L1
Synonyms:	CRTR1; LBP-9; LBP9
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

Fully Sequenced ORF: >SC304124 representing NM_014553.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)


```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAATACGACTCACTATAGGGCGGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGCTCTTCTGGCACACGCAGCCCGAGCACTACAACCAGCACAACTCCGGCAGCTACCTGCGTGATGTG
CTCGCTCTGCCCATCTTCAAGCAGGAGGAACCCAGCTGTCCCCGAGAACGAGGCCCGCTGCCACCC
CTGCAATATGTGTTGTGTGCTGCCACGTCCCAGCCGTGAAGCTGCATGAAGAGACGCTGACCTACCTC
AACCAAGGTCAGTCTTATGAAATCCGACTACTGGAGAATCGGAAGCTGGGAGACTTCAAGATCTGAAC
ACAAAATATGTCAAGAGCATCATCCGTGTGGTCTTCCATGACCGCCGGCTGCAGTATACGGAGCACCAG
CAGCTGGAGGGCTGGCGTGGAGTCGGCCAGGGGACCGGATCCTGGACATCGATATTCCACTGTCTGTT
GGTATCTTGACCCAGGGCCAGCCGACCCAGCTGAATGCAGTCGAGTTTTTGTGGACCCTGCGAAG
AGAGCTTCTGCATTCATTAGGTACTGCATCAGCACAGAATTCACCCCGGAAGCACGGGGCGAG
AAGGGAGTGCCCTTCGAGTCCAGATTGACACGTTTAAGCAGAACGAGAATGGGGAGTACACGGAGCAC
CTGCACTCAGCCAGCTGCCAGATCAAGGTGTTCAAGCCGAAGGGAGCCGATCGGAAACAGAAGACTGAC
CGGAGAAGATGGAGAAAAGAACTGCCAAGAGAAGGAGAAAATACCAGCGCTCCTATGAAACCACCATC
CTCACAGAGTGCTCTCCATGGCCGACGTGGCCTACCAGGTGAACAGCGCCCCGTCGCCAAGCTACAAT
GGTTCTCAAACAGCTTTGGCCTCGGCGAAGGCAACGCCTCTCCGACCCACCCGGTGGAGGCCCTGCC
GTGGGCAGTGACCACCTGCTCCCATCAGCTTCGATCCAGGATGCCCAGCAGTGGCTTACCAGCAACAGG
TTCTCGCAGTTCTGCCGGCTCTTGGCAGCTTCTCAGGTGCTGACTTGCTGAAGATGTCCCGAGATGAT
TTGGTCCAGATCTGTGGTCCCGCAGATGGGATCCGGCTTTCAACGCCATCAAAGCCGGAATGTGAGG
CCAAAGATGACCATTTATGTCTGTCCAGGAGCTGGAGCAGAATCGAGTGCCCTGCAGCAGAAGCGGGAC
GGCAGTGGAGACAGCAACCTGTCTGTACCACGCCATCTTCTGGAAGAGCTGACCACCTGGAGCTG
ATTGAGAAGATCGCCAACCTGTACAGCATCTCCCCCAGCACATCCACCGAGTCTACCGGCAGGGCCCC
ACGGGCATCCATGTGGTGGTGGAGCAACGAGATGGTGCAGAACTTCCAAGATGAATCCTGTTTTGTCTC
AGCACAATTAAGCTGAGAGCAATGATGGCTACCACATCATCCTGAAATGTGGACTTGA
ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
  
```

Restriction Sites:

SgfI-MluI

Plasmid Map:



ACCN: NM_014553

Insert Size: 1440 bp

OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none"> 1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_014553.2</u>
RefSeq Size:	9306 bp
RefSeq ORF:	1440 bp
Locus ID:	29842
UniProt ID:	<u>Q9NZI6</u>
Cytogenetics:	2q14.2
Protein Families:	Transcription Factors
MW:	54.6 kDa
Gene Summary:	Transcription factor that facilitates establishment and maintenance of pluripotency in embryonic stem cells (ESCs) (PubMed:25215486, PubMed:26906118). With KLF2, acts as the major effector of self-renewal that mediates induction of pluripotency downstream of LIF/STAT3 and Wnt/beta-catenin signaling (By similarity). Required for normal duct development in the salivary gland and kidney (By similarity). Coordinates the development of the kidney collecting ducts intercalated (IC) and principal (PC) cells, which regulate acid-base and salt-water homeostasis, respectively (By similarity). Regulates the expression of IC genes including subunits B1 and D2 of the V-ATPase complex, OXGR1, CA12, SLC4A1, AQP6 and IC-specific transcription factor FOXI1 (By similarity). Regulates also the expression of JAG1 and subsequent notch signaling in the collecting duct (By similarity). JAG1 initiates notch signaling in PCs but inhibits notch signaling in ICs (By similarity). Acts as a transcriptional suppressor that may suppress UBP1-mediated transcriptional activation (By similarity). Modulates the placental expression of CYP11A1 (PubMed:10644752).[UniProtKB/Swiss-Prot Function]