

Product datasheet for **SC304172**

ZNF623 (NM_014789) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	ZNF623 (NM_014789) Human Untagged Clone
Tag:	Tag Free
Symbol:	ZNF623
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

Fully Sequenced ORF: >SC304172 representing NM_014789.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCACGATCGCC
ATGATTTTGTGTCTTTTGTTCAGATTCTAATGTAGGAAGTGGTGAGAAGAAGGTGACTGAAGCCTGG
ATTTCTGAGGATGAAACTCACATAGGACGACGTACAGACAGACTCACGGTGATGGAGCTCCCTCTCCC
GAGTCTGAGGAAGTCCACGAGCCAGATTAGGGGAGCTCTTGGGAAATCCAGAAGGTGAGAGCCTGGGG
AGTTCCCTCTCAGGACAGGGGCTGCAAGCAGGTGACAGTGACCCATTGGAAGATCCAGACAGGAGAG
ACAGCTCAAGTGTGACCAAGTCAGGAAGAAACCATATTCTGAAGTCAAGCCTTCTTCTGCTTCAGAGA
GAGCTCATAGAGGGGAAGCCAATCCTTGCGATATCTGTGGCAAAACCTTCACGTTTAATTCGGACCTA
GTTAGGCATCGGATTTTCGATGCTGGGAGAAACCTTACACGTGCGATCAGTGTGGGAAAGGCTTTGGC
CAGAGCTCACACCTTATGGAGCATCAGAGAATTCACACTGGAGAGAGACTCTACGTCTGTAATGTGTGT
GGGAAAGACTTCATTCACTATTCAAGTCTCATTGAGCATCAGCGCTTCATTGAGGAGAAAGCCCTTC
AAATGTGCGCAGTGTGGGAAGGCGTTTTGTACAGTTCAGACCTGATTAGGCACCAGAGAGTTACACC
AGAGAGAGACCTTTGAATGCAAAGAGTGTGGGAAAGGCTTCAGTCAGAGCTCCTTACTTATTGCCAT
CAGAGGATTACACGGGAGAAAGGCCCTATGAGTGCAATGAATGTGGGAAATCCTTCATAAGGAGCTCG
AGCCTCATTGCCATTATCAGATCCACACAGAAGTGAACAGTATGAATGCAAAGATGTGGGAAGGCA
TTCCGTCATCGCTCAGACCTTATTGAACACCAGAGAATTCACACCGGAGAGAGACCTTTGAATGCAAT
GAGTGTGGGAAAGCCTTTATTGAGGATTCAGAGCTCATTGAGTGTGGGAAAGCCTTTCTTCTGAGTTCATACCTTATT
CACACTGGAGAGAACTCTATGAATGTAACGAGTGTGGGAAAGCCTTTCTTCTGAGTTCATACCTTATT
CGACACCAGAAATCCACACTGGAGAGAGAGTGTATGAATGTAAGGAATGTGGGAAAGCCTTTCTCCAG
AAAGCCCATCTCACTGAGCACCAGAAGATCCACTCTGGGACAGGCCCTTCGAATGTAAGACTGTGGG
AAAGCCTTCATCCAGAGTCCAAGCTGCTTCTGCACCAGATTATTCACACTGGAGAAAAGCCTATGTG
TGCAGTTATTGTGGGAAAGGCTTTATTGAGAGTCAAACCTTCTTCAACACCAGAAAATTCATACTGAA
GAGAAGCTCTATGAATGTAGTCAGTATGGGAGAGATTTAACTCAACTACAAACGTTAAAAATAATCAA
AGGTTTACCAAGAGGAGTCTCCTTGAGTAAGGCCCCATACATTTGGGTGAGAGGTCTGTAGATAAG
GGGGAACACACAGGTAACCTTAATAA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGCGC
  
```

Restriction Sites: SgfI-MluI

ACCN: NM_014789

Insert Size: 1611 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_014789.3</u>
RefSeq Size:	3947 bp
RefSeq ORF:	1611 bp
Locus ID:	9831
UniProt ID:	<u>O75123</u>
Cytogenetics:	8q24.3
Domains:	zf-C2H2
MW:	61.4 kDa
Gene Summary:	May be involved in transcriptional regulation.[UniProtKB/Swiss-Prot Function] Transcript Variant: This variant (1) encodes the longer isoform (1). COMPLETENESS: complete on the 3' end.