

Product datasheet for **SC304178**

IQSEC1 (NM_014869) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	IQSEC1 (NM_014869) Human Untagged Clone
Tag:	Tag Free
Symbol:	IQSEC1
Synonyms:	ARF-GEP100; ARFGEP100; BRAG2; GEP100; IDDSSBA
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)




[View online »](#)

Fully Sequenced ORF: >SC304178 representing NM_014869.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCCGCATCGCC
ATGTGGTGCCTGCACTGCAACTCGGAGAGGACCCAGTCCCTTCTGGAGCTAGAGCTTGACAGCGCGT
GAGGGCAGGCCCCCAGCAGTGAAGTGGACATCCCTGGACAGCCCTCAGCCTACCCCAAGGCCCC
TTGGTGCCCGTTCCAGCCTGAGCCGGATCACTACGAGCACACGTCAGTGGGAGCCTATGGGCTGTAC
TCGGGGCCCGCGGGCAACAGCAGCGCACGCGGAGGCCAAGCTGCAGCACTCGACCTCCATCCTGCGC
AAGCAGGCTGAGGAGGAGCCATCAAGCGCTCACGCTCACTCTCCGAGAGCTATGAGCTCTCCTCGGAC
CTGCAGGACAAGCAGGTGGAGATGCTAGAACGAAAGTATGGGGGGCGCTGGTAACCCGCCATGCGGCC
CGCACCATCCAGACGGGTTTCGCCAGTACCAGATGAACAAGAATTGAGCGCTTGCGCAGCTCCATG
TCAGAGAACCGCATGTCACGCCGATTGTGCTGTCCAACATGAGGATGCAGTTCTCTTTGAGGGGCT
GAGAAAGTGCACAGCTCCTACTTCGAGGGGAAGCAGGTCTCAGTACTAACGACGGCTCCAGCTGGGA
GCCCTGGTGTCCCCTGAGTGTGGTGACCTCAGCGAGCCACCACCCTCAAGTCTCCGGCCCCCTCAGT
GACTTTGCGGACGCCATACCGAGCTGGAGGACGCTTCTCTAGCAAGTGAAATCACTGGCCGAGTCC
ATCGACGATGCCCTCAACTGCCGACGCTGCACACTGAGGAGGCACCGGCCCTGGATGCGGCGGGGCC
CGGGACACCGAACCCAGACAGCCCTGCACGGCATGGACCACCGCAAAGTGGACGAGATGACGGCCTCG
TACAGTGTATGTCACCCTGTACATCGATGAGGAGGAGTGTGCGCCCTCTGCCCTCTCGCAGGCAGGG
GACCGGCCGTCCAGCACCGAGTCCGACCTGCGGCTACGGGCTGGGGGCGCAGCCCCAGACTACTGGCC
CTGGCCCAAAAGAGGACAAGGCTGACACGGACACGAGCTGCCGGAGCACGCCGTGCTGGAGCGGCAG
GAGCAGCGGCTGCGGGTGGAGCATCTGCCGCTGCTCACCATCGAGCCACCAGCGACAGCTCTGTGGAC
CTTAGTACCCTCGGAGCGGGGTCACTCAAGAGCGAGAGTGTACGAGCGCAGCTTGGCGGGCAG
CAGGGCAGTCCCAAGCATGGTCCCACAGCGGCCCCCAAGAGCCTCCCGGGAGGAGCCTGAGTTG
CGGCCCGGCCCCAGGCCCTGGACAGCCACTTGCCATCAATGGCTCAGCAACCGGACAGGAAG
TCTGAGTCGGAATACTCAGACGGTGAATGACAGCATCAACAGCACGTCCTCAACTCCAACGATACCATC
AACTGCAGCTCCGAGTCACTCGTCCCGTGCAGCCCTGCGGGAGCAGACGCTCAGCAAGCAGACCTACCAC
AAGGAGGCCCGCAACAGCTGGGACTCGCTGCCTTTAGCAACGATGTATCCGCAAGAGGCACTACCGC
ATCGCCCTGAACCTTTCAACAAGAAGCCTGAGAAGGGAGTCCAGTACCTCATCGAGCGTGGCTTTGTG
CCCAGACGCCCGTGGGGTGGCCACTTCTGCTGCAGCGCAAGGGCCTCAGCCGGCAGATGATCGGC
GAGTTCCTGGGAACCGGCAGAAGCAGTCAACCGTGACGTGCTCGACTGCGTGTGGACGAGATGGAC
TTCTCTACCATGGAGCTGGATGAGGCCCTCAGGAAATCCAGGCCACATCCGTGTCCAAGGGGAGGCT
CAGAAAGTGGAGCGGCTCATAGAGGCGTTCAGCCAGCGCTACTGCATCTGCAACCCTGGGGTGGTGCGG
CAATTCGGAAACCCAGACACCATTTTCATCCTGGCCTTCGCCATCATCCTGCTGAACACCGACATGTAC
AGCCCCAATGTCAAGCCCGAGCGGAAAATGAAGCTAGAGGACTTCATCAAGAACCTCCGAGGTGTGGAC
GATGGTGAAGACATTCGCCGTGAGATGCTGATGGGGATCTATGAACGGATCCGTAAGCGAGAGCTAAG
ACCAATGAGGACCATGTGTCCAGGTGCAGAAGGTGGAGAAGCTCATTGTGGGAAAAAGCCGATCGGA
TCCCTGCATCCCGGGCTCGGCTGTGTGCTCTCTGCCCCACCGTCCGTTGGTCTGCTACTGCCGGCTC
TTGAGGTTCCAGACCCAAACAAGCCCCAGAACTCGACTACACCAGCGAGAAATCTTCTGTTCAAC
GACCTCCTGGTGGTACCAAGATCTTCCAGAAGAAGAAGAACTCGGTGAGGTACAGCTTCCGACAGTCC
TTCTCCTGTACGGCATGCAGGTCTGCTCTTCGAGAACCAGTACTACCCCAATGGCATCCGGCTCACC
TCGTCTGTCCCGGAGCAGATATCAAAGTGTTAATAAACTTCAACGCCCCCAACCCTCAAGACCGGAAG
AAATTCACCGATGACCTGCGGGAGTCCATTGCGGAAGTCCAAGAGATGGAGAAGCACAGGATAGAGTCG
GAGCTCGAGAAGCAGAAAGGCGTGTGCGGCCAGCATGTCCAGTGTCTAGCCTCAAAAAGGAGTCG
GGCAACGGAACACTGAGCCGGGCTGCCTGGACGACAGCTATGCCAGCGGTGAGGGCCTCAAGCGCAGC
GCCCTCAGCAGCTCCCTGCGGGACCTCTCGGAAGCCGGGAAGCGAGGGCGTGCAGCAGTCCGGGATCG
CTAGAGAGCAATGTGAATTTACGCTTTTCGAGCCACTGCAGCCGTGAGTGTGCTCTCTAA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
  
```

Restriction Sites: SgfI-MluI

Plasmid Map:


ACCN: NM_014869

Insert Size: 2892 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_014869.6](#)

RefSeq Size: 5486 bp

RefSeq ORF: 2892 bp

Locus ID: 9922

UniProt ID: [Q6DN90](#)

Cytogenetics: 3p25.2-p25.1

Protein Pathways:	Endocytosis
MW:	108.3 kDa
Gene Summary:	<p>Guanine nucleotide exchange factor for ARF1 and ARF6 (PubMed:24058294). Guanine nucleotide exchange factor activity is enhanced by lipid binding (PubMed:24058294). Accelerates GTP binding by ARFs of all three classes. Guanine nucleotide exchange protein for ARF6, mediating internalisation of beta-1 integrin.[UniProtKB/Swiss-Prot Function]</p> <p>Transcript Variant: This variant (2) differs in the 5' UTR, lacks a portion of the 5' coding region and initiates translation at an alternate start codon, compared to variant 1. The encoded isoform (b) is shorter and has a distinct N-terminus, compared to isoform a. Sequence Note: This RefSeq record was created from transcript and genomic sequence data to make the sequence consistent with the reference genome assembly. The genomic coordinates used for the transcript record were based on transcript alignments.</p>