

Product datasheet for **SC304194**

SV2C (NM_014979) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	SV2C (NM_014979) Human Untagged Clone
Tag:	Tag Free
Symbol:	SV2C
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

Fully Sequenced ORF: >SC304194 representing NM_014979.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCCGGATCGCC
ATGGAAGACTCTTACAAGGATAGGACTTCACTGATGAAGGGTGCCAAGGACATTGCCAGAGAGGTGAAG
AAACAAACAGTAAAGAAGGTGAATCAAGCTGTGGACCGAGCCAGGATGAATACCCAGAGGTCCTAC
AGTCGGTTCCAAGATGAAGAAGATGATGATGACTACTACCCGGCTGGAGAAACCTATAATGGTGAGGCC
AACGATGACGAAGGCTCAAGTGAAGCCACTGAGGGGCATGATGAAGATGATGAGATCTATGAGGGGAG
TATCAGGGCATCCCCAGTATGAACCAAGCGAAGGACAGCATCGTGTCAAGTGGGCAGCCCAAGGGCGAT
GAGTACAAGGACCGGGGAGCTGGAATCAGAAAGGAGAGCTGACGAGGAAGAGTTAGCCAGCAGTAT
GAGCTGATAATCCAAGAATGCGGTATGGTCGTTTTTCAGTGGGCCCTTTCTTCGTCTGGGCATGGCT
CTTATGGCAGACGGTGTAGAGGTGTTTTCGTTGGCTTCGTGTTACCCAGTGTGAGACAGACCTCTGC
ATCCCAAATTCAGGATCTGGATGGCTAGGCAGCATAGTGTACCTCGGGATGATGGTGGGGCGTTCTTC
TGGGGAGGACTGGCAGACAAAGTGGGAAGGAAACAGTCTCTTCTGATTTGCATGTCTGTCAACGGATTC
TTTGCCTTCTTTCTTCAATTTGTCCAAGGTTATGGCTTCTTTCTTCTGTGCTTACTTTCTGGATTC
GGGATTGGAGGAGCCATACCCACTGTGTTCTCGTACTTTGCTGAAGTCTGGCCCGGAAAAGCGGGGC
GAACACTTGAGCTGGCTCTGCATGTTCTGGATGATCGGTGGCATCTACGCCTCTGCCATGGCCTGGGCC
ATCATCCCGCACTACGGGTGGAGCTTCAAGCATGGGATCGGCCTACCAGTTTACAGTTGGCGTGTGTTT
GTCATCGTCTGTGCACTCCCCTGTGTCTCCTCCGTGGTGGCCCTCACATTCATGCCTGAAAGCCACGA
TTCTTGTGGAGGTTGAAAACATGATGAAGCTGGATGATTCTGAAGTTAATTCATGACACCAACATG
AGAGCCCGGGGTCAAGCTGAGAAGGCTTCAAGGTAACAAAAAATAAACTCCTAAACAAATAGATGAG
CTGATTGAAATTGAGAGTGACACAGGAACATGGTATAGGAGGTGTTTTGTTCCGATCCGCACCGAGCTG
TACGGAATTTGGTTGACTTTTATGAGATGTTTCAACTACCCAGTCAGGGATAATAACAATAAAGCTTACA
ATTGTTTGGTTCACCCTGTCTTTGGTACTATGGATTATCCGTTTGGTCCCTGATGTCATTAACCT
CTGCAGTCCGATGAATATGCATTGCTAACCAGAAATGTGGAGAGAGATAAATATGCAATTTCACTATT
AACTTTACAATGGAAAATCAGATTCATACTGGAATGGAATACGACAATGGCAGATTCATAGGGGTCAAG
TTCAAATCTGTAACCTTCAAAGACTCTGTTTTTAAAGTCTGCACCTTTGAGGATGTAACCTCAGTGAAC
ACCTACTTCAAGAACTGCACATTTATTGACTGTTTTTACAACACAGATTTTGGCCATATAAATTC
ATTGACAGTGAATTTAAAACTGCTCGTTTTTACAACAAGACGGGATGTCAGATTACCTTTGATGAT
GACTATAGTGCCTACTGGATTTATTTTGTCAACTTTCTGGGGACATTGGCAGTATTGCCAGGGAACATT
GTGCTCGCTCTGCTGATGGACAGAATTGGGCGCTTAACAATGCTAGGTGGCTCTATGGTGTCTTCGGGG
ATCAGCTGTTTCTTCTTTGGTTCGGCACCAGTGAATCCATGATGATAGGCATGCTGTGTCTGTACAAT
GGATTGACCATCTCAGCCTGGAACCTCTTGTACGTGGTCACTGTGGAACGTACCCACAGACCGGAGG
GCAACAGGCTTTGGCTTCTTAAATGCGCTATGCAAGGCAGCAGCCGCTCGGGAACTTAATATTTGGC
TCTCTGGTCAGCATCACCAAATCAATCCCCATCTGCTGGCTTCTACTGTGCTCGTGTGGAGGACTC
GTTGGGCTGTGCCTGCCTGACACACGAACCCAGGTTCTGATGTAA
ACGGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
  
```

Restriction Sites: SgfI-MluI

ACCN: NM_014979

Insert Size: 2184 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_014979.3
RefSeq Size:	10801 bp
RefSeq ORF:	2184 bp
Locus ID:	22987
UniProt ID:	Q496J9
Cytogenetics:	5q13.3
Protein Families:	Secreted Protein, Transmembrane
Protein Pathways:	ECM-receptor interaction
MW:	82.3 kDa
Gene Summary:	Plays a role in the control of regulated secretion in neural and endocrine cells, enhancing selectively low-frequency neurotransmission. Positively regulates vesicle fusion by maintaining the readily releasable pool of secretory vesicles.[UniProtKB/Swiss-Prot Function] Transcript Variant: This variant (1) encodes the longer isoform (1). Sequence Note: This RefSeq record was created from transcript and genomic sequence data to make the sequence consistent with the reference genome assembly. The genomic coordinates used for the transcript record were based on transcript alignments.