

Product datasheet for **SC304217**

AVL9 (NM_015060) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	AVL9 (NM_015060) Human Untagged Clone
Tag:	Tag Free
Symbol:	AVL9
Synonyms:	KIAA0241
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

Fully Sequenced ORF: >SC304217 representing NM_015060.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGCACTG
GATCCGGTACCGAGGAGATCTGCCGCCCGCATCGCC
ATGGAGAAGGCCAGGAGAGGCGGGGATGGCGTCCCCGGGGGCCCGTACTGCACATCGTGGTGGTCGGA
TTTCACCACAAGAAGGGCTGCCAGGTTGAATTCCTTACCGCCCTGATTCCAGGAGATGGACATGAC
AGCCACACTTTACCTGAAGAATGGAAGTATTTGCCCTTCCTTGCCCTACCAGATGGCGCACACAACATAC
CAGGAAGATACTGTGTTTTTCACTTGCCACCCAGAAATGGAATGGAGCCACAGTATTTGGTATCTCT
TGCTATCGACAAATTGAAGCCAAGGCACTGAAAGTAAGGCAAGCAGATATCACCAGAGAGACTGTTGAG
AAAAGTGTCTGTGTTCTAAGCAAGCTGCCTCTGTATGGTTTACTTCAAGCAAACTTCAACTCATTACA
CATGCATATTTTGAAGAGAAGGATTTTCCCAATTTCTATTCTAAGGAGCTTTATGAACATATGAAT
AGTTCCTTGGGAGGTGCTTCATTAGAAGGATCCCAAGTATATCTTGGTCTGTACCTCGAGATCTTGTC
CTTCATTTTCGACACAAGGTCTTAATCCTATTTAAGCTAATCTTCTTGAAGGAGGTTCTTTTTTAT
ATTTCTCCAGTGAATAAATTGGTGGTGCAGTGTGCTGTTATCCCTTTTCCAGGCATGATTGAA
CATGGTCTCAGTACTGTTCTCAGTATAGACCCCGGAAAAGTATGTCTGAAGATGGTGGGCTTCAGGAA
AGTAACCCATGTGCAGATGATTTGTTTCTGCATCCACTGCTGATGTTTACATACCAACTTGGGAAGT
ATCAGGAAAGTCAATGGCAGGAAACCATGGAGAAGATGCTGCCATGAAGACTGAGGAGCCTTTGTTCCAA
GTGGAAGACAGCAGCAAGGGCAGGAACCAATGATACCAATCAATATTTGAAACCTCCATCTCGCCCA
TCTCCAGATTCTCAGAAAGTACTGGGAACTTTGGATCCTAGTGTCTTAGAGGACCCCACTTGAAG
GAAAGGGAACAGCTGGGATCAGACCAGACAAATTTGTTTCCAAAGGACTCTGTCCCTCAGAGAGTCTT
CCAATTACTGTACAACCTCAAGCTAATACGGGACAGGTAGTCTGATACCAGGGCTCATTTCCGGTTTG
GAAGAGGATCAGTATGGCATGCCCTGGCCATCTTCAAAAGGATATCTGTGTTGCCTTACATGGCA
TTGCAGCAGCATCATCTCTCTCCGATGTCACCGTTTCGGGGTTTGTGCTGGAGCTACTAACATCCTT
TTTCGACAAACAGAAACACCTCAGTGATGCCATTGTGGAAGTAGAAGAAGCTCTGATCCAGATCCATGAT
CCAGAAGTCAAGGAGCTGCTTAACCCAACCACTGCAGACCTAAGGTTTCGAGACTACCTAGTGAGGCAC
GTGACTGAGAATCGGGATGACGTCTTCTAGATGGCACGGGCTGGGAGGAGGTGACGAATGGATCCGG
GCCAGTTTGGGCTCTACATTATGCCCTGCTGGCTGCCACACTGCAATTAGACAATGAAAGATATTA
TCGGACTATGGGACAACCTTTGTTACAGCATGGAAGAATACTCACAACACTACAGGGTGTGGAACAGCAAC
AAGCATCCAGCACTTGCAGAAATAAATCCAAACCATCCATTTCAAGGCCAATACTCAGTATCAGACATG
AAGTTAAGGTTCTCATTCTGTTCAGAAATAGTGAACGTGGCAAAAAATTGGAACGTCATGGTCACA
ACTAGCCGGAATGTTGTACAAACAGGAAAAGCTGTTGGCCAGTCAGTTGGAGGAGCTTTTCCAGTGCA
AAGACAGCTATGTCTTATGGCTTTCCACTTTCACCACTTCCACCTCCCAAGTCTCACTGAGCCACCA
GATGAGAAGCCTTGA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
  
```

Restriction Sites: Sgfl-MluI

ACCN: NM_015060

Insert Size: 1947 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_015060.2</u>
RefSeq Size:	7065 bp
RefSeq ORF:	1947 bp
Locus ID:	23080
UniProt ID:	<u>Q8NBF6</u>
Cytogenetics:	7p14.3
MW:	71.9 kDa
Gene Summary:	Functions in cell migration.[UniProtKB/Swiss-Prot Function]