

Product datasheet for **SC304233**

N4BP3 (NM_015111) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	N4BP3 (NM_015111) Human Untagged Clone
Tag:	Tag Free
Symbol:	N4BP3
Synonyms:	LZTS4
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

Fully Sequenced ORF: >SC304233 representing NM_015111.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGGCCACAGCCCCAGGCCCTGCTGGCATTGCCATGGGCAGCGTGGGCAGCCTGTTGGAACGGCAGGAC
TTCTCCCTGAAGAGCTCGGGCGGCACCTTGCCGGGTCTCGGGGCTCCCGCCAGCCTGATGGGCTCCTC
CGGAAGGGCTTGGGCCAGCGTGAGTTCTCAGCTACCTGCACCTCCCAAGAAGGACAGCAAGAGCACC
AAGAACACCAAGCGGGCCCTCGGAACGAGCCTGCCGACTATGCCACCCTCTACTACCGGAACATTCT
CGCGCGGGTGACTTCAGCAAGACCTCGCTGCCAGAACGGGGTCTGTTGACAAGTGCCGCAATTCGCCCC
TCAGTGTTCAAGCCTACGGCGGGCAACGGGAAAGGCTTCTATCCATGCAAAGTCTGGGTCCCACAAA
GGCCAGAAGCTGTGGCGCAGCAATGGCAGCCTGCACACGCTGGCCTGCCACCGCCCTGAGCCCGGG
CCCCGGGCCAGCCAGGCCCGGGCACAGTGCTGCACGCCCTCAGCCTAGATGAGGGCGGCCCTGAGCCC
GAGCCCAGCCTGTCCGACTCTCCAGTGGGGTAGTTTTGGTGCAGTCTGGTACTGGCCCTAGCCCC
TTCAGTCTCTCCCTTGGCCACCTTAACCACCTCGGGGGTCCCTGGACCGGGCCTCTCAAGGACCAAG
GAGGCTGGGCCACCAGCTGTGCTGAGCTGCCTGCCCGAGCCACCACCCCTACGAGTTCTCTGTCTCC
TCTGCCGAGGAAATGGGAGCCGTGCTGCCCGAGACCTGTGAGGAGCTCAAGAGGGGCTTGGCGATGAG
GACGGCTCCAACCCCTTCACGCAGGTGCTGGAGGAGCGCCAGCGGCTGTGGCTGGCTGAGCTGAAGCGC
CTGTATGTGGAGCGGCTGCACGAGGTGACCCAGAAGGCTGAGCGCAGCGAGCGCAACCTCCAGCTGCAG
CTGTTTATGGCTCAGCAGGAGCAGCGGCCCTGCGCAAGGAGCTGCGGGCTCAGCAGGGCCTGGCTCCG
GAGCCTCGGGCCCCCGGCACCTCCAGAGGCTGACCCAGTGCACGACCAGAGGAGGAAGCCCGATGG
GAGGTGTGCCAGAAGACAGCAGAGATTAGCCTCTTGAAGCAGCAGTGCCTGAAGCCAGGCGGAAGT
GCCCAGAAGCTGGCGGAGATCTTCACTGTAAGACACAATTCGGGGCAGCCGGGCACAAGCCAGGCT
CAGGACGCAGAGCTGGTCCGGCTGCGCGAGGCTGTGCGCAGCCTGCAGGAGCAGGCCCTCGGGAGGAA
GCCCCAGGCAGCTGTGAGACTGATGACTGCAAGAGCAGGGGCTGCTAGGGGAGGCAGGAGGCAGCGAG
GCCAGAGACAGTGTGAGCAGCTGCGGGCTGAGCTGCTGCAGGAGCGACTTCGGGGCCAGGAGCAGGCG
CTGCGCTTTGAGCAGGAGCGGGGACTTGGCAGGAGGAGAAGGAGCGCGTGTGCGCTACCAGCGGGAG
ATCCAGGAGGGTACATGGACATGTACCGCCGAACAGGCACTGGAGCAGGAAGTGCGGGCACTGCGG
GAGCCCCCACACCTGGAGTCCCCGGCTCGAGTCTCCAAGATCTGA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
  
```

Restriction Sites: SgfI-MluI

ACCN: NM_015111

Insert Size: 1635 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_015111.1</u>
RefSeq Size:	5987 bp
RefSeq ORF:	1635 bp
Locus ID:	23138
UniProt ID:	<u>O15049</u>
Cytogenetics:	5q35.3
MW:	60.5 kDa
Gene Summary:	Plays a role in axon and dendrite arborization during cranial nerve development. May also be important for neural crest migration and early development of other anterior structures including eye, brain and cranial cartilage.[UniProtKB/Swiss-Prot Function]