

Product datasheet for **SC304268**

CSTF2T (NM_015235) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	CSTF2T (NM_015235) Human Untagged Clone
Tag:	Tag Free
Symbol:	CSTF2T
Synonyms:	CstF-64T
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

Fully Sequenced ORF: >SC304268 representing NM_015235.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCACGATCGCC
ATGTCGAGTTTGGCGGTGAGAGACCCGGCAATGGATCGATCACTGCGTTCCGTGTTCTGTTGGGAACATT
CCATATGAGGCAACTGAGGAGCAGTTAAAGGACATTTTCTCGGAGGTTGGTTCTGTTGTGAGTTTCCGG
CTGGTATACGATAGAGAGACGGGAAAACCAAGGGCTATGGCTTCTGCGAATACCAAGACCAGGAGACC
GCGCTTAGTGCCATGCGGAACCTCAATGGGCGGGAGTTCACTGGGAGAGCGCTTCGGGTGGACAATGCT
GCCAGTAAAAGAATAAGGAGGAGTTAAAGAGCCTTGGGCTGCAGCGCCATTATTGACTACCCCTAT
GGGGATCCCATCGATCCAGAAGATGCCCTGAATCGATTACCAGAGCAGTAGCCAGTCTCCCCCGGAG
CAGATGTTTGAGCTGATGAAGCAGATGAAGCTCTGTGTCCAAACAGCCACCAGGAAGCTCGAAACATG
TTACTTCAAAATCCAACTGGCTTATGCACTGTTGCAGGCACAAGTAGTGATGAGAATCATGGATCCA
GAGATTGCTCTGAAAATTCTGCATCGGAAGATACATGTACACCACTGATCCAGGCAAACTCAGTCT
GTGCTGTCTCTGGCCCTGGCCCTGGCCCTGGGCTCTGCCAGGACCTAATGTTCTGCTGAAC
CAGCAGAAATCCTCCAGCTCCTCAGCCTCAGCATTTGGCTAGAAGACCTGTGAAGGACATTCTCCTCTG
ATGCAGACTCCTATCCAGGTGGAATTCAGCTCCAGGGCCAATACCAGTGCAGTTCGCGGAGCTGGT
CCTGGTTCCTTAATCCTGGAGGAGCAATGCAGCCCCAACTTGGAATGCCAGGGGTTGGCCAGTGCCT
TTAGAGCGGGGACAAGTGCAGATGTGAGATCCTAGAGCTCCTATACCTCGCGGACCCGTGACTCCTGGT
GGTCTGCCTCCTCGAGGACTGTTAGGAGATGCTCCAAATGACCCACGTGGAGGGACTTTGCTTTCAGTC
ACTGGAGAAGTGGAGCCCAGAGGTTATCTGGGTCCACCCCATCAGGGTCCCCCATGCATCATGCCTCT
GGTCATGACACTCGTGGCCCTTCTCACATGAGATGAGGGGAGGGCCATTAGGAGATCCCAGACTGCTA
ATTGGAGAGCCCAGAGGCCCATGATAGATCAAAGGGGTCTACCTATGGATGGTAGAGGTGGTAGAGAT
TCTCGAGCGATGGAGACTCGTGCCATGGAACCTGAGGTCTTAGAGACACGTGTAATGGAGAGGAGAGGA
ATGGAGACCTGTGCGATGGAACCAAGAGGGATGGAAGCAAGGGGCATGGATGCAAGAGGATTGGAGATG
AGGGGCCCTGTCCCAGTTCAGAGGCCCTATGACTGGTGAATTCAGGGTCTGGTCCCATTAAATATA
GGGGCAGGTGGCCCTCCTCAGGGACCCAGACAGGTCCCAGGCATTTAGGGGTGGGGAATCCTGGAGCT
GGTATGCAGGGTACAGGCATACAAGGAACAGGCATGCAGGGAGCAGGCATACAAGGAGGAGGGATGCAG
GGGGCAGGCATACAAGGAGTCAATACAAGGAGGAGGTACAAGGAGGAGGTACAAGGGGGCAAGC
AAGCAAGGTGGAAGCCAGCCTAGCAGTTTTAGTCTGGGAGAGCCAGGTCACTCCACAGGATCAGGAG
AAGGCAGCTTTGATCATGCAGGTTCTTCACTGACTGCAGATCAGATTGCCATGCTGCCCTGAGCAA
AGGCAGAGTATCCTGATTTTAAAGGAACAAATCCAGAAATCCACTGGAGCGTCTGA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
  
```

Restriction Sites: Sgfl-MluI

ACCN: NM_015235

Insert Size: 1851 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_015235.2</u>
RefSeq Size:	4127 bp
RefSeq ORF:	1851 bp
Locus ID:	23283
UniProt ID:	<u>Q9H0L4</u>
Cytogenetics:	10q21.1
MW:	64.4 kDa
Gene Summary:	May play a significant role in AAUAAA-independent mRNA polyadenylation in germ cells. Directly involved in the binding to pre-mRNAs (By similarity).[UniProtKB/Swiss-Prot Function]