

Product datasheet for **SC304309**

FAM149A (NM_015398) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	FAM149A (NM_015398) Human Untagged Clone
Tag:	Tag Free
Symbol:	FAM149A
Synonyms:	MST119; MSTP119
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

Fully Sequenced ORF: >SC304309 representing NM_015398.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGCACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGTTATTTGAAGGGAAGTGAACCTCAGACCCAGAGTCTGCTGGCCGAATGCGGGGAGTGGACAAGA
AGATCCCTCCATTTGAGAGTATTAGGAAGACAGCTGATCCTGCCCACTGACAAAGGCGTCCAGCATTTT
CAGGGCAGCACTCTGCCTCCGCAGTCCACAGACCCCGCTCAGTGCCTGCGGACACAGCAGCAACATC
AGAGAGTTGTGCAATTTCTGGCTCTCAAATAGTCCAGCAGCACTCTCAGCCTCTGCCCTGCCAGGCCCT
GATGACACAGGGTTGCTGACCTAACGGCACGTTCCCTGGAAGAAGAGGTTTACCATGTGGATGGA
AAGATTGAGGAGTATTTGCTTTTACAGAAAAGAGGATGATGATGAATGTCTTGAACAAAAACCAGCT
CAGCCCGTAGGAAATGGCGCAAACTCGGACTTCTCTGTTTCCCGCGTGACTGTGTCAAAGATGCC
GTGGCAGCAGAAGTGTGATCACGTCTGGACAAATATGGTAGAACTTTTGAAGAGCTGATTAGAAAA
CACTGGGAACTACACTCACAGAAGGGAAAAACAGAGAGAAACATTGAAAGTGGCTGGAACAGATTT
CCGCACGTCTCTGTTCCACACGCTCACGCCGATGGAGCCAGTGGCCCCCGTCCGGACACGCCGAGGCT
CAGGGCATCTCCCTGGCTTCTCGTCTGAACCCGCCCCAGATTCATCACTTCTCCAGCAGTTTTATAGT
GACATGAATGGTGTCTGACAATTCAAGCAAAACCGCTTCAGCGGAGACCTGCCTATTTTGTGACAGA
ACGCAAGATGAGAAGGAGGACAAAGCATCGGGTGGAGGGGAGGTGCTCTCTCCCGCACCCGACAGA
CTGGGACGGGCCCTCAGACACTCATGGATTATCACCTTCTGCAAGAAAACACCAAGTGCCTGGAGGCTG
CCTTCTCTTGTCTCAGATTCACAGAGACTAAAACTCCCAACATCTATAGTGACGAAGTTCTTGGGGA
ACAAAAGTCCGACTGGCGTGGACCACATGGCTTCCCCACTGGTTCAAACGTACGGAGCAGGTTCCCC
CCGCTAGTCACGGAGACCAGGGGGCAGAATACAGCAGTTCTGGATGCCGCTTGTCTTACAGAGGA
AGGCATCTACAAAACCGTGTGTTGAGTGCCATGCCTGACGGTACAGAACGATCGCGTCTTCGAGAAAGA
ACAGCCACCCTGGAACGGTTGTCAAGGCCAGCACAACCCACAGTTCCGGCAGTCGGATACGCCCTCGA
AAAAGTTCAATTGACACAAATGGAATTTGCTGCTCACACATGGACAGGTCAAAGTATTTTACAGGTTCA
CAATATGTGCCTAAATCTTTTACAGAGACAACCTTGACTTTCAAGAGGAGATTCCAAGTGACATCTGA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGCGC
  
```

Restriction Sites: SgfI-MluI

ACCN: NM_015398

Insert Size: 1449 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_015398.2](#)

RefSeq Size: 2721 bp

RefSeq ORF: 1449 bp

Locus ID: 25854

UniProt ID: [A5PLN7](#)

Cytogenetics: 4q35.1

MW: 53.3 kDa