

Product datasheet for **SC304314**

ZNF473 (NM_015428) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	ZNF473 (NM_015428) Human Untagged Clone
Tag:	Tag Free
Symbol:	ZNF473
Synonyms:	ZFP100; ZN473
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

Fully Sequenced ORF: >SC304314 representing NM_015428.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGCACTG
GATCCGGTACCGAGGAGATCTGCCGCCCGATCGCC
ATGGCTGAGGAATTTGTGACCCCAAGGATGTCGGCATGGACTTCACCTTGGGAGACTGGGAGCAGCTC
GGGCTGGAACAGGGGGACACGTTCTGGGACACAGCGTTGGACAATTGCCAGGACCTCTTCTGCTGGAC
CCCCAAGACCCAACTGACCTCCACCCAGATGGCAGTGAAGATCTGGAGCCTCTGGCAGGAGGAAGC
CCAGAAGCAACAAGCCCTGATGTGACTGAGACCAAGAAGCTCTCTCTGATGGAGGATTTCTTGAAGAA
GGATTCTCCAGGAGATTATAGAGATGTTATCCAAGGATGGCTTCTGGAAGTCCAATTTCCGAGAAGCC
TGTATAGAGGACACCTGGTTAGATAGTTTCTAGGCGATCCAGAAAGTCTTCTGAGGTCTGATATTGCC
ACCAACGGGGAAAGTCCACGGAATGCAAGAGTCATGAATTAAGAGAGGACTCAGTCTGTGTCCACC
GTTTCCACGGGAGAAGATTCCATGGTGCATAATGTTTCTGAAAAGACCCTCACACCAGCTAAGTCTAAG
GAATATAGGGGTGAGTTTTTCTCTACTCCGACCACAGCCAGCAGGATTCTGTTCAAGGAGGGAGAAA
CCATATCAATGTAGTGAATGTGGGAAAAGCTTCAGTGGGAGTTACCGTCTTACCCAGCACTGGATCACT
CATACTAGGGAGAAACCACTGTCCATCAAGAGTGTGAGCAAGGTTTTGACCGGAATGCTCCCTTTCT
GTGTATCCGAAAAGTCAACGCGGTACAAATTCTATGTGTGAATGAATATGGGACAACCTTTAGTCAG
AGTACATACCTGTGGCATCAGAAAAGTCACTGGAGAAAAACCATGTAAGAGTCAAGATAGTGACCAC
CCACCCAGTCATGACACACAGCCTGGTGGAGCATCAGAAAAGTCAACAGATAGTAAGTCTACAAGTGT
AACGAATGCGGCAAGGCTTTTACCCGGATCTTCCACCTTACTCGGACCCAGAAGATCCACACTCGGAAA
CGCTATGAGTGTTCAGTGGCAGGCGACCTTCAACTTGAGAAAACACCTCATCAACATCAGAAAAGT
CACGCTGCAAAAAGTACCTCTGAGTGTGAGGAGTGTGGGAAAGATTTTAGGCACAGTTCGCTGCTCATT
GAACACAGGCTCTTCTGCTGGAGAGGAGCCTTATAAGTGTAAACGAAGTGGGAAATCCTTCAGGCAT
AACTCTACCTAAAGATCCATCAGAGGGTTCACAGTGGAGAGAAGCCTTACAAATGCAGTGAAGTGTGGG
AAGGCCTTCCACCGGCACACTCACCTTAATGAACATCGGCGAATTCATACAGGCTACAGACCCCAAAA
TGTGAGGAATGCGTCAGGAGTTTCAAGCGGCCCTCACATCTGATGCGACATCAGGCCATTACACCGCA
GAAAAGCCCTATAGCTGTGCTGAATGCAAGGAGACTTTCAAGGATAACAATCGCCTTGTGCAACACCAG
AAAATGCACACTGTCAAAACCCCATATGAATGTCAGGAGTGGGAGAACGCTTCATTGCGGCTCAACC
CTGAAGTGCCACGAGAGTGTTCAGCCAGAGAAAAACAAGGATTTTGTGAGTGGGAAAGATCTTGAT
CAGAACCCAGAACAGAAAGAGAAGTGTCTTAAGTGTAAACAAATGTGAGAAAACCTTAGCTGCAGCAAA
TACCTAACTCAGCAGAGAGGATTCACACCAGGGGAGTGAAGCCCTTTGAATGTGACCAAGTGTGGGAAA
GCCTTTGGCCAAAGTACTCGGCTCATTACCATCAAAGAATCCACTCTAGAGTGAGGCTGTATAAATGG
GGTGAGCAAGGGAAAGCCATCAGCAGTGCTCCCTTATCAAACTTCAGTCTTCCACACAAGGAGCAC
CCTTTTAAATGTAAAGATGCGGAAAGACCTTCAGCCACAGTGCACACCTCTCAAAACATCAGTTAAT
CACGCTGGAGAGAATCCCTTTAAATGTAGTAAGTGTGACAGAGTCTTACCCAGAGAAATACCTTGTT
CAGCATGAGCGAACTCATGCCAGAAAGAAGCCGTTGGTGTGAACGAATGCGGGAAGACGTTCCGTCAG
AGCTCATGCCTTTCTAAGCATCAGAGAATTCACCTCAGGTGAGAAGCCCTATGTATGTGATTACTCGGG
AAGGCCTTCGGCTGAGTGTGAGCTTGTCCGCCACAGAGAATTCACACTGGAGAAAAGCCTTATGTT
TGTGAGGAATGCGGGAAGCCTTCACCCAGAGCTCATGCCTTTCTATTACCCGAGAGTTTCACTGGG
GAGAAGCCCTACAGATGTGGTGAATGTGGGAAAGCCTTTGCCAGAAAGCAATCTAACACAGCACCAG
AGAATTACACAGGGGAGAAGCCTTACTCTGTATGTGTGGCAAAGCTTTTGTCTCAGTGCCCAT
CTCAACCAGCACCTGAGAGTTTACACCCAGGAGACACTTTATCAGTGTCAACGTTGCCAGAAAGCCTTT
CGGTGCCACTCGAGCCTCAGCGCCATCAGCGTGTACACAACAAGCAGCAATACTGCCTGTAG
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
  
```

Restriction Sites: SgfI-MluI
ACCN: NM_015428
Insert Size: 2616 bp

OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_015428.3</u>
RefSeq Size:	4773 bp
RefSeq ORF:	2616 bp
Locus ID:	25888
UniProt ID:	<u>Q8WTR7</u>
Cytogenetics:	19q13.33
MW:	100.2 kDa
Gene Summary:	This gene encodes a member of the Krueppel C2H2-type zinc-finger family of proteins. The encoded protein, a component of the U7 snRNP complex, plays a role in histone 3'-end pre-mRNA processing and may be required for cell cycle progression to S phase. Expression level and methylation status of this gene may be correlated with bone mineral density. [provided by RefSeq, Jul 2016] Transcript Variant: This variant (1) represents the longest transcript and encodes the longer isoform (1). Variants 1 and 2 encode the same isoform.