

Product datasheet for **SC304332**

SAMD4A (NM_015589) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	SAMD4A (NM_015589) Human Untagged Clone
Tag:	Tag Free
Symbol:	SAMD4A
Synonyms:	SAMD4; SMAUG; SMAUG1; SMG; SMGA
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

Fully Sequenced ORF: >SC304332 representing NM_015589.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCCGGATCGCC
ATGATGTTTCGCGACCAGGTCTGGGGTCTGGCGGGCTGGTTTAAGGGCTGGAACGAGTGCAGCAGACT
GTTGCGCTGCTGCTGCTCAAGCGCTGAGCCAGACCAGGCCGCTTCTCCAGCTCTGCCTGGAG
CACTCGCTGGCCGACTGCGCCGAGCTGCAGTCTCGAACGCGAGGCCAACAGCCCCGGAATCATTAA
CAATGGCAACAGGAATCCAAGGATAAAGTGATTTCCCTCTGTTAACTCATCTGCCTTTGCTGAAGCCA
GAAACCTCGACGCGAAAGTAGAATATATGAAACTGCTGCCAAAATCCTGGCTCACTCTATTGAACAC
AACCAGCACATTGAGGAGAGCAGGCAGCTGCTGCTATGCTTTGATACATCCAGCCACTTCGTTAGAA
GACCGTAGTGCTTAGCCATGTGGCTGAATCACTTGGAGGACCGCAGCTCGACCAGCTTTGGTGGCCAG
AACCGAGGCCGCTCAGACTCTGTGGATTATGGACAGACACTACTATCACAAAGACAGAACTCTGAT
GACAAGCTCAATGGGTGGCAGAACTCTCGGGATTCTGGGATTTGCATCAATGCCTCCAAGTGGCAGGAC
AAAAGCATGGGGTGTGAGAAATGGCCATGTGCCCTCTACTCCTCCTCATCTGTCCCACCAATCAAT
ACGATTGGAACAGCACAAAGTACAATTCTCTCAGGCCAGGCACACCAGCCCTTTGAAACGATCTGTG
TCCCTTACCCACCCATGAATGTGCCAAACCAGCCTCTAGGACATGGATGGATGTCTCATGAGGACTTA
CGAGCTAGAGGACCCAGTGCCTCCCATCCGATCATGCCCCCTGTCTCCACAAAGCAGTGTAGCCTCT
TCAGGAAGTGGTGGGAGTGAACACTTAGAAGATCAGACCAGTCTCGTAACACATTTCAAGAAGAAGGT
AGTGGGATGAAAGATGTTCCAGCCTGGCTGAAAAGCCTCCGCCTGCACAAATATGCCGCGCTTTTCTCC
CAGATGACCTATGAGGAGATGATGGCCCTACCCGAGTCCAGCTGGAGGCGCAGAATGTTACAAAGGT
GCAAGACACAAAATTGTCATCAGTATTCAGAAGCTCAAAGAAAGACAAAATCTCCTGAAGCTTTGGAA
AGGGACATCATCGAGGGGGCAGCCTCCGCATCCCGCTCCAGAACTGCACCAGATGATCCTGACTCCG
ATCAAGGCCTACAGCTCCCCGAGCACCACCCCGAGGCTCGCCGCGGGAGCCCCAGGCCCGCGTCTAG
CCCTCACTGATGGGCCCGAGAGCCAGAGCCCGACTGCAAAGATGGGGCCGAGCCACTGGCGCCACG
GCCACCCCTCGGCCGGGGCAGCGGGGGCTCCAGCCGACCCAGCTGAGCAGCTGCGATGGGGAGCTG
GCCGTCGCCCCCTGCCAGAGGGGGACCTCCCCGGCAGTTCACACGCGTATGGGGAAAGTGTGCACA
CAGCTCTTGGTCTCCAGACCTGATGAGGAAAATAAAGTTCCTATTTACAGCTCATAGACAAGTGTCTA
ATTCATGAGGCATTTACAGAGACACAGAAAAAAGATTGTTGTCATGGAAACAGCAGGTGCAGAAGCTC
TTTCGGTCTTTCCCTCGAAAACCTTTCTAGACATATCAGGATATCGACAGCAAAGAAATCGAGGCTTT
GGCAATCCAACCTCCCTCCGACGGCTGGCTCTGTGGCGGTGGCATGGGCAGACGGAACCCGCGCCAG
TACCAGATCCCCTCTCGGAACGTCCCTTCCGCCGCTGGGCCTTTGGGACCAAGTGGATTTCGTGACG
TCCAACAGCGCAACACCACAGCTACCCCAACCATCATGAAACAAGGAAGACAGAACCTGTGGTTTGCC
AACCCCGGGGGCAGCAATAGCATGCCAAGCCGACCCACAGCTCAGTCCAGAGGACCCGCTCGTGCCC
GTGCACACTTCCCACAGAACATGCTGATGTTCCAGCAGCCAGAATTCAGCTTCCCCTGACCGAACCT
GACATCAACAACAGGCTGGAGTCGTTGTGCCTCAGTATGACCGAACACGCCCTGGGAGACGGGGTTGAC
CGGACCTCCACCATCAG
AGCGGACCGACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGAT
ATCCTGGATTACAAGGATGACGACGATAAGGTTTAA
  
```

Restriction Sites: SgfI-RsrII

ACCN: NM_015589

Insert Size: 2157 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_015589.5
RefSeq Size:	6839 bp
RefSeq ORF:	2157 bp
Locus ID:	23034
UniProt ID:	Q9UPU9
Cytogenetics:	14q22.2
MW:	79.4 kDa
Gene Summary:	<p>Sterile alpha motifs (SAMs) in proteins such as SAMD4A are part of an RNA-binding domain that functions as a posttranscriptional regulator by binding to an RNA sequence motif known as the Smaug recognition element, which was named after the Drosophila Smaug protein (Baez and Boccaccio, 2005 [PubMed 16221671]).[supplied by OMIM, Mar 2008]</p> <p>Transcript Variant: This variant (1) represents the longest transcript and encodes the longest isoform (1).</p>