

Product datasheet for **SC304335**

LRIT1 (NM_015613) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	LRIT1 (NM_015613) Human Untagged Clone
Tag:	Tag Free
Symbol:	LRIT1
Synonyms:	FIGLER9; LRRC21; PAL
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)




[View online »](#)

Fully Sequenced ORF: >SC304335 representing NM_015613.
Blue=Insert sequence Red=Cloning site Green=Tag(s)

```
GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGAGGGTGGCATTAGGCATGCTGGCTCTTGGCCCTTGCGTGGCCCCCAGGCCGGGGCTTCTGC
CCCTCTCAATGCAGCTGCAGCTCCATATCATGGGTGATGGCAGCAAGGCCAGGACAGTAGTGTGCAAC
GACCCCGACATGACCCCTGCCCGCGCTCCATCCCCCGACACCTCCAGACTGCGCCTGGAGCGGACG
GCCATACGCAGGGTTCTGGCGAGGCCTTCAGGCCCTGGGCCGCTGGAGCAGCTGTGGCTGCCTTAC
AACGCCCTCAGCGAGCTCAACGCCCTCATGCTGCGGGGCTGCGACGCCTGCGGGAGCTGCGGCTGCC
GGGAACCGCTGGCCGCTTCCCTGGCGGGCTCAGGGACGCCCAAGCTGCGGCTGCTGGACTG
CAGGCCAACCGCTCTCGGCTGTGCCGCTGAGGCCGCGCGCTTCTGGAGAACCTCACCTTCTCGAC
CTCTCCAGCAACCAGCTGATGAGGCTCCCGCAGGAGCTCATCGTCTCTGGGCTCACCTGGAGACCGGT
ATCTTCTCTCCCGCCACCACCCAGCGGGTCTAGGGCTACAGGACAACCCCTGGGCATGTGACTGC
GACTCTATGACTGGTTCATCTTTGGATGGCTGGGCCCAAATTGGCCTTCATTGAGACTGAACTG
AGATGTGCCAGCCACGTAGCCTGGCCGAGTGGCCTTCAGCCAGCTTGAAGTGAAGTGGCAGGGC
CCAGAGCTCCATCCAGGAGTGGCCAGCATCAGGTCCCTTTGGGTGGCACAGCACTGCTACGCTGTGGA
GCTACTGGAGTCCCTGGGCCGAGATGAGCTGGAGGAGGCCAATGGCAGGCCCTTAATGGTACAGTG
CACCAGGAAGTCTCCAGTGACGGCACGAGCTGGACTCTGCTGGCCGCTGCAGTGTCCACCTTGAC
TCCGGAGACTACATCTGCCAAGCAAGAAGTCTCTGGGAGCCTCTGAAACTGTTATCTCCTTGATTGTC
ACTGAGCCACCGACTTCCACAGAACACAGTGGGAGCCAGGAGCACTATGGGCAAGGACAGGTGGCGGA
GGGGAAGCTGCTGTTACAACAACAAGCTGGTGGCCAGGCATGTCCCCAGATTCCCAAACCCGCTGTC
CTGGCCACTGGACCCTCTGTGCCAGCACAAAGGAGGAGCTGACCCTTGAGCACTCCAGATGGATGCC
CTGGGAGAGCTCTCTGATGGGCGGGCAGGACCCTCAGAGGCACGAATGGTGAAGTCTGTGAAGTGGTG
GGGACACTTACCACAGCGTGTCTTGGTGTGGAAGGCACCCAGGCTAAGAACAACAAGTGCCTTCACT
GTCCTCTACGCGGTCTTGGGCAGCACAGCATGCGGCGGGTGATTGTGCAGCCTGGGAAGACCAGAGTG
ACCATCACTGGGCTGTTGCCAAGACCAAGTATGTGGCGTGTGTCTGTGTGCAGGGCCTGGTGCCTCCG
AAGGAGCAGTGTGTTATTTCTCCACCAATGAAGTGGTGGATGCTGAGAACACTCAGCAGCTTATCAAT
GTGGTGGTATCAGTGTGGCCATCGTATTGCCCTGCCTCTCACGCTGCTGTCTGCTGCAGTGTCTT
CAGAAGCGCTGCCAAAAGTGTCTCAACAAGGACTCCACTGAGGCCACAGTTACCTACGTCAACCTAGAG
AGACTGGGCTACAGCGAGGACGGCTTGGAGGAGCTGTCCCGGCACAGTGTGAGCGAGGCTGACAGGCTT
CTCTCAGCTCGTCCAGTGTGACTTTCAGGCCTTTGGAGTCAAAGGGGGCAGGAGAATCAATGAGTAC
TTCTGCTGA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
```

Restriction Sites: SgfI-MluI

Plasmid Map:


ACCN: NM_015613

Insert Size: 1872 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_015613.2](#)

RefSeq Size: 2301 bp

RefSeq ORF: 1872 bp

Locus ID: 26103

UniProt ID: [Q9P2V4](#)

Cytogenetics:	10q23.1
Protein Families:	Transmembrane
MW:	68 kDa
Gene Summary:	Possible role in phototransduction.[UniProtKB/Swiss-Prot Function]