

Product datasheet for **SC304336**

PYGO1 (NM_015617) Human Untagged Clone

Product data:


Product Type:	Expression Plasmids
Product Name:	PYGO1 (NM_015617) Human Untagged Clone
Tag:	Tag Free
Symbol:	PYGO1
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>SC304336 representing NM_015617. Blue=Insert sequence Red=Cloning site Green=Tag(s)

```
GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGCCCGCCGAGAAGCTCCAGCTCCCGCTTACAAAGTTTCTCGCATGGTGGTATAGTGGACTGGAT
GGGTTAGGAGGACCAGGTGTACAACCTAGGAAGCCAGATAAGAAAAAGCGCAAGGCAAATACACAGGGA
CCTTCTTTCCCTCCATTGTCTGAGTATGCTCCACCACCGAATCCAACCTGACCATCTAGTGGCTGCT
AATCCATTTGATGACAATAAATACTATTTCTATAAACCACCTACCTTCGTCAAATCCATATCTTGGC
CCTGGTTATCCTGGCTTTGGAGGCTATAGTACATTCAGAATGCCACCTCACGTTCCCCCAAGAATGTCT
TCCCATACTGTGGTCTTACTCACTCAGGAACCCAGCCACACCCATTTCTCAGAATCCTCTGGGCATG
GGTTTTAATCGACCTCATGCTTTTAACTTTGGGCCACATGATAATTCAAGTTTCGGTAATCCATCTTAT
AATAATGCACTAAGTCAGAATGTCAACATGCCTAATCAACATTTTAGACAAAATCCTGCTGAAAATTTT
AGTCAAATTCCTCCACAGAATGCTAGCCAAGTTTCTAACCCTGATTTGGCATCTAATTTTGTCTGGA
AATAATTCAAATTTTACTTCTCGTTAGAATCTAATCATTCTTTTATTCTCCCCAAACACTTTTGGT
CAAGCAAAGCACCACCCCAAAACAAGACTTTACTCAAGGAGCAACCAAAACACTAATCAAATTTCC
TCTGCTCATCCACCTCACTTGAATATGGATGACACAGTGAATCAGAGTAATATTGAATAAAAAATGTT
AATCGAAACAATGCAGTAAATCAGGAGAACAGCCGTTCAAGTAGCACTGAAGCCACAACAATAACCT
GCAAATGGGACGCAGAATAAGCCACGACAACCAAGAGGTGCAGCAGATGCCTGCACCACAGAAAAAAGC
AATAAATCCTCTCTTACCCAAACCGTCATGGCCATTCGCTTCTGACCCAGTGTATCCTTGTGGAATT
TGTACAAACGAGGTGAACGATGATCAGGATGCCATCTTATGTGAGGCCTTGTGCAGAAATGGTTTCAT
CGGATCTGTACTGGAATGACTGAAACAGCTTATGGCCTCTTAACTGCAGAAGCATCTGCAGTATGGGGC
TGTGATACCTGTATGGCTGACAAAGATGTCCAGTTAATGCGTACTAGAGAACTTTTGGTCCATCTGCA
GTGGGCAGTGATGCTTAA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
```

Restriction Sites: SgfI-MluI



Plasmid Map:



ACCN: NM_015617

Insert Size: 1260 bp

OTI Disclaimer: Due to the inherent nature of this plasmid, standard methods to replicate additional amounts of DNA in E. coli are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at custsupport@origene.com or by calling 301.340.3188 option 3 for pricing and delivery.

The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_015617.2</u>
RefSeq Size:	8399 bp
RefSeq ORF:	1260 bp
Locus ID:	26108
UniProt ID:	<u>Q9Y3Y4</u>
Cytogenetics:	15q21.3
Protein Families:	Druggable Genome
MW:	45.1 kDa
Gene Summary:	<p>Involved in signal transduction through the Wnt pathway.[UniProtKB/Swiss-Prot Function]</p> <p>Transcript Variant: This variant (1) encodes isoform 1. Variants 1 and 2 encode isoforms that are the same length, but have distinct protein sequences. Sequence Note: This RefSeq record was created from transcript and genomic sequence data to make the sequence consistent with the reference genome assembly. The genomic coordinates used for the transcript record were based on transcript alignments.</p>