

Product datasheet for **SC304357**

Cytokeratin 2 (KRT76) (NM_015848) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Cytokeratin 2 (KRT76) (NM_015848) Human Untagged Clone
Tag:	Tag Free
Symbol:	KRT76
Synonyms:	HUMCYT2A; KRT2B; KRT2P
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)




[View online »](#)

Fully Sequenced ORF: >SC304357 representing NM_015848.
Blue=Insert sequence Red=Cloning site Green=Tag(s)

```
GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCACGATCGCC
ATGAACAGACAAGTTTGCAAGAAATCCTTCAGTGGCAGGAGCCAGGGTTTCTCTGGCCGCTCTGCTGTG
GTCTCCGGCAGCAGCAGGATGAGCTGTGTGGCCCGCTCTGGGGAGCTGGTGGAGGGCCGTGGCTTC
AGGAGCGGAGCAGGCAGCTTTGGCAGTCGCAGCCCTCTACAACCTGGGCAGCAACAAGAGCATCTCCATC
AGCGTGGCAGCTGGCAGCTCCCGGGCTGGAGGCTTTGGGGAGGGCGGAGCAGCTGTGGCTTTGCAGGT
GGCTATGGAGGTGGCTTTGGGGCAGCTATGGTGGTGGCTTTGGTGGTGGCAGAGGAGTAGGTAGTGGT
TTTGGAGGGGCTGGTGGCTTTGGTGGAGCTGGTGGCTTTGGTGGTCTGGTGTATTTGGTGGCCTGGC
AGCTTTGGTGGTCTGGTGGCTTTGGCCCTGGGGCTTTCTGGGGAATTCAGGAAGTGATTGTAAC
CAGAGTCTCTGCAGCCCTCAATGTGGAGATCGACCCAGATTGGCAAGTAAAGGCCAGGAGCGG
GAACAGATCAAGACCCTCAACAACAAGTTGCCTCCTTCATCGACAAGTGCGGTTCTGGAACAGCAG
ACAAGGTCTGGAGACCAAGTGGAACTGCTCCAGCAGCAGACCACAGGCTCAGGGCCAGCAGCTG
GAGCCTTGTTTGAATCCTACATCAGCTTCTATGCAAGCAGCTAGATCACTTCTAGGGGAGAGGGGG
AACCTGGAGGGAGAGCTGAAAAGCATGCAGGACCTGGTGGAAAGACTTCAAAAAGAAGTATGAAGTGAA
ATCAACAAACGCACTGCCGAGAGAATGAGTTTGTGGGGCTCAAGAAGGATGTGGATGCGGCTTTTCATG
ACAAGGTGGAGCTGCAGGCCAAAGTGGACAGCCTGACAGATGAAGTCAAGTCTCTGAGGACCCTCTAT
GAGATGGAGCTGTCCAGATGCAAGCCATGCCAGTGACACGCTGTGGTCTGTCCATGGACAACAAC
CGCTGCCTGGACTGGGAGCATCATTGCCGAGGTCCGCGCCAGTATGAGGAGATTGCCAGAGGAGC
AAGTCTGAAGCTGAGGCCCTGTACCAGACCAAGCTTGGGGAGCTGCAGACCACAGCTGGCAGGCATGGG
GATGACCTGAGGAACACCAAGAGTGAGATCATGGAGCTCAACAGGATGATCCAGAGGCTACGGGCTGAG
ATTGAAAATGTCAAGAAGCAGAATGCCAACCTGCAGACGCAATTGCAGAGGCTGAGCAGCGTGGAGAG
ATGGCCCTCAAGGACGCAATGCCAAGCTCCAAGACTTGCAAGCTGCCCTACAGAAGGCTAAGGATGAC
CTGGCTCGGCTCCTGCGTACTACCAGGAGCTGATGAACGTCAAGCTGGCCCTGGATGTGGAGATTGCC
ACCTACCGCAAGCTGCTGGAGGGAGAGGAGTGCAGGATGTCTGGAGAATGTGAGAGTGGCGTGTGCATT
TCAGTGGTCAGCAATGTCACCAGCACAAGTGGCAGCTCTGGCAGTAGCCGTGGAGTTTTTGGAGGGGTC
AGTGGCAGTGGCAGTGGTGGCTACAAAGGCGGCAGCAGCAGCAGCAGCAGCAGTGGCTATGGAGTCAGT
GGCGGCAGCGGCAGTGGCTATGGAGGGTCAAGTGGCAGCACTGGGGCAGGGGTAGCAGCGGGAGC
TACCAGAGCAGCAGTAGTGGGAGCAGGCTCGGCGGTGCAGGTAGCATCTCCGTGAGCCACAGTGGAAATG
GGCTCCAGCTCTGGCAGCATCCAGACTTCTGGAGGAAGTGGCTACAAGTCTGGTGGTGGAGGCAGCACC
AGTATCCGCTTCTCCAGACCACGAGCTCAAGCCAGCATAGCTCCACCAAGTGA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGCGC
```

Restriction Sites: SgfI-MluI

Plasmid Map:


ACCN: NM_015848

Insert Size: 1917 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_015848.4](#)

RefSeq Size: 2529 bp

RefSeq ORF: 1917 bp

Locus ID: 51350

UniProt ID: [Q01546](#)

Cytogenetics: 12q13.13

MW: 65.8 kDa

Gene Summary: Keratins are intermediate filament proteins responsible for the structural integrity of epithelial cells and are subdivided into epithelial keratins and hair keratins. The type II keratins are clustered in a region of chromosome 12q13. [provided by RefSeq, Jun 2009]