

## Product datasheet for **SC304402**

### **FBXO40 (NM\_016298) Human Untagged Clone**

#### **Product data:**

Product Type:	Expression Plasmids
Product Name:	FBXO40 (NM_016298) Human Untagged Clone
Tag:	Tag Free
Symbol:	FBXO40
Synonyms:	FBX40
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)




[View online »](#)

Fully Sequenced ORF: >SC304402 representing NM\_016298.  
Blue=Insert sequence Red=Cloning site Green=Tag(s)

```
GCTCGTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCCGCATCGCC
ATGGGGAAAGCCCGCAGATCCCCGCCAGGGCACCACAGGCATTGTGAGGGATGCTTCAACCGCCACTGC
CACATTCCTGTGGAACCAACACCTCCTGCCTGGTAATAAGCTGCCACCTGCTCTGTGGTCCACCTTC
CACATGTGCAAAGAGGCAGAGCACCAGCTCCTCTGCCCTTAGAGCAGGTTCCGTGCCTCAACTCCGAA
TATGGCTGCCCTCTGTCCATGTCCCGCCACAACTGGCCAAGCACCTGCAGGTGTGCCCGCCAGCGTG
GTCTGCTGCTCCATGGAGTGAACCGCTGGCCAAATGTGGACTCTGAAACCACCCTTCATGAAAACATC
ATGAAAGAGACCCCGAGTGAGGAGTGTGGACACAGCCCTGGCCCTGCAGGATCAGAAGTCTCTTC
AGATCCTTGAATAAGTGGAACTTTCCAGAACTAGAGAGGCTACTGAGGAGGAACCAACTATGAAT
GGTGAACCGAGTGGAGGAAATGGGAGGAGCAGTGGTGGAGTGGATATCGTTTGGTACCACATGGT
CTGTCAGCAACTAATGGGAGATGGCAGAGCTAAGTCAAGAAGAACGGGAGGTGTAGCCAAAACAAA
GAAGGGATGGACCTGGTCAAGTTTGGCCAGTGGGAAAATATTTTCAGCAAAGAGCACGCAGCCTCTGCT
TTAACAAATTCATCAGCGAGCTGTGAGAGCAAGAACAAGAATGACTCCGAGAAAGAACAGATTTCCAGT
GGCCATAACATGGTAGAAGGAGAGGGCGCTCCCAAAAAGAAAGAACACAGGAAAATCAGAAGCAGCAG
GACGTTCTGACAGCCATGGAACACAGGGCTTGGCCCTTGGCAGGATGGTGTCTGAAAGACTGAAA
ACAGCTGTGGATGCAAAGGACTATAACATGTATCTAGTGACAATGGGCGGATGCTGATACACTTTGGT
CAGATGCCTGCTGTACACCCAAGGAGAGAGACTTTGTTTATGGCAAGCTGGAGGCTCAGGAAGTTAAG
ACTGTTTACACCTTCAAAGTTCCTGTGAGCTACTGTGAAAGCGAGCTCGACTTGGAGATGCCATGTTG
AGTTGTAAGCCAAGTGAACACAAGGCAGTGGATACTTCAGATTTGGGGATCACTGTGGAGGACCTGCC
AAATCAGATCTCATCAAGACCACCCTCCAGTGTGCTTTGGAAAGAGAAGTCAAAGGCCACGTCATCTCT
GAATCCAGAAGCATTGATGGACTGTTTATGGATTTTGGCACACAACATACAACCTTTGAGCCAGAACAG
TTTTCTCTGGGACAGTCTGGTGTGACCTAACCGCTGCCACCCAGGGGGACTCCACGTGGAGCTCCAC
AGCGAGTGTGTGACCAGGAGACACAACAAAGCAGCTCTGCCTTCACTTTCACTTGAACAAATTTCTT
AGGAGGGATGAGTTCCCTGCACTTCAAGAATGTCCACACAGACATTCACTCATGTCTCAATGGCTGG
TTCCAGCATCGATGCCCCCTCGCCTACTTGGGATGTACATTTGTTCAAACCAATTTCCGTCCCCAGGG
CAAAAGGCAAAAGTAATCTATAGCCAGGAGCTCAAGACCTTTGCCATTAAGCCGGAGGTTGCTCCAGAG
CTGAGCGAGGGAAGGAAGAACAACCATCTTTGGGTGATGGAGGAAAAAGCCAGAATTCTTAACCAGC
CTGCCCTGGAGATTTTGAAGTACATTGCTGGTCTTGGACAGCGTCAGCCTGGCCAGCTCTCCAG
GTGCTGTGCTGATGAGGAATATCTGTGCCACTTTGTTACAAGAGAGAGGAATGGTCTTTTGAATGG
AAGAAAAAGAGGTATCCCATGGAGGCACCTCCTGGAGAGTCCACAGAGAGATCTGGCAGTTCAGCAGC
CTCTTCTCAAAAATCAAGAGCTGGGAGTTAATGAAGTCACTCCATGTCTGAGCACCTGAAGTCTGT
CCTTTCAACATTGTAGAGCACAACAACTGACCCGATTCTTTGACTAGCATGTGTCAGCCCCGTGAGCAG
GCCCGAGAGAGCTTAGTCTCCACCTTTAGAATCAGACCACGAGGAAGATACGTCTCTAA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
```

Restriction Sites: SgfI-MluI

**Plasmid Map:**


**ACCN:** NM\_016298

**Insert Size:** 2130 bp

**OTI Disclaimer:** Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

**OTI Annotation:** This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

**Reconstitution Method:**

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

**RefSeq:** [NM\\_016298.3](#)

**RefSeq Size:** 5722 bp

**RefSeq ORF:** 2130 bp

**Locus ID:** 51725

**UniProt ID:** [Q9UH90](#)

**Cytogenetics:** 3q13.33

**Protein Families:** Druggable Genome

**MW:** 79.8 kDa

**Gene Summary:** Members of the F-box protein family, such as FBXO40, are characterized by an approximately 40-amino acid F-box motif. SCF complexes, formed by SKP1 (MIM 601434), cullin (see CUL1; MIM 603134), and F-box proteins, act as protein-ubiquitin ligases. F-box proteins interact with SKP1 through the F box, and they interact with ubiquitination targets through other protein interaction domains (Jin et al., 2004 [PubMed 15520277]).[supplied by OMIM, Mar 2008]