

Product datasheet for **SC304419**

ELOA2 (NM_016427) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	ELOA2 (NM_016427) Human Untagged Clone
Tag:	Tag Free
Symbol:	ELOA2
Synonyms:	HsT832; TCEB3B; TCEB3L
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)




[View online »](#)

Fully Sequenced ORF: >SC304419 representing NM_016427.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCACGATCGCC
ATGGCGGCAGGGTCCACTACGCTGCACGCAGTGGAGAAGCTGCAGGTGCGTCTGGCCACTAAGACGGAG
CCGAAAAAGCTAGAGAAATATTTGCAGAAACTCTCCGCTTGCCATGACGGCAGACATCTGGCGGAG
ACTGGAATCAGAAAGACGGTGAAGCGCCTGCGGAAGCACCAGCACGTGGGCGACTTTGCCAGAGACTTA
GCGGCCCGGTGGAAGAAGCTGGTCTCGTGGACCAGAAACACCCGCTGGCCACAGGACCTGAGGAG
AGCGCTTCCCGACAGCGCTTCGGGGAGGCTCTCAGGACCAGGAAAAGGCTGGGGCTTCCAGAAAAAC
GCGACGGCCCCAGGAGCCATCTCACAGCCTGAGCACAGACGGACAGCACGAGAACCTCCGGGG
CAACAGAGACCTACCCGAGGTCTCACAGTCGCGAGCCAGAGCTGAGAGAAAGTCCCCAGAATAGCC
CCAGCTGATTCGGCCGCTATCGGGCTCTCAACGCGCACAGCTCCCTCCGGATGCCGAGGGCCCT
GAGCCCCTGCGCCGGGAAGCAACCCGGAAGAGGCCACACTCACGCGGCTCAGGGCGGGCTCTGCTG
TGTCAGGCTGCCAGGGCAACCCAGGGAAAGCCGTTGTGAGCCACAGCAAGGGGCACAAATCGTCT
CGCCAGGAAAAACGCCCTTGTGTGCCAGGAGATTGGCACTCCCTACTTTGATCAGGGAGAAATCA
TGCGGGGCTGCTTAAGAGAGGAAACCCCAAGGATGCCCTCCTGGGAAGTGCCAGGGACAGGCAGCCT
TCGGACTTCAAGACAGACAAGGAAGGGGGCAAGCTGGCAGCGGCCAGCGTGTCCCTGCCTTGGAGGAG
GCTCCAGACAGTACCAGAAGAGGCCTCAGCACAGTCACTCGAACAGAAGAGGCCAGTCTAGACGGC
CGGGACCCAGGAAATGGGACACACGGCTGTGCGCCGAGGAGAAAGAGCAGCTTTCACACAGCCGAGAG
ACTCAAGAGGGGAAGCCACCGACTGCTCATTGGACAGAACGTCCTGAGCTCCCTCTCTGAGGTGGAG
GAGGTAGATATGGCTGAGGAATTCGAGCAGCCACTCTGTATGTAAAAATACCTCACCTACGATCAG
TTGCGGAAGCAAAAAGAAAAGACTGAAAAATCTGCCACCACTGCACTTGGAGATAAACAAAAGGAAAGCA
AACGAATCCAAGGGCACTCGTGAGTCTGGGATTGCGCTAAGAAATTGCCTCCTGTCCAGGAAAGCCAG
TCAGAGAGGCTGCAGGCGCCGGCGCTGATTCGCGGGCCGAAAAACGGTGCCAGCCATGTCTTCTCA
GAGCTCTGGACCTCTCAGAGGCCTGGATGCAGGCCAACTACGATCCGCTTTCGGATTCTGACTCCATG
ACCTCCCAGGCAAAGCCAGAAGCACTCTTCCACAAAGTTCCGGGAGGAAGCTGCTTCCCTGGACGC
AGAGTGAATGCTAAGATGCCGTGTACTCGGGCTCCAGGCCTGCCTGCCAGCTCCAGGTGCCGACGCTG
CGCCAGCAGTGTGCCAGGTGCTTAGAAACAATCCGGACGCCCTCAGCGACGTGGGAGAGGTCCCTAC
TGGTTCTTGAACCTGTTCTGGAAGGTGGAGGCCGATCAGCTGTATCGCAGAAAGAAAGACAATCAC
GCACTCGTTAGAGAGACAGACGAATTACGGAGGAATCATTGTTTCCAGGACTTCAAGGAAGAAAAGCCA
CAGGAAAAACAAACTTGGAGGGAGCAGTACCTGCGGCTTCCGGACGCCCCAGAGCAGCGCTGAGAGTA
ATGACAACGAATATCCGATCTGCACGTGGAACAACCCCAACGGCAGAGAGGCAAAGATGATCTGTTTC
AAATCTGTGGCCAAGACGCTTATGATACTTCAAGGAGGCAAGAGAAGTCTGCAGGAGACGCTGACCCC
GAAAATGGGGAGATCAAGCCAGCCTCAAGCCCAGGGAAGCAGCCACACTCCCTCCAGCCAGAGCAGC
AGCGGCGGTGGCAGAGACAGCAGCAGCAGCATCTTCGCTGGCTCCCTGAGAAGCGGGCCAAACCCCTGC
CTGAGCAGCAGCAATGAGCACGCGGCGCCCGGGCCAAAACCCGAAACAGGCTGCCAAGAAAGTGGCC
CCGCTGATGGCCAAGGCAATTCGAGACTACAAGAGAAGATTCTCCCGACGTAA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
  
```

Restriction Sites: SgfI-MluI

Plasmid Map:


ACCN: NM_016427

Insert Size: 2262 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_016427.2](#)

RefSeq Size: 3062 bp

RefSeq ORF: 2262 bp

Locus ID: 51224

UniProt ID: [Q8IYF1](#)

Cytogenetics: 18q21.1

Protein Families: Transcription Factors

MW: 83.9 kDa

Gene Summary: This gene encodes the transcriptionally active subunit of the SIII (or elongin) transcription elongation factor complex, which also includes two regulatory subunits, elongins B and C. This complex acts to increase the rate of RNA chain elongation by RNA polymerase II by suppressing transient pausing of the polymerase at many sites along the DNA template. Whereas a related protein with similar function, elongin A, is ubiquitously expressed, the encoded protein is specifically expressed in the testis, suggesting it may have a role in spermatogenesis. [provided by RefSeq, Jul 2008]