

## Product datasheet for **SC304444**

### CRNKL1 (NM\_016652) Human Untagged Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	CRNKL1 (NM_016652) Human Untagged Clone
Tag:	Tag Free
Symbol:	CRNKL1
Synonyms:	CLF; Clf1; CRN; HCRN; MSTP021; SYF3
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)




[View online »](#)

**Fully Sequenced ORF:** >SC304444 representing NM\_016652.  
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCCGGATCGCC
ATGACGGCAACAGTGGAGAACCTAACGTTCCAGAAGGACACTTTGGGAAATGCTGTTGACAAAAACACA
TCCAGATTGGAGCTCCGTTCTTACTCCCTTGCAGGGCGCCACGGCTCGACAGAGCCACTCGTCTCGCC
TGGTCATCCAGTCCGGAGGCTGACCTGGGGATGCGCCCTGGACGCCCTGCACCGCAGTCTTGGGTC
GCTGCTAGTCAGCACGGCGTTACCCACCTCATCCGCAGCTCCAGGACGCCGCACTCCACGCGCTGCCGC
AAGGAGGACGCGCAGCCTGGTCACCATGGAAACGGTGCCGCTCCGTACGGCTCAAGCGAGAGGTCAA
AGGTCAGTCTGCAAGTCCGCTCCGTTCCGAGATCTGTTTGTCTGAAAGCTTTGTGGTTTCG
GTGAGCTCTCAGAGCCGATTTCTAGCGTCCGTGCCGGGACAGGTGTCCAGAGTTCGACTGCTGCAGAC
ATGGCGGCTCCACCGCGCCGGGAAGCAGCGGATTCCCAAAGTGGCCAAGGTGAAAAACAAGCCCG
GCTGAGGTACAGATAACTGCTGAACAACCTTAAGAGAGGCTAAAGAAAGAGAATTGAGCTTCTCCA
CCTCCACCTCAACAGAAGATCACAGATGAAGAAGAATTAATGATTATAAACTAAGGAAAAGGAAGACT
TTTGAAGATAATAAAGAAAAACAGGACTGTGATTAGTAACTGGATAAAAATACGCACAATGGGAAGAA
AGCCTAAAGGAGATTCAAAGGGCTCGATCCATATACGAGCGTGCTTTAGATGTAGACTACCGAAATATT
ACACTCTGGCTGAAATACGCAGAAATGGAAATGAAGAATCGCCAAGTCAACCATGCTCGAAATATCTGG
GACCGGGCCATAACAACGCTGCCTCGAGTTAATCAGTTCTGGTACAAGTACACGTACATGGAGGAAATG
TTGGGAAACGTTGCCGGTGCCCGCAGGTGTTTGAGCGCTGGATGGAGTGGCAGCCTGAGGAGCAAGCC
TGGCACTCCTACATCACTTTGAGCTGAGATACAAAGAGGTGGATCGGGCCCGCACCATTATGAGCGA
TTTGTCTCGTGCACCTGATGTTAAGAAGTGGATCAAGTATGCCCGCTTTGAAGAAAAACATGCTTAT
TTTCCCATGCACGGAAAGTGTATGAGAGAGCTGTGGAATCTTTGGAGATGAACATATGGATGCAGAC
CTTTATGTTGCCTTTGCCAAGTTTGAAGAAAAACAGAAAGAGTTTGAAGGGTACGAGTGAATTAACAAG
TATGCCCTGGACAGAATTTCAAACAAGATGCCAAGAAGTCTTTAAAAATTATACCATCTTTGAGAAG
AAGTTTGGTGATAGGCGGGTATTGAAGATATCATTGTGAGCAAACGGAGATTCCAGTACGAAGAAGAA
GTGAAGGCGAATCCACACAATTATGATGCATGGTTTGATTACTTGCCTTGGTAGAAAGTACGCAGAA
GCTGAAGCCGTGAGAGAAGTCTATGAAAGGGCCATTGCCAATGTCCACCCATTGAGGAGAAAGGCAC
TGAAGCGCTACATTTATCTTTGGATCACTATGCACTCTATGAAGAATTGGAGGCAAAGGATCCTGAG
AGGACAAGACAGGTGTATCAAGCCTCTTTGGAAGTAACTCCTCACAAAAAGTTCACATTTGCCAAAATG
TGGATACTGTATGCACAGTTTGAATACGACAGAAAGTCTGTCTATTAGCCAGAAGAGCATTGGGAACT
TCCATAGGCAAAATGTCCAAAGAACAATTTAAAGTTTACATAGAATTGGAGCTACAGCTTCGAGAA
TTTGACAGATGCCGGAAGCTTTATGAAAAGTTCCTGGAATTTGGACCTGAAAATTGTACCTCATGGATT
AAATTCGCTGAATTAGAGACAATCCTTGGTGATATTGACAGAGCACGGGCAATCTATGAATTAGCCATC
AGTCAGCCACGTTTAGACATGCCAGAGGTGCTTTGGAAATCATATATTGATTTTGAATTTGAGCAGGAA
GAAACAGAAAGAACACGAAACCTTTACCGCGGTTGCTTCAACGGACGCAGCATGTCAAGGTATGGATC
AGCTTTGCTCAGTTTGAATGTCTTTCAGGAAAAGAAGGAAAGTTTACTAAATGCAGACAAATTTATGAA
GAAGCTAACAAAACCATGCGAAACTGTGAAGAAAAGGAAGAGAGACTTATGCTGCTGGAATCTTGGCGA
AAGAAGAGAAGAAAGTCCAGACTGATGATGGTCTGATGCAGGCTGGGAAGAATACTTTGATTACATC
TTTCCAGAAGATGCTGCCAACCAACCTAACCTCAAACCTCGCCATGGCCAAACTGTGGAAAGAACAG
CAGCAGGAAAAGGAGGATGCTGAGCACCATCCAGATGAGGACGTCGATGAGAGTGAATCTGA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
  
```

**Restriction Sites:** SgfI-MluI

**Plasmid Map:**


**ACCN:** NM\_016652

**Insert Size:** 2547 bp

**OTI Disclaimer:** Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

**OTI Annotation:** This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

**Reconstitution Method:**

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

**RefSeq:** [NM\\_016652.5](#)

**RefSeq Size:** 4413 bp

**RefSeq ORF:** 2547 bp

**Locus ID:** 51340

**UniProt ID:** [Q9BZJ0](#)

**Cytogenetics:** 20p11.23

**Protein Pathways:** Spliceosome

**MW:** 100.5 kDa

**Gene Summary:** The crooked neck (crn) gene of Drosophila is essential for embryogenesis and is thought to be involved in cell cycle progression and pre-mRNA splicing. A protein encoded by this human locus has been found to localize to pre-mRNA splicing complexes in the nucleus and is necessary for pre-mRNA splicing. Alternatively spliced transcript variants have been described. [provided by RefSeq, Jul 2013]  
Transcript Variant: This variant (1) represents the longest transcript and encodes the longest isoform (a).