

Product datasheet for **SC304460**

DEC1 (NM_017418) Human Untagged Clone

Product data:

Product Type: Expression Plasmids
Product Name: DEC1 (NM_017418) Human Untagged Clone
Tag: Tag Free
Symbol: DEC1
Synonyms: CTS9
Mammalian Cell Selection: Neomycin
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Fully Sequenced ORF: >SC304460 representing NM_017418.
Blue=Insert sequence **Red**=Cloning site **Green**=Tag(s)


```
GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGCCGGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGACAATGAATGTTCTGGAGGCTGGGAAGTGAAGAGCATTGTGCCAGCACCTGGTGAGGGCCTTCTT
GCCGTGTACACATGATGGTTTTACTGATGCCCTGCACAGAGAGAGGTCTGTAAAGTGGCAAGCAGGA
GTCTGCTACAATGGAGGAAAGGATTTTGTGTATCTCTTGCCAGGCCCAAGGCTGCAGAGGGAATTGCA
GATTGA
ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
```

Restriction Sites: Sgfl-Mlul



[View online »](#)

Plasmid Map:



ACCN: NM_017418

Insert Size: 213 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_017418.2](#)

RefSeq Size: 1251 bp

RefSeq ORF: 213 bp

Locus ID: 50514

Cytogenetics: 9q33.1

MW: 7.5 kDa

Gene Summary: The function of this gene is not known. This gene is located in a region commonly deleted in esophageal squamous cell carcinomas. Gene expression is reduced or absent in these carcinomas and thus this is a candidate tumor suppressor gene for esophageal squamous cell carcinomas. [provided by RefSeq, Jul 2008]