

Product datasheet for SC304516

STAP2 (NM_017720) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	STAP2 (NM_017720) Human Untagged Clone
Tag:	Tag Free
Symbol:	STAP2
Synonyms:	BKS
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>SC304516 representing NM_017720. Blue=Insert sequence Red=Cloning site Green=Tag(s)


```
GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTACTATAGGGCGCCGGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGGCCTCTGCCCTGAGGCCACCCCGTGTCCCAAGCCTAAGGGTGTCTGCCTCACACTACTATGAG
AGCTTTCTAGAGAAGAAGGGGCCCTGTGACCGGGATTACAAGAAGTTCTGGGCAGGCCTGCAGGGTCTC
ACCATTTATTTCTACAATAGCAATCGGGACTTCCAGCACGTGGAGAAGCTCAACTGGGAGCATTGAG
AAACTCACAGATGAGATTCCTGGGAAGCTCACGTGACCCTGGCACCCACTTCAGCCTGATTCTCCGG
GATCAGGAGATCAAGTTCAAGGTAGAGACCTTGGAGTGTGGAAATGTGAAAGGCTTCATCTTAACG
GTGGTGGAGCTCCGTGTCCGACCGACTTGACCCTGCTTCTGGGCACCTATACATGATGTCTGAAGTC
TTGGCCAAAGAGGAGGCCGCGCCGTGCACTGGAGACACCCTCGTGCTTCTGAAGGTGAGCCGGCTGGAG
GCACAACCTGCTCCTGGAGCGCTACCCCGAGTGGCGGAACCTGCTGCTGCGGCCACGCGGGACGGCGCC
GACGGCGTGTCCGTACCCACGCGGCAGATGCACAACGGGACGCACGTGGTCCGGCATTACAAGGTGAAG
CGGGAGGGCCCCAAGTACGTGATCGATGTGGAACAGCCGTTCTCTTGACCTCCCTGGACGCCGTGGTC
AACTATTTCTGTGTCGATACCAAAAAGGCGCTGGTGCCATTCTGTTAGACGAGGACTACGAGAAGGTG
CTAGGCTACGTGGAAGCCGATAAGGAGAATGGCGAGAATGTGTGGGTGGCGCCCTCCGCTCCGGGCCCA
GGTCTGCACCCCTGCACAGGTGGCCCCAAGCCGCTGTACCTGCGTCTAGCCAGGACAAGCTGCCCCCA
CTGCCCCACTACCGAACCAAGGAGAACTACGTGACCCCATTTGGAGATGGCCAGCTGTTGACTAT
GAGAACCAAGATGTGGCTTCTCTAGTTGGCCAGTCATCCTGAAGCCAAAGAAGTTGCCAAAGCCCTCT
GCCAAGCTTCAAAGCCACCCGTTGGACCAAGCCAGTAGAGAAGGGGTTTACCATGTTGCCAGGCT
GGTCTCGAACTCTGACCTCAAGTGATCCACCCACCTCAGCCTCCCAAAGTGTGGGATTACCGGCGTG
AGCCACCACACCTGGCCTCATCTGTCTTCTTCCAGAGCCCAAAGTCTTTAATGGTGGCTTGGGCAGG
AAGCTGCCAGTCAGTTCAGCCAGCCTCTTCCCCACAGCCGGGCTGGCAGACATGACGGCAGAGCTA
CAGAAGAAGCTGGAGAAGAGGGCGGCACTGGAGCACTGA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGCGC
```



View online »

Restriction Sites: SgfI-MluI

Plasmid Map:



ACCN: NM_017720

Insert Size: 1350 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_017720.2](#)

RefSeq Size: 1557 bp

RefSeq ORF: 1350 bp

Locus ID: 55620

UniProt ID: [Q9UGK3](#)

Cytogenetics: 19p13.3

MW: 49.6 kDa

Gene Summary: This gene encodes the substrate of breast tumor kinase, an Src-type non-receptor tyrosine kinase. The encoded protein possesses domains and several tyrosine phosphorylation sites characteristic of adaptor proteins that mediate the interactions linking proteins involved in signal transduction pathways. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Jul 2008]
Transcript Variant: This variant (1) represents the longer transcript and encodes the longer isoform (1).