

Product datasheet for **SC304552**

DET1 (NM_017996) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	DET1 (NM_017996) Human Untagged Clone
Tag:	Tag Free
Symbol:	DET1
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

Fully Sequenced ORF: >SC304552 representing NM_017996.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)


```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGGTGGGACAGATACTAAAAAGAGATGTGATAATGGATCATCATGTTTCTACCATCAAGCCTCGAAGA
ATCCAAAACCAAATGTCAATTCACCGCTTGAACGCGCGGATCAGTTCAGGCAAGGCAGGTACCCAC
TGGCACC AAGTCCGAGTGTCCATCAGAATGTCTTCCCAACTTCACAGTTGTCAACGTTGAAAAGCCT
CCTTGTTTCTTGCATAAATTCTCACCTGATGGACGCTACTTTATTGCTTTTTCTCAGACCAGACATCT
CTTGAAATCTATGAGTACCAGGGCTGCCAGGCAGCAGAGGACCTACTGCAGGGATACGAAGGAGAAATC
CTGTCCAATGGCAATGACCAGCGGTGAGTGAATATCCGGGGCCGGCTCTTTGAACGCTTTTTTGTCTG
CTGCACATTACCAATGTTGCGGCAATGGTGAGCACCTGAACCGGGAGTGTAGTCTTCTACTGATGAC
TGCCGCTGTGTCATCGTGGGCTCAGCTGCCTACCTCCAGATGAGCCTCACCTCCATTTTTGAGGTA
TATCGAACAGTGAATCAGTGACCCCAACCCACGGTCCCCTCTAGAAGACTATTCCTCCATATCATT
GACCTTACACCGCGCTTATGTGATACACGCAGTTCAAGTGTGACAAGGTGGTCTTGTACACAAAC
CAAGGGCTGACTTGTACAAAACATCCCTGGCCATCTTGTCTGTGCAACAACAGACCATCCATGTCTTC
CAGGTGACTCCTGAAGGCACTTTATTGATGTGCGGACCATTGGCCGCTTTTGTATGAGGATGACCTG
CTCACTGTGTCAGCTGTTTTCCCTGAGGTACAGCGGGACAGTCAGACAGGCATGGCCAATCCCTTTAGG
GATCCTTTCATCAATCCCTCAAACACCGGTTGCTGGTATATTTGTGGCGCCGGGAGAACAGGATGGT
AGTGCAATGGCCAAGAGGCGCTTCTCCAGTATTTGACCAACTGCGGCAGCTGCGAATGTGAAAAATG
CAGCTTCTGGATGAAAACCACTGTTTATCAAGTACACTAGTGAGGATGTAGTAACACTGCGAGTCACA
GATCCATCACAGGCATCTTCTTTGGTGTACAATATGGTGACGACAGAGGTGATTGCTGTGTTTGAG
AATACATCAGATGAGCTTTTGGAGCTCTTTGAGAATTCTGTGACCTTTTTCGTAATGCTACCTGCAC
AGTGAAGTTCAGTTTCCCTGCTCAGCTTCTAGCAACAATTTTGAAGGCAGATCCAGCGCCGTTCAA
GACACTATTATAAATGCCAAGTATGGAGGGCACACAGAGGCAGTACGCCGGCTGCTGGGTGAGTCCCC
ATCAGTGTCTCAGTCTTACAGCGGTAGCCCTATCTGGATTTGTCTCTTTCAGTTATGATGACAAGTGG
GTATCTGTGATGGAGCGGCCAAGACTTGTGGAGATCACCCAATCAGGTTCTATGCCCGGGACTCGGGC
CTGCTCAAGTTTGGATCCAGGCGGGTATTGGGCCGCCCATCAACCACACAGTGCAGCAGCTTGT
GCCTTACCTTTCACCTTTGAGCCTTTCGCTATTTCTGTGACAGGACTAATGCTGAGTATGTTGTC
AACTTCCATATGCGACACTGCTGCAGTAG
ACGCGTACGCGGCCGCTCGAGCAGAAACTATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTAAACGGCCGGC
  
```

Restriction Sites:

Sgfl-MluI

Plasmid Map:



ACCN:	NM_017996
Insert Size:	1686 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none"> 1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_017996.4</u>
RefSeq Size:	2382 bp
RefSeq ORF:	1686 bp
Locus ID:	55070
UniProt ID:	<u>Q7L5Y6</u>
Cytogenetics:	15q26.1
Protein Pathways:	Ubiquitin mediated proteolysis
MW:	65.1 kDa
Gene Summary:	<p>Component of the E3 ubiquitin ligase DCX DET1-COP1 complex, which is required for ubiquitination and subsequent degradation of target proteins. The complex is involved in JUN ubiquitination and degradation.[UniProtKB/Swiss-Prot Function]</p> <p>Transcript Variant: This variant (1) represents the longest transcript and encodes the longer isoform (1).</p>