

Product datasheet for **SC304585**

MBNL3 (NM_018388) Human Untagged Clone

Product data:


Product Type:	Expression Plasmids
Product Name:	MBNL3 (NM_018388) Human Untagged Clone
Tag:	Tag Free
Symbol:	MBNL3
Synonyms:	CHCR; MBLX; MBLX39; MBXL
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>SC304585 representing NM_018388. Blue=Insert sequence Red=Cloning site Green=Tag(s)

```
GCTCGTTTGTGAACCGTCAGAATTTGTAAACGACTACTATAGGGCGCCGGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGACGGCTGTCAATGTTGCCCTGATTCGTGATACCAAGTGGCTGACTTTAGAAGTCTGTAGAGAATTT
CAGAGAGGAACCTGCTCTCGAGCTGATGCAGATTGCAAGTTTGCCATCCACCAAGAGTTTGCCATGTG
GAAAATGGTCGTGTGGTGGCCTGTTTTGATTCTCTAAAGGGTCGGTGTACCCGAGAGAAGTCAAGTAC
CTTCACCTCCTCCACACTTAAAAACGCAGCTGGAGATTAATGGGCGGAACAATCTGATTCACAGAAG
ACTGCCGAGCCATGTTGCCAGCAGATGCAGCTTATGCTCCAAAACGCTCAAATGTCATCACTTGGT
TCTTTTCTATGACTCCATCAATTCAGCTAATCCTCCATGGCTTTCAATCCTTACATACCACATCCT
GGGATGGGCTCGTTCTGCAGAACTTGTACCAATACACCTGTTCTGATTCCTGAAAACCCACCTCTT
GCAATGCCAGGAGCTGTTGGCCAAAACCTGATGCGTTCAGATAAACTGGAGGTTTGCCGAGAATTTAG
CGTGGAATTTGTACCCGTGGGGAGAATGATTGCCGCTATGCTCACCTACTGATGCTTCCATGATTGAA
GCGAGTGATAACTGTGACAATCTGCATGGATTACATCAAAGGTCGATGCTCGCGGGAGAAATGCAAG
TACTTTTATCCTCCTGCACACTTGAAGCCAGACTCAAGGCGCTCATCATCAGATGAACCATTCAGCT
GCCTCTGCCATGGCCCTGCAGCCTGGTACACTGCAACTGATACCAAGAGATCAGCACTGGAAAAGCCC
AATGGTGCCACCCGGTCTTTAATCCACTGTTTTCCACTGCCAACAGGCTCTGACTAACCTGCAGCTC
CCACAGCCGGCATTATCCCTGCAGGGCAATACTGTGCATGGCACCCGCTTCAAATATTGTGCCATG
ATGCACGGTGCTACACCTACCCTGTGTCTGCAGCAACAACACCTGCCACCAGCGTTCCGTTCCGTGCA
CCAACACAGGAATCAGCTGAAATTCGTA
ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGCGCCGC
```

Restriction Sites: SgfI-MluI



[View online »](#)

Plasmid Map:


ACCN: NM_018388

Insert Size: 1065 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_018388.3](#)

RefSeq Size: 11502 bp

RefSeq ORF: 1065 bp

Locus ID: 55796

UniProt ID: [Q9NUK0](#)

Cytogenetics: Xq26.2

MW: 38.5 kDa

Gene Summary: This gene encodes a member of the muscleblind-like family of proteins. The encoded protein may function in regulation of alternative splicing and may play a role in the pathophysiology of myotonic dystrophy. Alternatively spliced transcript variants have been described. [provided by RefSeq, Dec 2009]
Transcript Variant: This variant (1, also known as variant G) encodes the longest isoform (1, also known as isoform G). Sequence Note: This RefSeq record was created from transcript and genomic sequence data to make the sequence consistent with the reference genome assembly. The genomic coordinates used for the transcript record were based on transcript alignments.