

## Product datasheet for **SC304596**

### **ARHGAP15 (NM\_018460) Human Untagged Clone**

#### **Product data:**

Product Type:	Expression Plasmids
Product Name:	ARHGAP15 (NM_018460) Human Untagged Clone
Tag:	Tag Free
Symbol:	ARHGAP15
Synonyms:	BM046
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)


**Fully Sequenced ORF:** >SC304596 representing NM\_018460.  
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAATACGACTCACTATAGGGCGGCCGGGAATTCGTGCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGCAGAAATCTACAAATCTGATACTCCGTGGAAACACTGAATTCTACCCGCCAAGGCACAGGAGCT
GTGCAAAATGAGAATCAAAAATGCCAACAGCCACCATGACAGGCTCAGCCAAAGTAAATCCATGATCCTC
ACCGATGTCGGGAAGGCTACTGAACCTATATCCAGACACAGAAGGAATCATTACAGCATATCTTGAAA
GATGTCATTCTCCATTGGAACAACTGATGGTTGAAAAAGAAGTTATCTGCAAAAAGCTAAAATTGCA
GATGGAGGAAAGAACTAAGGAAAACTGGTCTACTTCCTGGATTGTTCTTTCTAGTCGAAGAATTGAA
TTTTACAAAGAATCCAAGCAACAGGCTCTGTCCAATATGAAAACCTGGGCACAAACCAGAAAGTGGAT
TTGTGTGGAGCACACATTGAATGGGCCAAGGAAAAATCGAGCAGAAAGAATGTCTTTAGATCACAACA
GTATCAGGAAATGAGTTCCTTCTACAGTCAGATATTGACTTCATCATATTGGATTGGTCCACGCTATC
AAAAATGCAATTGACAGATTGCCAAAGGATTCAAGTTGTCCATCAAGAAACCTGGAATTATCAAATC
CAAAGATCCTCTAGCACTGAATTGCTAAGTCACTACGACAGTGATATAAAAGAACAGAAACCAGAGCAC
AGAAAATCTTTAATGTTCCAGACTGCATCACAGTCTCCGATACAAGCGACAAAAATCGAGTTAAAAGC
AGATTAAGAAGTTTATTACCCGAAGACCTTCCCTGAAAACCTCTGCAAGAAAAAGGACTTATTAAGAT
CAAATTTTGGCTCTCATCTGCACAAAGTGTGTGAACGTGAAAATTCACAGTTCCTGGTGTGTAAG
CAATGCATTGAAGCTGTTGAGAAAAGAGGTCTAGATGTTGATGGAATATATCGAGTTAGTGGCAATCTG
GCAACAATACAGAAGTTAAGATTTATTGTCAACCAAGAAGAGAAGCTGAATTTGGACGACAGCCAGTGG
GAGGACATCCACGTTGTCACCGGAGCACTGAAGATGTTTTCCGGGAGCTGCCTGAGCCGCTCTCCCT
TACAGTTTCTTTGAGCAGTTTGTGGAAGCGATCAAAAAGCAAGACAACAACAAGAATTGAAGCTGTA
AAATCTCTGTACAAAACCTCCCTCCGCAAATCGTGACACCATGAAAGTCTCTTTGGACATCTAAT
AAGATAGTGGCCAAAGCCTCCAAGAACCTCATGTCCACGCAAGCTTGGGGATTGTATTTGGACCTACC
CTTCTGCGAGCTGAAAATGAAACAGGAAACATGGCGATCCACATGGTCTACCAGAACCAGATAGCTGAG
CTCATGCTGAGTGAGTACAGTAAGATCTTCGGCTCAGAGGAAGCTGA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
  
```

**Restriction Sites:** SgfI-MluI

**Plasmid Map:**



**ACCN:** NM\_018460

**Insert Size:** 1428 bp

<b>OTI Disclaimer:</b>	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
<b>OTI Annotation:</b>	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<u><a href="#">NM_018460.3</a></u>
<b>RefSeq Size:</b>	1774 bp
<b>RefSeq ORF:</b>	1428 bp
<b>Locus ID:</b>	55843
<b>UniProt ID:</b>	<u><a href="#">Q53QZ3</a></u>
<b>Cytogenetics:</b>	2q22.2-q22.3
<b>MW:</b>	54.5 kDa
<b>Gene Summary:</b>	RHO GTPases (see ARHA; MIM 165390) regulate diverse biologic processes, and their activity is regulated by RHO GTPase-activating proteins (GAPs), such as ARHGAP15 (Seoh et al., 2003 [PubMed 12650940]).[supplied by OMIM, Mar 2008]