

Product datasheet for **SC304643**

UNC93A (NM_018974) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	UNC93A (NM_018974) Human Untagged Clone
Tag:	Tag Free
Symbol:	UNC93A
Synonyms:	dj366N23.1; dj366N23.2; Unc-93A
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>SC304643 representing NM_018974. Blue=Insert sequence Red=Cloning site Green=Tag(s)

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCC**CGATCGCC**
ATGGACAGAAGTCTAAGGAACGTCCTTGTGGTTTCCTTTGGGTTCTGCTTCTTTACAGCCTATGGA
GGTCTGCAGAGCCTGCAGAGCAGCCTGTACAGCGAGGAGGGCCTGGGTGTACAGCGCTCAGCACCTC
TATGGAGGCATGCTCCTGTCTCCATGTTCTCCACCGCTCCTCATCGAGAGGCTGGGTGCAAGGGG
ACCATCATCTCTCCATGTGTGGCTACGTGGCCTTCTCCGTGGGCAACTTCTCGCCAGCTGGTACACT
TTGATCCCCACCTCCATACTGCTGGGACTCGGGCGCGCCCGCTGTGGTCTGCACAGTGCACATACCTC
ACGATCACGGGAAACACACATGCAGAGAAGCGGGAAAGCGTGGCAAAGACATGGTGAACCAAGTATTTT
GGCATCTTCTTCTCATATTCCAGTCATCCGGTGTGTGGGGCAACTTGATCTCATCGCTGGTATTTGGC
CAGACTCCAGCCAAGAGACCCTTCCAGAAGAGCAGCTCACGTCTGTGGGGCCAGTGACTGCCTGATG
GCCACCACAACCAACAGCACCCAGAGGCCCTCCAGCAGCTGGTCTACACCCTCCTGGGCATCTAC
ACTGGGAGTGGTGTCTGGCTGTCTGATGATAGCTGCGTTCTCCAACCCATACGAGATGTTACAGCGG
GAAAGTGAAGGAGAGAAGAAATCAGTACCTTCTGGTCCACTTTACTGTCGACTTTCAAGCTATATAGA
GATAAACGTCTGTGCCTCTTAATTCTGCTGCCGTGTACAGTGGATTGCAGCAAGGATTCCTCTCCAGC
GAATACACAAGGTCTATGTCACCTGCACCCTGGGCATCCAGTTCGTCGGCTACGTGATGATCTGCTTC
TCGGCCACTGACGCGCTGTGCTCCGTGTTGTATGGAAGGTCTCGCAGTACACGGGCAGGGCTGTGCTG
TACGTGCTGGGCGGGTGACCCACGTGTCTGCATGATTGCCCTACTGCTGTGGAGACCTCGTGCTGAC
CATCTGGCAGTGTCTTCGTATTCTCTGCCTGTGGGGCGTGGCAGATGCCGTCTGGCAGACACAAAAC
AATGCTCTCTACGGCTTCTGTTTGAGAAGAGCAAGGAAGCTGCCTTCGCCAATTACCGCTGTGGGAG
GCCCTGGGCTTCGTATTGCCTTCGGGTACAGCATGTTTTTGTGCGTGCACGTCAAGCTCTACATTCTG
CTGGGGTCTCTGAGCCTGACCATGGTGGCGTATGGGCTTGTGGAGTGCCTGGAGTCCAAGAACCAGATC
AGACCCACGCTCCAGGACAGGTCAACCAGGCAGAGGATGAAGAAATACAAACAAAAT**GTGA**
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTAAACGGCCGGC



[View online »](#)

Restriction Sites:	Sgfl-Mlul
ACCN:	NM_018974
Insert Size:	1374 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_018974.3</u>
RefSeq Size:	2132 bp
RefSeq ORF:	1374 bp
Locus ID:	54346
UniProt ID:	<u>Q86WB7</u>
Cytogenetics:	6q27
Protein Families:	Transmembrane
MW:	50.3 kDa