

Product datasheet for **SC304648**

FBXO42 (NM_018994) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	FBXO42 (NM_018994) Human Untagged Clone
Tag:	Tag Free
Symbol:	FBXO42
Synonyms:	Fbx42; JFK
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)




[View online »](#)

Fully Sequenced ORF: >SC304648 representing NM_018994.
Blue=Insert sequence Red=Cloning site Green=Tag(s)

```
GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCCGCATCGCC
ATGGCCAGTCCTCGGACAGTGAAGATGACAGTTTCATGGCTGTGGACCAGGAAGAACTGTGCTGGAA
GGGACAATGGATCAAGATGAGGAGCCCCACCCAGTATTGGAGCTGAGGAGACTAGACATAATAGGTCC
ATGTCGGAGCTGCCAGAAGAGTTTTGGAGTATATCCTGTCTTTCTCTACCGTATCAGGAACACAAA
ACTGCGGCCCTTGTCTGCAAACAGTGGTATCGACTTATCAAAGGTGTAGCCCATCAGTGTTATCATGGT
TTCATGAAGGCTGTCCAGGAAGGAAACATTCAGTGGGAGAGCCGTACCTATCCTTATCCTGGAACCCCA
ATCACTCAGCGCTTCTCGCACAGTGCATGCTATTATGATGCTAATCAGTCTATGTATGTGTTGGAGGC
TGTACCCAGAGCAGCTGCAATGCTGCTTCAATGACCTCTGGAGACTTGACCTAAACAGCAAAGAGTGG
ATCCGACCTTTGGCTTCAGGGTCTATCCTTCCCCAAAGCTGGAGCAACTCTGGTCGTGTACAAGGAC
TTGCTAGTGCTGTTTGGTGGCTGGACGCGGCCAAGCCCTTATCCCCTACACCAGCCAGAGAGATTCTTT
GATGAAATACACACTTACTCACCTCTAAAAATTGGTGGAACTGCATTGTGACAACCCATGGGCCACCT
CCCATGGCTGGCCACTCCTCCTGTGTGATAGATGATAAAATGATTGTCTTTGGTGCTCTTAGGATCC
CGGCAATGAGCAATGATGTCTGGGTCCTTGACCTTGAGCAGTGGGCGTGGTCCAAGCCGAACATCTCT
GGCCCCAGTCCTCATCTCGAGGTGGCCAACTCAGATTGTATAGATGATGCAACTATCTTAATCCTC
GGAGGGTGTGGCGGTCCCAATGCTCTATTCAAGGATGCTTGGTGTGTCACATGCATTCTGGTCCTTGG
GCCTGGCAGCCACTCAAGGTAGAAAAATGAAGAGCATGGGGCCCCAGAAGTGGTGGCCATCCAGCTTGC
CGGGTGGACAGTGTGTGGTGGTCTTCAGCCAGGCTCCTAGTGGGAGAGCCCCACTCAGCCCCAGTTTG
AACTCTCGCCATCACCTATCAGTGCCACTCCTCAGCTCTCGTTCCTGAAACCCGAGAGTACCGCTCT
CAGTCTCCAGTAAGAAGCATGGATGAAGTCCCTTGTGTTAACGGCCGCTGGGGAACACTGAGACCCAGG
GCTCAAAGGCAGACTCCTCAGGTTCCCGGGAAGGGAGCCTTCCCCAGCCAGAGGAGACGCTCTCCT
ATCCTCAATGGTGGGAGTTTGTCTCCAGGAACGGCAGCTGTGGTGGCTCTTCTTTGGACAGTCTGTA
CAGGCCATATCTCCAAGTACTCCATCTGCTCCTGAAGGATACGACCTGAAAATAGGACTTCTTTGGCC
CCCCGACGAGGATCACTACCAGATCAGAAAGATCTGAGATTAGGATCCATAGATCTGAATTGGGATCTG
AAACCCGCTTCCAGTAGTAATCCCATGGATGGCATGGACAATAGGACAGTTGGGGGAAGTATGAGACAC
CCTCCTGAACAGACAAAATGGTGTGCATACCCACCTCACGTGGCCAGTGCCCTTGACGGGCGCTCTCC
CCAGGTGCCCTGCGTCGGAGTCTGGAAGCCATCAAAGCGATGTCCTCAAAGGCCCTCGGCCTCTGCA
GCTAAGTCTCCTCTTGGGTCTTCTCCAGGCTCTCCTGGGAGCCAGAGTTTGAGCAGTGGAGAAACA
GTGCCATCCCTCGCCAGGGCCTGCCAAGGAGATGGACATTCTTACCTCCCATTGCTCGCCGCTG
GGCCACCACCCTCCACAGTCCCTAAATGTTGGCAAACCCCTATACCAGAGTATGAACTGCAAGCCCATG
CAGATGTACGTGCTGGACATTAAGACACCAAGGAGAAGGGGCGGGTCAAATGGAAGTATTTAATAGC
AGTTCTGTGGTTGGACCTCCTGAAACCAGCCTGCATACCGTGGTACAAGGCAGGGGTGAACATCATA
TTTGGAGGACTCATGGACAAGAAACAGAATGTGAAGTACTATCAAAAACAAACGCTTGTACTTTGTA
CGAGCAAAGAGATAA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
```

Restriction Sites: SgfI-MluI

Plasmid Map:


ACCN: NM_018994

Insert Size: 2154 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_018994.2](#)

RefSeq Size: 6253 bp

RefSeq ORF: 2154 bp

Locus ID: 54455

UniProt ID: [Q6P3S6](#)

Cytogenetics:	1p36.13
Protein Families:	Druggable Genome
MW:	77.8 kDa
Gene Summary:	Members of the F-box protein family, such as FBXO42, are characterized by an approximately 40-amino acid F-box motif. SCF complexes, formed by SKP1 (SKP1A; MIM 601434), cullin (see CUL1; MIM 603134), and F-box proteins, act as protein-ubiquitin ligases. F-box proteins interact with SKP1 through the F box, and they interact with ubiquitination targets through other protein interaction domains (Jin et al., 2004 [PubMed 15520277]).[supplied by OMIM, Dec 2010]