

Product datasheet for **SC304650**

CSGLCAT (CHPF2) (NM_019015) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	CSGLCAT (CHPF2) (NM_019015) Human Untagged Clone
Tag:	Tag Free
Symbol:	CSGLCAT
Synonyms:	chPF-2; ChSy-3; CSGLCA-T; CSGLcAT
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

Fully Sequenced ORF: >SC304650 representing NM_019015.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCCGCATCGCC
ATGCGACTGAGCTCCCTGTTGGCTCTGCTGCGGCCAGCGCTTCCCCTCATCTTAGGGCTGTCTCTGGGG
TGCAGCTGAGCCTCTGCGGGTTTCTGGATCCAGGGGAGGGAGAAGATCCCTGTGTCGAGGCTGTA
GGGGAGCGAGGAGGGCCACAGAATCCAGATTCAGAGCTCGGCTAGACCAAAGTGATGAAGACTTCAA
CCCCGATTGTCCCCTACTACAGGGACCCCAACAAGCCCTACAAGAAGGTGCTCAGGACTCGGTACATC
CAGACAGAGCTGGGCTCCCGTGAGCGGTTGCTGGTGGCTGCTGACCTCCCGAGCTACACTGTCCACT
TTGGCCGTGGCTGTGAACCGTACGGTGGCCATCACTTCCCTCGTTACTCTACTTCACTGGGCAGCGG
GGGGCCCCGGCTCCAGCAGGGATGCAGGTGGTGTCTCATGGGGATGAGCGGCCCGCTGGCTCATGTCA
GAGACCTGCGCCACCTCACACACACTTTGGGGCGACTACGACTGGTCTTTCATCATGCAGGATGAC
ACATATGTGCAGGCCCGCCCTGGCAGCCCTTGTGGCCACCTCAGCATCAACCAAGACCTGTACTTA
GGCCGGGCAGAGGAGTTCAATTGGCGCAGGCGAGCAGGCCGGTACTGTCATGGGGCTTTGGCTACCTG
TTGTCACGGAGTCTCTGCTTCTGCTGCGGCCACATCTGGATGGCTGCCGAGGAGACATTCTCAGTGCC
CGTCTGACGAGTGGCTTGGACGCTGCCTCATTGACTCTCTGGGCGTCCGCTGTGTCTCACAGCACCAG
GGGCAGAGTATCGCTCATTGAACTGGCCAAAAATAGGGACCCTGAGAAGGAAGGGAGCTCGGCTTTC
CTGAGTGCCTTCGCCGTGACCCCTGTCTCCGAAGGTACCCTCATGTACCGGCTCCACAAACGCTTACG
GCTCTGGAGTTGGAGCGGGCTTACAGTGAATAGAACAAGTGCAGGCTCAGATCCGGAACCTGACCGTG
CTGACCCCGAAGGGGAGGCAGGGCTGAGCTGGCCGTTGGGCTCCCTGCTCCTTTCACACCACACTCT
CGCTTTGAGGTGCTGGGCTGGACTACTCACAGAGCAGCACACCTTCTCCTGTGCAGATGGGGCTCCC
AAGTGCCCACTACAGGGGGTACGAGGGGAGCTGCGGTGATGCGTTGGAGACTGCCCTGGAGCAGCTC
AATCGGCCGATACAGCCCGCCTGCGCTTCCAGAAGCAGGACTGCTCAACGGTATCGGGCTTCGAC
CCAGCACGGGCGATGGAGTACACCCTGGACCTGCTGTTGGAATGTGTGACACAGCGTGGGACCGGCGG
GCCCTGGCTCGCAGGGTACGCTGCTGCGGCCACTGAGCCGGGTGGAATCCTACCTATGCCCTATGTC
ACTGAGGCCACCCGAGTGCAGCTGGTGTGCCACTCCTGGTGGCTGAAGCTGCTGCAGCCCCGGCTTTC
CTCGAGGCCTTTGCAGCCAATGTCTGGAGCCACGAGAACATGCATTGCTCACCCGTGTTGCTGGTCTAC
GGGCCACGAGAAGGTGGCCGTGGAGCTCCAGACCCATTTCTTGGGGTGAAGGCTGCAGCAGCGGAGTTA
GAGCGACGGTACCCTGGGACGAGGCTGGCTGGCTCGCTGTGCGAGCAGAGGCCCTTCCAGGTGCGA
TCATGGACGTGGTCTCGAAGAAGCACCCCTGTGGACACTCTTCTTCTTACCACCGTGTGGACAAGG
CCTGGGCCGAAGTCTCAACCGCTGTCGATGAATGCCATCTCTGGCTGGCAGGCTTCTTTCCAGTC
CATTTCCAGGAGTTCAATCCTGCCTGTACCACAGAGTACACCCAGGGCCCCCGGGGCTGGCCCT
GACCCCCCTCCCTCCTGGTGTGACCCCTCCCGGGGGCTCCTATAGGGGGGAGATTTGACCGGCAG
GCTTCTGCGGAGGGCTGTTCTACAACGCTGACTACCTGGCGGCCGAGCCCGGCTGGCAGGTGAACTG
GCAGGCCAGGAAGAGGAGGAAGCCCTGGAGGGGCTGGAGGTGATGGATGTTTTCTCCGGTTCTCAGGG
CTCCACCTCTTTCGGGCCGTAGAGCCAGGGCTGGTGCAGAAGTTCTCCTGCGAGACTGCAGCCCAGG
CTCAGTGAAGAACTCTACCACCGCTGCCCTCAGCAACCTGGAGGGGCTAGGGGGCCGTGCCAGCTG
GCTATGGCTCTCTTTGAGCAGGAGCAGGCCAATAGCACTTAG
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
  
```

Restriction Sites: SgfI-MluI

ACCN: NM_019015

Insert Size: 2319 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none"> 1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_019015.2</u>
RefSeq Size:	4001 bp
RefSeq ORF:	2319 bp
Locus ID:	54480
UniProt ID:	<u>Q9P2E5</u>
Cytogenetics:	7q36.1
Protein Pathways:	Chondroitin sulfate biosynthesis, Metabolic pathways
MW:	85.9 kDa
Gene Summary:	<p>Transfers glucuronic acid (GlcUA) from UDP-GlcUA to N-acetylgalactosamine residues on the non-reducing end of the elongating chondroitin polymer. Has no N-acetylgalactosaminyltransferase activity.[UniProtKB/Swiss-Prot Function]</p> <p>Transcript Variant: This variant (1) represents the longer transcript and encodes the longer isoform (1). Sequence Note: This RefSeq record was created from transcript and genomic sequence data to make the sequence consistent with the reference genome assembly. The genomic coordinates used for the transcript record were based on transcript alignments.</p>