

Product datasheet for **SC304671**

UGT1A3 (NM_019093) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	UGT1A3 (NM_019093) Human Untagged Clone
Tag:	Tag Free
Symbol:	UGT1A3
Synonyms:	UDPGT; UDPGT 1-3; UGT-1C; UGT1-03; UGT1.3; UGT1A3S; UGT1C
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

Fully Sequenced ORF: >SC304671 representing NM_019093.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)


```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGGCCACAGGACTCCAGGTTCCCCTGCCGTGGCTGCCACAGGACTGCTGCTTCTCCTCAGTGTCCAG
CCCTGGGCTGAGAGTGGAAAGGTGTTGGTGGTGCCCATTTGATGGCAGCCACTGGCTCAGCATGCGGGAG
GTCTTGGCGGAGCTCCATGCCAGAGGCCACCAGGCAGTGGTCTCACCCAGAGGTGAATATGCACATC
AAAGAAGAGAACTTTTTACCCTGACAACCTATGCCATTTCTGGACCCAGGATGAATTTGATCGCCAT
GTGCTGGGCCACACTCAACTGTACTTTGAAACAGAACATTTTCTGAAGAAATTTTTCAGAAGTATGGCA
ATGTTGAACAATATGTCTTTGGTCTATCATAGGTCTTGTGTGGAGCTACTACATAATGAGGCCCTGATC
AGGCACCTGAATGCTACTTCTTTGATGTGTTTTAACAGACCCCGTTAACCTCTGCGCGGCAGTGTG
GCTAAGTACCTGTCGATTCTACTGTGTTTTTTGAGGAACATTCATGTGATTTAGACTTTAAGGGC
ACACAGTGTCAAACCCCTCTCTATATTCCTAGATTACTAACAAACCAATTCAGACCACATGACATTC
ATGCAAAGGGTCAAGAACATGCTCTACCCTCTGGCCCTGTCTACATTTGCCATGCTTTTTCTGCTCCT
TATGCAAAGCCTTGCCCTGAGCTTTTTTCAGAGAGAGGTGTGAGTGGTGGATATTCTCAGTCATGCATCT
GTGTGGCTGTTCCGAGGGGACTTTGTGATGGACTACCCAGGCCAATCATGCCAACATGGTCTTCATT
GGGGGCATCAACTGTGCCAACAGGAAGCCACTATCTCAGGAATTTGAAGCCTACATTAATGCTTCTGGA
GAACATGGAATTTGGTTTTCTTTGGGATCAATGGTCTCAGAAATTCAGAGAAGAAAGCTATGGCA
ATTGCTGATGCTTTGGGCAAAATCCCTCAGACAGTCTGTGGCGGTACACTGGAACCCGACCATCGAAT
CTTGCGAACAACACGATACTTGTAAAGTGGTACCCCAAAACGATCTGCTTGGTCACCCGATGACCCGT
GCCTTTATCACCCATGCTGGTCCCATGGTGTATGAAAGCATATGCAATGGCGTTCCCATGGTGATG
ATGCCCTTGTGGTGGTATCAGATGGACAATGCAAAGCGCATGGAGACTAAGGGAGCTGGAGTACCCTG
AATGTTCTGGAATGACTTCTGAAGATTTAGAAAATGCTCTAAAAGCAGTCATCAATGACAAAAGTTAC
AAGGAGAACATCATGCGCCTCTCCAGCCTTACAAAGACCGCCCGGTGGAGCCGCTGGACCTGGCCGTG
TTCTGGGTGGAGTTGTGATGAGGCACAAGGGCGGCCACACCTGCGCCCGCAGCCACGACCTCACC
TGGTACCAGTACCATTCTTGGACGTGATTGGTTTCTCTTGGCCGTCGTGCTGACAGTGGCCTTCATC
ACCTTTAAATGTTGTGCTTATGGCTACCGGAAATGCTTGGGGAAAAAAGGGCGAGTTAAGAAAGCCAC
AAATCCAAGACCCATTGA
ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGCGC
  
```

Restriction Sites:

Sgfl-MluI

Plasmid Map:



ACCN:	NM_019093
Insert Size:	1605 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none"> 1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_019093.2</u>
RefSeq Size:	2345 bp
RefSeq ORF:	1605 bp
Locus ID:	54659
UniProt ID:	<u>P35503</u>
Cytogenetics:	2q37.1
Protein Families:	Transmembrane
Protein Pathways:	Androgen and estrogen metabolism, Ascorbate and aldarate metabolism, Drug metabolism - cytochrome P450, Drug metabolism - other enzymes, Metabolic pathways, Metabolism of xenobiotics by cytochrome P450, Pentose and glucuronate interconversions, Porphyrin and chlorophyll metabolism, Retinol metabolism, Starch and sucrose metabolism
MW:	60.3 kDa

Gene Summary:

This gene encodes a UDP-glucuronosyltransferase, an enzyme of the glucuronidation pathway that transforms small lipophilic molecules, such as steroids, bilirubin, hormones, and drugs, into water-soluble, excretable metabolites. This gene is part of a complex locus that encodes several UDP-glucuronosyltransferases. The locus includes thirteen unique alternate first exons followed by four common exons. Four of the alternate first exons are considered pseudogenes. Each of the remaining nine 5' exons may be spliced to the four common exons, resulting in nine proteins with different N-termini and identical C-termini. Each first exon encodes the substrate binding site, and is regulated by its own promoter. Substrates of this enzyme include estrone, 2-hydroxyestrone, and metabolites of benzo alpha-pyrene. [provided by RefSeq, Jul 2008]