

Product datasheet for **SC304683**

IRGC (NM_019612) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	IRGC (NM_019612) Human Untagged Clone
Tag:	Tag Free
Symbol:	IRGC
Synonyms:	CINEMA; IFGGE; ligp5; IRGC1; R30953_1
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

Fully Sequenced ORF: >SC304683 representing NM_019612.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAATTTTGAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCACGATCGCC
ATGGCTACTTCAAAGTTGCCCGTGGTGCCTGGGAGGAGGAAAACACCATCCTTATGGCCAAGGAAAGG
CTGGAGGCCCTGCGCACAGCCTTTGAGTCGGGTGACCTCCCCAGGCCGCTCTACCTCCAGGAGCTG
CTGGCCTCCACGAAAGCATCCGCCTGGAGGTGGGCGTACGGGCGAGTCGGGCGCGGCAAGTCTCTC
CTCATCAATGCCCTGCGTGGCCTGGAGGCCGAGGACCCTGGCGCGGCTCTCACGGGCGTCATGGAGACC
ACGATGAACCGTCGCCCTATCCACACCCACAGTTCCTGACGTGACCCTCTGGGACCTGCCAGGAGCC
GGCTCTCCAGGCTGCCCGGCTGACAAGTACCTAAAGCAGGTAGACTTCAGCCGCTATGACTTCTTCTG
CTGGTCTCCCCCGCGCTGCGGGGCGCTCGAGACCCGCTGGCCGCTGAGATCCTGTGCCAGGGCAAG
AAGTTCTACTTTGTGCGACCAAGGTGGACGAGGACCTGGCGGCCACGCGACCCAGCGGCCGTCGGC
TTCAGAGAGGCCGCTGTCTGCAGGAGATCCGAGACCACTGTGCCGAGCGGCTGCGGGAGGCCGGCGTG
GCTGACCCTCGCATCTTCTGGTGTCAACCTCTCGCCGCGCCGCTACGACTTTCCACGCTGGTGTCC
ACCTGGGAGCACGACCTGCCCTCCACCGGCGCCACGCTGGCCTGCTGTGCTCCCCGACATCTCGCTG
GAGGCCTTGCAGAAGAAGAAGGCCATGCTTCAAGAGCAAGTCTCAAGACCGCCCTGGTGTGGGCGTC
ATCCAGGCCCTGCCGGTCCCAGGGCTGGCGGCCGCTACGATGATGCGTTGCTCATCCACTCACTGCGT
GGCTACCACCGCAGCTTTGGTCTGGACGACGACTCGCTGGCCAAGCTGGCCGAGCAGGTGGGCAACAG
GCAGGTGACCTGCGCTCGGTATCCGCTCCCCACTGGCCAACGAGGTCTCGCCTGAGACTGTCTGCGG
CTCTATTCAGTCGTCGACGCGCCATGCGGGTGGCCCGCGCCTTTGAGAGGGGATCCCTGTGTTT
GGGACGCTGGTGGCTGGCGGCATCAGCTTTGGCGCTGTCTACACCATGCTCCAGGGCTGCCTCAACGAG
ATGGCTGAGGACGCCAGCGTGTCCGATCAAGGCCCTGGAGGATGACGAGCCGACGCCGAGGTGAGC
TTGGAAGTGGCCAGTGACAATGGCGTGGAAAAGGGGGCTCCGGGAGGGAGGTGGGAGGAAGCCCA
CTCTAACCTGCAGGAAGCTCGGCCTCCTTCTAAGTACATTCTGGACAGCTGGAAGAAACAGACTCA
GAAGAGAAATAA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGCGC
  
```

Restriction Sites: SgfI-MluI

ACCN: NM_019612

Insert Size: 1392 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_019612.3](#)

RefSeq Size: 1619 bp

RefSeq ORF: 1392 bp

Locus ID: 56269

UniProt ID: [Q6NXR0](#)

Cytogenetics: 19q13.31

MW: 50.3 kDa